



•

LOCATIONS

Aidacare South Australia

ADELAIDE

71 Port Road Thebarton SA 5031 TEL: (08) 8118 5600 FAX: (08) 8352 8799



5

P 1300 133 120 E des@aidacare.com.au W www.aidacare.com.au

2

DES Enquiries:

CONTENTS

Bathroom & toilet	7 (7
Bedroom	12
Chairs - Therapeutic	13 (j)
Ramps	15 J
Crutches & Walkers	16 🝞
Wheelchair & Seating	17 (7)
Mobility Scooters	29



Aidacare & DES

Aidacare, in conjunction with our supply partners Pride, Permobil, Sunrise Medical, OttoBock and Invacare are proud to be recognized suppliers of a large number of products listed on the current DES (Domiciliary Equipment Service) contract that provides equipment and home modification services assisting people to be independent at home and in the community.

Each product has a specific order form which must be completed and submitted with the application. Along with a prescription detailing each person's particular needs.

In addition to the items in this catalogue, Aidacare is also able to supply other subsidised equipment including:

- Special purpose wheelchairs
- Pressure cushions and specialised seating components
- Mobile and static commodes
- Pressure care mattresses
- Other bedroom, bathroom and Daily Living equipment.

Our electronic catalogue is available for download at www.aidacare.com.au

For the full range of our equipment please browse the complete online Aidacare Catalogue at www.aidacare.com.au.

For additional information please contact our friendly staff on: 1300 133 120, specialised equipment professionals

OUR SERVICES

Aidacare offer a broad range of services to DES prescribers and clients including:

- Access to an extensive range of demonstration stock for trial
- Specialist sales staff to assist with trials and support prescription of customised and specialised equipment
- Clinical specialists to provide additional advice and support for healthcare
- Clinical education and support
- Service & repair field technicians
- General customer service and support including the booking of equipment trials

OUR CASE CONSULTING SPECIALISTS

Helen Smith

Helen is an Occupational Therapist with many years of experience in the field of Mobility and Seating. She is passionate about achieving accurate seating and positioning outcomes for clients and has been involved in the prescription of both powered and manual mobility and complex seating here and abroad. In her role at Aidacare, Helen will be working with both adult and children in wheelchair and seating prescriptions and specialised equipment solutions.

Michael Mifsud

Michael is an experienced seating and wheelchair technician who specialises in Complex Rehab equipment. Michael is skilled in the fabrication of customised seating (foam on ply) and the setup of complex seating systems whether it be modular or customised. He has been involved in power mobility programming and electronics and the repair and diagnostics of power mobility electronics.

Linh Yeomans

Linh is a Physiotherapist with many years of experience in the prescription of manual and powered mobility and customised complex and modular seating systems for both adults and children. Linh has also been involved in the prescription of paediatric mobility and positioning and general specialised paediatric equipment such as standing frames, strollers, walkers and sleep systems. Linh manages the Complex Rehab department at Aidacare SA, however, also carries out wheelchair and seating assessments for both adults and children. Helen Smith Equipment Specialist - Complex Rehab

Occupational Therapist Mobile: 0419 497 298 Email: helen.smth@aidacare.com.au

Linh Yeomans

Manager Complex Rehab Physiotherapist (B.Physio) Mobile: 0409 643 098 Email: linh.yeomans@aidacare.com.au

Michael Mifsud Service Techniciam Wheelchair Techinican Mobile: (08) 8118 5600





DES Order & Approval Process

Aidacare offers a service to Prescriber's that are based on "Partnering to provide the best solution". We listen, respond and work with therapists to provide the highest possible customer service solution for the client.



DES Catalogue

Referral is received by Prescriber

- Prescriber carries out an assessment to identify client's needs and DES funding eligibility
- Together the client and prescriber consider possible equipment solutions or modifications required

If client is approved for DES funding, the prescriber contacts DES for product/equipment availability. If no suitable product match is found, the prescriber contacts Aidacare to arrange an in home trial or the client is seen in the Aidacare seating clinic. Trials are only started once funding has been approved



www.aidacare.com.au

Aidacare Admin:

P 1300 133 120 F 1300 787 052 E des@aidacare.com.au

SOLUTION

Once possible equipment choices have been trialled, the client along with the prescriber and the Aidacare Complex Rehab Specialist will decide what equipment will work best to meet the client needs

ORDER STEP 1

Aidacare provides the prescriber with detailed product information and a formal quote

ORDER STEP 2

Prescriber forwards to DES: • Aidacare quote

Product specifications

Appropriate prescriber application forms

DES APPROVAL GRANTED

DELIVERY

When all equipment is ready for delivery, an Aidacare Complex Rehab team member will contact the prescriber. Depending on what equipment or modifications are supplied, these will need to either go to the DES warehouse for 'barcoding' or the work (such as wheelchair and seating modifications) may be carried out or supplied directly to the client in consultation with the prescriber as required.



DES Enquiries:



CUSTOMISED AND SPECIALISED EQUIPMENT



CUSTOMISED SEATING & MOBILITY SOLUTIONS

Aidacare Complex Rehab Specialists are dedicated to provide specialised advice, service and equipment based solutions for our clients.

This specialised field is supported by our experienced team of specialists, Customer Care and Technical service staff with practical hands on experience and professional training including qualified Therapists and Rehab Engineers.

The key focus areas of our Complex Rehab team includes:

- Scripted Manual Wheelchairs
- Scripted Power Wheelchairs
- Seating & Positioning Products
- Pressure Management Cushions

The Aidacare Complex Rehab Specialists in conjunction with therapists assess and recommend the clinically appropriate product solution based on the individual client needs. We are not limited to a specific range or style, enabling us to recommend the right option. Aidacare represents the full scope of leading international industry brands and models.





DES Enquiries:



HAND HELD SHOWER HOSE PUSH ON SINGLE

- Flexible shower head with single push on attachment
- Gentle flow shower head for sensitive skin
- High pressure clamp for attaching to tap

Code	Description	Length	Warranty
BTS114500	Hand Held Shower Push On Single	2 m	12 Months



SHOWER COMMODE TRANSIT

- Attendant Propelled with Swing-up Arms and Sliding Foot Plate
- 4 Stainless Steel Locking Castors
- Lower Height at seat front
- Stainless Steel Frame
- Available in Open and Closed Seats

Code	Description	Frame Width	Seat Width	Seat Depth	Seat Height	Max User Weight
BTC070600SP1	Shower Commode Transit - Low (Open Seat)	510 mm	440 mm	430 mm	500 mm	120 kg
BTC070600SP2	Shower Commode Transit - Low (Closed Seat)	510 mm	440 mm	430 mm	500 mm	120 kg
BTC070600SP3	Shower Commode Transit - Std (Open Seat)	510 mm	440 mm	430 mm	590 mm	120 kg
BTC070600SP4	Shower Commode Transit - Std (Closed Seat)	510 mm	440 mm	430 mm	590 mm	120 kg



SHOWER COMMODE SELF PROPELLING

- Self Propelled Rear Wheels with Swing-up Arms and Swing Away Foot Rests
- 2 Stainless Steel Locking Castors
- Stainless Steel Frame for Corrosion Resistances
- Includes Bed Pan Carrier for K-Care Shower commodes
- Available in Open and Closed Seats

Code	Description	Frame Width	Seat Width	Seat Depth	Seat Height	Max User Weight
BTC077900SP1	Shower Commode Self Propelled - Std (Open Seat)	510 mm	440 mm	430 mm	590 mm	120 kg
BTC077900SP2	Shower Commode Self Propelled - Std (Closed Seat)	510 mm	440 mm	430 mm	590 mm	120 kg
BTC077900SP3	Shower Commode Self Propelled - Low (Open Seat)	510 mm	440 mm	430 mm	500 mm	120 kg
BTC077900SP4	Shower Commode Self Propelled - Low (Closed Seat)	510 mm	440 mm	430 mm	500 mm	120 kg
BTC080500SP1	Shower Commode Folding (Open Seat)	510 mm	440 mm	430 mm	590 mm	120 kg
BTC080500SP2	Shower Commode Folding (Closed Seat)	510 mm	440 mm	430 mm	590 mm	120 kg





SHOWER COMMODE SWING AWAY FOOTPLATES

- Attendant Propelled with Swing -up Arms and Swing-away Foot Rests
- 4 Stainless Steel Locking Castors
- Stainless Steel Frame
- Includes: Bed Pan Carrier, Stainless Steel
- Available in Open and Closed Seats

Code	Description	Frame Width	Seat Width	Seat Depth	Seat Height	Max User Weight
BTC070900SP1	Shower Commode - Swing Away (Open Seat)	510 mm	440 mm	430 mm	590 mm	120 kg
BTC070900SP2	Shower Commode - Swing Away (Closed Seat)	510 mm	440 mm	430 mm	590 mm	120 kg



SHOWER COMMODE TRANSIT FOLDING

- Attendant Propelled with Swing-up Arms and Sliding Foot Plate
- 4 Stainless Steel Locking Castors
- Lower Height at seat front
- Stainless Steel Frame
- Includes: Bed Pan Carrier, Stainless Steel
- Available in Open and Closed Seats

Code	Description	Frame Width	Seat Width	Seat Depth	Seat Height	Max User Weight
BTC071400SP1	Shower Commode Transit Folding (Open Seat)	510mm	440mm	430mm	590mm	120 kg
BTC071400SP2	Shower Commode Transit Folding (Closed Seat)	510mm	440mm	430mm	590mm	120 kg



SHOWER COMMODE BARIATRIC

- Attendant Propelled with Padded Back Rest
- Removable Padded Arms and Swing-away Footrest
- 4 Stainless Steel Frame
- Includes: Bed Pan Carrier, Stainless Steel
- Available in Hand Made Open and Closed Seats

Code	Description	Frame Width	Seat Width	Seat Depth	Seat Height	Max User Weight.
BTC062100SP1	Shower Commode Transit Bariatric (Open Seat)	610 mm	550 mm	440 mm	590 mm	190 kg
BTC062100SP2	Shower Commode Transit Bariatric (Closed Seat)	610 mm	550 mm	440 mm	590 mm	190 kg
BTC062200SP1	Shower Commode Bariatric (Open Seat)	660 mm	600 mm	450 mm	590 mm	210 kg
BTC062200SP2	Shower Commode Bariatric (Closed Seat)	660 mm	600 mm	450 mm	590 mm	210 kg





SHOWER COMMODE TRANSIT TILT IN SPACE

- Attendant Propelled
- Gas Lift Tilt-in-Space with Swing-Up Arms
- Padded Calf Rest and Swing-away Foot Rests
- Removable Padded Arms and Swing-away Footrest
- 4 Stainless Steel Frame
- Includes Bed Pan Carrier, Stainless Steel
- Available in Hand Made Open and Closed Seats

Code	Description	Frame Width	Seat Width	Seat Depth	Seat Height	Max User Weight
BTC071100SP1	Shower Commode Transit Tilt in Space (Open Front)	510 mm	440 mm	430 mm	590 mm	135 kg
BTC071100SP2	Shower Commode Transit Tilt in Space (Closed Front)	510 mm	440 mm	430 mm	590 mm	135 kg



SHOWER CHAIR WITH ADJUSTABLE ARMS

- Adjustable Height with Clip-On Plastic Seat and Arms
- Zinc Plated Frame
- Fitted with Non Slip Rubber Tips
- Shower Chair Arms Assist Transfers

Code	Description	Frame Width	Seat Width	Seat Depth	Seat Height	Max User Weight
BTS114900	Shower Chair with Adjustable Arms	450 mm	440 mm	430 mm	420-550 mm	110 kg



SHOWER CHAIR BARIATRIC

- Reinforced wide frame designed for bariatric use
- Fixed Height frame for extra stability
- Fixed with non-slip rubber tips
- Padded seat, armrests and backrest for user comfort and support
- Zinc treated anti corrosive frame
- Warranty: 12 months

Code	Description	Frame Width	Seat Width	Seat Depth	Seat Height	Max User Weight
BTS105700	Shower Chair Bariatric	680 mm	550 mm	430 mm	520 mm	300 kg



DES Enquirie



SHOWER STOOL WITH ARMS

- Shower stool arms assist transfers
- Fully adjustable seat height and angle to suit user
- Fitted with non-slip rubber tips
- Moulded clip on seat for quick assembly, cleaning or replacement
- Frame options Zinc treated (anti corrosive) or Aluminum (non-corrosive)
- Warranty: 12 months

Code	Description	Frame Width	Seat Width	Seat Depth	Seat Height	Max User Weight
BTS115200	Shower Stool with arms	460 mm	350 mm	300 mm	490-610 mm	110 kg



OVER TOILET AID EXTRA WIDE

- Extra wide and reinforced frame accommodating the larger user with seat width of 500 mm
- Fully adjustable seat height to suit user
- Fitted with non-slip rubber tips
- Removable lift back seat for easy cleaning, comes standard with splash guard
- Zinc treated anti corrosive frame
- Warranty: 12 months

Code	Description	Seat Width	Seat Depth	Seat Height	Max User Weight
BTT146300	Over Toilet Aid - Extra wide	500mm	420mm	480-620mm	135kg



OVER TOILET AID ADJUSTABLE

- Heavy Duty
- Adjustable Height with Clip-on Plastic Seat Flap and Splashguard
- Zinc Plated Frame

Code	Description	Seat Width	Seat Depth	Seat Height	Max User Weight
BTT153200	Over Toilet Aid Adjustable	445 mm	420 mm	480-620 mm	110 kg







RAISED TOILET SEAT WITH OPEN LID

- Provides a raised seat and support to assist with toileting and transfers
- Fully adjustable seat height to suit user
- Fitted with non-slip rubber tips
- Removeable lift back seat for easy cleaning
- A standard splash guard provided
- Frame options Zinc treated (anti corrosive) or Aluminium (non-corrosive)

Code	Description	Height	Warranty	Max User Weight
BTT149800	Raised Toilet Seat with Open lid 50	50 mm	12 Months	190 kg
BTT150200	Raised Toilet Seat with Open lid 100	100 mm	12 Months	190 kg





DES Enquiries:



BEDSIDE COMMODE CHAIR

- Heavy Duty
- Adjustable Height
- Lift-up Padded Seat
- Padded Back
- Padded Arm Rest

- Lift-up Toilet Seat
- Fixed Bowl Carrier
- Autoclavable Bowl with carry handle and lid
- Grey Frame and Slate Upholstery

Code	Description	Colour	Seat Width	Seat Depth	Seat Height Range	Max User Weight
BEC046070	Bedside Commode Chair	Champagne / Vanilla Frame	440 mm	440 mm	480-605 mm	110 kg



OVER BED TABLE

- Laminate Top and Frame
- "C" Base Thermoform Top
- Spring Loaded
- Tilt Top ability

Code	Description	Top Length	Top Width	Height Range
BEA006700	Over Bed Table	855 mm	400 mm	710-1050 mm





CHAIRS



CHAIR AIR BED DELUXE

- Air Comfort Deluxe V2 with Seat Tilt
- Back Recline and Raising Leg Rest
- Full (Four) adjustable Air Bladder System
- Suitable for client 162cm 188cm tall
- Covered in Multi-Stretch PU and PVC material that is Anti-Microbial and Fire Retardant

Code	Description	Seat Width	Back Width	Back Height	Max. User Weight
CHP198300	Chair Air Bed Deluxe Small	460 mm	450 mm	660 mm	90 kg
CHP198200*	Chair Air Bed Deluxe Large	510 mm	500 mm	760 mm	180 kg



SIT TO STAND RECLINER

- Designed to be the ultimate in lift chair, comfort, style and performance
- State of the art lift system will recline further back that any other lift chair
- Equipped with emergency backup system should the power fail

Code	Description	Colour	Seat Width	Back Width	Back Height	Max. User Weight
CHP197102	Power Recliner Pride C5	Chocolate	560 mm	710 mm	660-710 mm	170 kg
CHP197104	Power Recliner Pride C5	Blue	560 mm	710 mm	660-710 mm	170 kg



KITCHEN STOOL ADJUSTABLE WITH ARMS

- Black Frame
- Ideal for the use in the kitchen or similar work area where a standard chair is not suitable
- Fully adjustable seat height to suit individual user requirements
- Padded seat and backrest for added comfort

Code	Description	Seat Width	Seat Depth	Back Rest Height		Max User Weight
CHP213500	Kitchen Stool with Adjustable Arms	360 mm	300 mm	400 mm	520-640mm	100 kg



DES Enquiries:



HI LITE CHAIR HIGH BACK

- High back lightweight orthopedic chair
- Lightweight and sturdy design for function and comfort
- Contoured passed backrest for lumbar and neck support
- Fully adjustable seat height and angle
- Rubbers stoppers for stability
- Easy to clean vinyl upholstery

Code	Description	Colour	Seat Width	Seat Depth	Backrest Height	Seat to Floor	Max User Weight
CHP204900	Hi Lite Chair Geriatric	Vanilla Frame/ Champagne	490 mm	480 mm	700 mm	430 - 580 mm	120 kg
CHP204500	Hi Lite Chair Geriatric	Grey Frame/ Slate	490 mm	480 mm	700mm	430 - 580 mm	120 kg
CHP205000	Hi Lite Chair Geriatric	Vanilla Frame/ Mulberry	490 mm	480 mm	700mm	430 - 580 mm	120 kg



LOW BACK

- Low back multifunction orthopedic chair
- Fully adjustable seat height
- Curved padded backrest provides postural support
- Rubber stoppers for stability
- Easy to clean vinyl upholstery
- Powered coasted steel frame

Code	Description	Colour	Seat Width	Back Rest Height	Seat to Floor	Max User Weight
CHP208500	Utility Chair Geriatric LB	Vanilla Frame/ Champagne	440 mm	450 mm	430 - 580 mm	120 kg
CHP207900	Utility Chair Geriatric LB	Grey Frame/ Slate	440 mm	450 mm	430 - 580 mm	120 kg



RAMPS



DEPAC RAMP WITH EDGE

- Fiberglass ramps with concertina folding action & transport
- 5cm Edge Barrier designed for safety & security to prevent steering off the edge of the ramp

Code	Description	Length	Width	Weight	Safe Working Load
RAM656120	Depac Ramp with Edge	2000 mm	820 mm	13.7 kg	300 kg



PERSONAL & SENIOR DEPAC RAMP

- Fiberglass ramp with concertina folding action for storage & transport
- The Personal is lightweight portable ramp, allows wasy access in domestic situaitions including climbing curbs. Features a strap for storage on back of the wheelchair and folds up to allow for storage and transporations
- The senior is remarkably lightweight and strong and folds easy to assist in storage and transportation. Available in a large range of size.

Code	Description	Length	Width	Maximum Height	Weight	Safe Working Load
RAM656200	Personal Depac Ramp	700 mm	745 mm	175 mm	5 kg	300 kg
RAM656360	Senior Depac Ramp	900 mm	750 mm	220 mm	8 kg	300 kg
RAM656310	Senior Depac Ramp	1200 mm	750 mm	300 mm	8 kg	300 kg
RAM656320	Senior Depac Ramp	1350 mm	750 mm	340 mm	8.5 kg	300 kg
RAM656330	Senior Depac Ramp	1650 mm	750 mm	412 mm	11 kg	300 kg







AXILLARY CRUTCH

- Patented clip-lock adjustable handle and push button foot pieces for fast, precise measuring and fitting
- Lightweight aluminium frame
- High density rubber cushions and hand grips providing long wear and greater comfort

Code	Description	User Height	Crutch Height	Ferrule Tip	Max User Weight
WAC692200	Axillary Crutch - Adult / Medium	1550 - 1755 mm	1120 - 1320 mm	22 mm	136 kg
WAC692400	Axillary Crutch - Adult / Large	1780 - 1980 mm	1320 - 1520 mm	22 mm	136 kg



ALUMINIUM FOLDING WALKING FRAME

- Provides assistance for those who have difficulty walking
- Folds flat for easy transportation and storage
- Lightweight aluminium frame ideal for indoor use
- Height adjustable to suit user requirements
- Supplied with four rubber tips

Code	Description	Width	Handle Height	Ferrule Tip	Warranty	Product Weight	Max User Weight
WAF706700	Aluminium Folding Walking Frame	610 mm	710-787mm	25 mm	12 Months	2 kg	127 kg



SEAT WALKER WITH HAND BRAKES

- Provides assistance to those who need confidence walking and require rest en route
- Folds away for easy storage or transportation
- Height adjustable handles to suit user
- Handbrakes are lockable for stability
- Padded backrest for seating comfort
- Lightweight aluminium construction

Code	Description	Width Between arms	Seat Height	Handle Height	Max User Weight
WAF705800	Seat Walker - Small Frame 7" Wheels	435 mm	475 mm	740-840 mm	120 kg



- 2 1300 133 120 5 des@aidacare.com.au
- W www.aidacare.com.au



BREEZY BASIX

- Folding lightweight chair
- Available in self-propelled or transit
- Adjustable push handle height (+/- 3")
- Available with folding backrest
- Tension adjustable upholstery on backrest
- Swing away legrests as standard
- Drum break option available (16, 18 & 20")
- 41 46 cm adjustable seat depth
- Castor angle adjustment 0, 3 & 6 degrees

Code	Description	Colour	Seat Width	Seat Depth	Back Height	Frame Weight	Product Weight	Max. User Weight
MWS449101	Basix 14"	Sliver	360 mm	410 - 460 cm	410 - 435 mm	8.4 kg	14.9 kg	125 kg
MWS449106	Basix 16"	Sliver	400 mm	410 - 460 cm	410 - 435 mm	8.4 kg	14.9 kg	125 kg
MWS449113	Basix 18"	Sliver	460 mm	410 - 460 cm	410 - 435 mm	8.4 kg	14.9 kg	125 kg
MWS449120	Basix 20"	Sliver	520 mm	410 - 460 cm	410 - 435 mm	8.4 kg	14.9 kg	125 kg



OSD ASPIRE EVOKE LIGHTWEIGHT

- Lightweight aluminium frame (15kg) for ease of transportation
- Quick release rear wheels to reduce weight of chair for storage and transportation
- Rear axle adjustability to vary the centre of gravity to suit amputees
- Angle adjustable front castors to counter-balance adjustments in rear wheel mount location
- Height adjustable attendant push-handles to better suit attendants/carers
- Angle & tension adjustable backrest & depth adjustable seating for a more customised seating position and personalised comfort
- Angle & Height adjustable backrest canes for maximum postural support & comfort
- Swing-in/out height adjustable & removable footrests provide accurate postural support

Code	Description	Colour	Seat Width	Seat Depth	Back Height	Warranty	Frame Weight	Product Weight	Max User Weight
MWS440100	W/ Chair Folding - OSD Evoke LW 350	Silver	350 mm	400 mm	400 mm	24 months	7.4 kg	14.7 kg	130 kg
MWS440200	W/ Chair Folding - OSD Evoke LW 400	Silver	400 mm	400 mm	400 mm	24 months	7.6 kg	14.7 kg	130 kg
MWS440300	W/ Chair Folding - OSD Evoke LW 450	Silver	450 mm	400 mm	400 mm	24 months	7.9 kg	14.7 kg	130 kg
MWS440500	W/ Chair Folding - OSD Evoke LW 500	Silver	500 mm	400 mm	400 mm	24 months	8.1 kg	14.7 kg	130 kg

Accessories	
Code	Description
WCA858370	Elevating Legrest - Black Left
WCA858390	Elevating Legrest - Black Right
WCA858440	Stump Support - Black Left
WCA858450	Stump Support - Black Right
MSP42102	Armrests (Height & Depth Adj.) - Pair
WCA858460	Anti Tip Bars





OTTO BOCK START INTRO M1

- Aluminum tube frame
- Tip-assist integrated in the frame
- Individual footrests, angle adjustable, plastic footplate
- Swing-away footrest holders
- Desk side panel with armrest, different

height settings possible

- Seat and back upholstery of black nylon
- Rear wheels with quick-release axles
- Wheel lock
- Frame colour: silver metallic
- 24" rear wheel with PU tyre

Code	Description	Colour	Seat Width	Seat Depth	Back Height	Frame Weight	Product Weight	Max. User Weight
MWC427350	Start 45.5	Silver	455 mm	400 - 460 mm	425 mm	10.5 kg	15.5 kg	125 kg
MWC427340	Start 40.5	Silver	405 mm	400 - 460 mm	425 mm	10.5 kg	15.5 kg	125 kg
MWC427352	Start 50.5	Silver	505 mm	400 - 460 mm	425 mm	10.5 kg	15.5 kg	125 kg



ACTION 3 JUNIOR

- Lightweight, foldable paediatric wheelchair designed for children aged between 3 and 15 years
- Matches the individual needs of the child and can grow as they grow
- Large range of options to accompany the child in their development and to

provide the necessary clinical support

- All models of Action 3 junior come in attractive, modern styles, designed to appeal to kids
- Comes in Self Propelled and Attendant Propelled versions

Description	Colour	Seat Width	Seat Depth	Back Height	Warranty	Max. User Weight
Action 3 Junior	Dark Red, Sapphire Blue	205 – 380 mm	200 – 350 mm	350 - 430 mm	24 months	80 kg







OTTO BOCK START M6 JUNIOR

The Start Junior is a folding wheelchair just for kids. It features swing in / swing out legrests that assist in transfers as well as large 20" rear wheels to encourage self propelling. This chair is customizable with a large number of options available including elevating legrests, anti-tip bars and a tray.

- Aluminium tube frame
- All main component groups are attached with screws
- Standardized caster forks for all wheel sizes
- Durable upholstery

- Small folding size with folding fixture
- Children's colours
- Optimized wheel lock
- Footplate for short lower leg lengths

Description	Seat Width	Seat Depth	Back Height	Warranty	Product Weight	Max. User Weight
Start M6 Junior	280-380 mm	310 - 370 mm	300 - 375 mm	24 months	14.5 kg	90 kg



LIFESTREAM XF

The Litestream XF delivers a high level of customisability in an economic lightweight folding package. With a unique rigid reinforced frame, adjustable rear axle mount and dozens of performance enhancing features and options, chair setup can be optimised to meet the unique requirements of the individual user. A Heavy Duty package is available in frame widths 20, 21 & 22".

- Customisable frame widths 12-22"
- Adjustable rear axle and castor forks
- Height adjustable backrest
- Swingaway tubular armrests
- 70-degree swing-in/swing-out legrests
- 24" quick release spoke rear wheels with pneumatic tyres
- 8" polyurethane castors

Description	Seat Width	Seat Depth	Back Height	Warranty	Product Weight	Max. User Weight
Litestream XF	355 - 530 mm	355- 508 mm	305-480 mm	24 months	12 kg	113 kg / 159 kg (HD Model)



LITESTREAM JUNIOR

From the classroom to the playground, the Litestream® Junior is a superb manual wheelchair for even the most active child. The Litestream Junior accepts a wide range of accessories and options making it a highly adaptable and versatile rehab solution.

Description	Colour	Seat Width	Seat Depth	Warranty	Frame Weight	Product Weight	Max. User Weight
Litestream Short Frame	Blue	254 mm	254 mm	24 months	7.4 kg	10.66 kg	56.7 kg
Litestream Medium Frame	Blue	330 mm	330 mm	24 months	7.6 kg	10.66 kg	68.04 kg
Litestream Long Frame	Silver	406.4mm	431.8 mm	24 months	7.6 kg	10.66 kg	79.4 kg







CROSSFIRE T6

The Invacare® Top End® Crossfire™ Series of everyday wheelchairs offers a range of frame materials and frame styles to suit a rider's individual needs and preferences. From the highly adjustable Invacare® Top End® Crossfire™ T7A Wheelchair to the Invacare® Top End® Crossfire™ T6 Wheelchair with a fixed rear seat-to-floor height, you can find the ideal chair for the most energy efficient ride. The Crossfire Series Wheelchairs are the perfect fit for anyone, from the new rider to the more experienced, who wants a custom, super light chair built for performance.

Description	Seat Width	Seat Depth	Warranty	Product Weight	Max. User Weight
Crossfire T6	300-500 mm	300-500 mm	24 months	8 kg	113 kg



SOLARA 3G TILT IN SPACE

- Mono-tube, open frame design making vehicle transfer easy
- Choice of fold down or fixed back
- 4.25" of center-of-gravity adjustment
- Full offerings of wheel and tire options

Description	Seat Width	Seat Depth	Back Height	Floor to Seat Height	Warranty	Product Weight	Max. User Weight
Solara 3G	300- 610 mm	300 - 560 mm	430 - 610mm	317-483mm	24 months	8 kg	136 kg



FREEDOM NXT

The NXT is the most rigid, folding, paediatric tilt-in-space manual wheelchair on the market. Its strong and durable, yet compact and lightweight. The NXT ensure both the child and carer remain comfortable by providing a smoother ride and a user friendly solution. When the growth spurt hits and the frame won't fit, the NXT adjusts easily for a good fit today and tomorrow.

Description	Seat Width	Seat Depth	Warranty	Max. User Weight
Freedom NXT	254 - 508 mm	305 - 406 mm	24 months	102 kg



E des@aidacare.com.au W www.aidacare.com.au







POWERED WHEELCHAIRS / SCOOTERS

Power wheelchairs are customisable and have more options for seating. For users who require positioning devices for trunk control and stability, a power wheelchair can be customized with seat and back cushions with many positioning features. Power wheelchairs also have options for electric tilt and recline, which are features that can be beneficial for people who have limited ability to reposition themselves in the wheelchair. Tilt and recline features provide comfort and pressure relief from sitting for long periods of time.

A main difference between a power wheelchair and electric scooter is the turning radius. Power wheelchairs typically have a smaller turning radius, which means that they are more maneuverable in small spaces. Power wheelchairs also have more options with center wheel drives or front wheel drives for a small turning radius in tight spaces.

Both electric wheelchairs and electric scooters allow a person to travel for long distances. This is an important benefit for people with limited mobility as it provides opportunity to go farther in their community than they might in a manual wheelchair.

Power wheelchairs also offer the benefit of traveling over a variety of terrain through the options of tire types and wheel positions. The most popular position for wheels in power wheelchairs are rear wheel drive, where the large wheel is in the back and the smaller casters are in front, which provides stability in driving both indoor and outdoor. Alternatively, center wheel drives or front wheel drives, where the large wheels are placed in the center or the front of the power wheelchair change the center of gravity for the wheelchair and can help in maneuvering over obstacles and curbs.



A200

The A200 is a compact rear wheel drive chair designed for all-round basic indoor and outdoor performance. A200 is easy to disassemble and fits into most cars for transportation. The standard seat frame is compatible with most modular seating options and is adjustable in width, depth and back angle. As an alternative to the standard seat, you may prefer to choose the optional Contour Seat.

- Adjustable seat frame in widths 38-48cm, depths 38-50cm
- Tension adjustable backrest upholstery and standard seat cushion
- Height adjustable armrests
- Swing away removable legrests
- Basic seating package

Description	Seat Width	Seat Depth	Range	Max Speed	Turning Radius	Product Weight	Max. User Weight
A200	380-480 mm	360-500 mm	22 km	6km/h	700 mm	66 kg	100 kg



DES Enquiries:

1300 133 120

des@aidacare.com.au W www.aidacare.com.au





TDX SP

The TDX SP (Total Driving eXperience) takes all the best features of the TDX series and improves upon them with enhanced performance, superior ride quality, quieter chassis and elegantly simple designs. With built-in technologies such as True Centre Wheel drive, Enhanced SureStep, Traction Control Design, Quiet Stability Lock, MK6i Electronics and powerful 4 pole motors, the TDX SP is designed to meet the needs of rehab clients who demand performance and style

- Traction control design
- Faster & more powerful
- Quiet stability lock

- Elegant simplicity
- Seating options

Description	Seat Width	Seat Depth	Turning Radius	Max Speed	Max. User Weight
TDX SP	300 - 600 mm	300 - 550 mm	560 mm	9.6 km/h	136 kg



Q6 EDGE

The latest powered in the Quantum range, the q6 Edge supersedes the Q600 offering superior drive and functionality similar to the Q6000 but at more economical price. Fully compatible with the Tru-balance seating system and equipped with ATX suspension the Q6 Edge offers comprehensive seating options with high levels of stability.

- Exceptional performance and value
- Full rehab capability
- Standard 4-pole motors
- Low seat-to-floor heights
- Adult model pictured

Description	Seat Width	Seat Depth	Floor to Seat Height	Max Speed	Warranty	Max. User Weight
Q6 Edge - Adult	250 -560 mm	250-560 mm	400-450 mm	10km/h	24 months	136 kg
Q6 Edge - Paediatric	250-560 mm	250-560 mm	440-490 mm	9.7km/h	24 months	136 kg



TORQUE SP

The Invacare® Storm Series® Torque™ SP power wheelchair with 4-pole motors and MK Series™ power wheelchair electronics offers many of the features and options found on the other Storm Series models while maintaining a truly affordable price. Built on the rugged Storm Series H-frame design, the Torque SP power wheelchair offers a proven, durable construction and rubber element independent rear suspension for a softer ride. Rehab has equipped the unit with powerful new GLIDE Motors for added power and torque.

- MKIV ™ RII non-expandable power wheelchair electronics standard
- TT ® TrueTrack technology (optional) keeps the power wheelchair on a true forward path and reduces driver fatigue
- 4-pole motors provide true gear

reduction and higher torque for many rehab needs

 Rubber element rear wheel suspension absorbs significant jolts and impacts. Upgrade to Monroe® independent rear suspension and/or shock forks for additional vibration dampening

Description	Seat Width	Seat Depth	Floor to Seat Height	Max Speed	Warranty	Max. User Weight
Torque SP	300 - 550 mm	300 - 550 mm	420-530 mm	9km/h	24 months	136 kg

P 1300 133 120 E des@aidacare.com.au

W www.aidacare.com.au





R44 LIGHTNING

The new R44 rear-wheel drive power chair features a high-forward pivot long trail arm, a state-of-the-art design that provides dramatically improved obstacle climbing performance and a smoother ride. The R44 also has outstanding tight space manoeuvrability, a sliding battery tray for easy access and accepts rehab seating and electronics. With its stylish, sleek design and a standard speed of up to 9.65 km/h (6 mph), the R44 is a superior power chair for even the most active clients.

- Narrow 64.1 cm (25") width and compact 88 cm length
- Outstanding tight space manoeuvrability
- Maximum speed up to 9.65 km/h (6 mph)
- High-performance, high-efficiency motor package
- Accepts rehab seating
- Swing-away joystick
- Accommodates full range of Pride® seating and accessories
- Sliding battery tray for easy access
- Available in paediatric sizes

Description	Seat Width	Seat Depth	Max Speed	Warranty	Max. User Weight
R44 Lightning	460 mm	410 mm	6.5 km/h	24 months	136 kg



TDX SP CG (TILT) & MCG TILT & RECLINE

The TDX SP is the first model in the next generation of TDX (Total Driving eXperience) power wheelchairs. The TDX SP takes all the best features of the TDX series and improves upon them with enhanced performance, superior ride quality, quieter chassis, and an elegantly simple design.

- The available Invacare FormulaCG Powered Seating which features up to 7" of weight shift tilt
- Rehab (ASBA) seat that adapts for growth in both width and depth
- Semi-recline captain[™] seat is available in 16", 18", 20" and 22" widths
- Optional heavy-duty package

Description	Seat Width	Seat Depth	Turning Radius	Max Speed	Max. User Weight
TDX SP	300 - 600 mm	300 - 550 mm	56 cm	9 km/h	136 kg



Q6000Z

The Quantum 6000 Rehab Power Wheelchair incorporates aggressive new features, powerful 4-pole Hammer motors, and patented Mid-Wheel 6TM design platform to deliver unmatched performance and tight-space maneuverability while delivering limitless rehab capability.

Description	Seat Width	Seat Depth	Floor to Seat Height	Max Speed	Warranty	Max. User Weight
Q6000Z	250-610 mm	250-610 mm	460-510 mm	18 km/h	24 months	136 kg



DES Enquiries:





QUICKIE XPLORE

Quickie Xplore is an excellent choice for an active user spending time outdoors. The directional stability at high speeds in all types of environments, helps ensure that the end user's optimal Center of Gravity will be maintained without any adjustments required.

Description	Seat Width	Seat Depth	Floor to Seat Height	Max Speed	Warranty	Max. User Weight
Quickie Xplore	250-610 mm	250-560 mm	400 mm	9.7 km/h	24 months	181 kg



Subscribe to our monthly solutions e-bulletin by visiting www.aidacare.com.au to recieve the latest equipment and industry news.

- New product releases and reviews
- Industry news & events
- Aidacare professional training (APT) sessions
- Aidacare news and views

You can also like our new facebook page!



E des@aidacare.com.au W www.aidacare.com.au





High Profile



Low Profile

ROHO SINGLE VALVE HIGH PROFILE/ LOW PROFILE CUSHION

- The Roho High Profile Wheelchair Cushion The LOW PROFILE Cushion is features 4" (10 cm) interconnected air cells to provide an adjustable, low deformation, counter-pressure environment that assist in the healing of ischemic ulcers. Appropriate for the high risk client.
- The LOW PROFILE ROHO cushion is an adjustable skin protection and positioning wheelchair cushion that has 2" (5cm) high cells to provide excellent protection against skin breakdown. Also available with dual valves for side-to-side or frontto-back positioning and support.
- recommended for clients who are at moderate risk to high risk -or skin breakdown, have no or partial sensation, are able to do independent weight shifts, and can propel their wheelchair with their feet.

Please consult with Aidacare's Complex Rehab Specialists for this item.

All models include

- Cushion Cover
- Pump and patch kit
- Operating Instructions
- 36 Month Warranty

AidacareCode	Supplier Code	Profile	Cushion Size(width x depth)	Cells	Max User Weight
PTC543700	1R1010C	HP	18.5 x 18.5″	10 x 10	Unlimited
PTC543800	1R1010LPC	LP	18.5 x 18.5″	10 x 10	Unlimited
PTC543900	1R1011C	HP	18.5 x 20.25"	10 x 11	Unlimited
PTC543910	1R1011LPC	LP	18.5 x 20.25"	10 x 11	Unlimited
PTC544000	1R109C	HP	18.5 x 16.75"	10 x 9	Unlimited
PTC546000	1R109LPC	LP	18.5 x 16.75"	10 x 9	Unlimited
PTC544200	1R1110C	HP	20.25 x 18.5"	11 x 10	Unlimited
PTC544250	1R1110LPC	LP	20.25 x 18.5"	11 x 10	Unlimited
PTC544300	1R1111C	HP	20.25 × 20.25"	11 x 11	Unlimited
PTC544350	1R1111LPC	LP	20.25 × 20.25"	11 x 11	Unlimited
PTC544400	1R119C	HP	20.25 x 16.75"	11 x 9	Unlimited
PTC544410	1R119LPC	LP	20.25 x 16.75"	11 x 9	Unlimited
PTC544500	1R1210C	HP	22 x 18.5"	12 x 10	Unlimited
PTC544550	1R1210LPC	LP	22 x 18.5″	12 x 10	Unlimited
PTC544560	1R1211C	HP	22 × 20.25"	12 x 11	Unlimited
PTC544590	1R1211LPC	LP	22 × 20.25″	12 x 11	Unlimited
PTC544600	1R1212LPC	LP	22 × 22″	12 x 12	Unlimited
PTC544420	1R129C	HP	22 x 16.75"	12 x 9	Unlimited
PTC544430	1R129LPC	LP	22 x 16.75"	12 x 9	Unlimited
PTC545900	1R910C	HP	16.75 x 18.5″	9 x 10	Unlimited
PTC546000	1R910LPC	LP	16.75 x 18.5″	9 x 10	Unlimited
PTC546100	1R911C	HP	16.75 × 20.25"	9 x 11	Unlimited
PTC546200	1R911LPC	LP	16.75 x 20.25"	9 x 11	Unlimited
PTC546600	1R99C	HP	16.75 × 16.75"	9 x 9	Unlimited
PTC546700	1R99LPC	LP	16.75 × 16.75"	9 x 9	Unlimited





High Profile Dual Valve



Low Profile Dual Valve With Cover



High Profile

Please consult with Aidacare's Complex Rehab Specialists for this item.

ROHO® - DUAL VALVE HIGH PROFILE/ LOW PROFILE

- Two separate compartments of the cushion allow positioning adjustments to address a pelvic obliquity (side-to-side) or sacral sitting (front to back).
- The Dual Compartment offers slightly less postural control and flexibility than the Quadtro Select range but retains the pressure and shear risk management standard for ROHOs
- The High Profile Cushion features 4" (10 cm) interconnected air cells for higher pressure disbursement
- The Low Profile Cushion features 2" (5cm) high cells for greater stability
- Appropriate for clients with high risk of ischemic ulcers

Aidacare Code	Supplier Code	Profile	Cushion Size (width x depth)	Cells	Max User Weight
PTC547300	2R1010C	HP	18.5 x 18.5″	10 x 10	Unlimited
PTC547400	2R1010LPC	LP	18.5 x 18.5″	10 x 10	Unlimited
PTC547290	2R1011C	HP	18.5 x 20.25"	10 x 11	Unlimited
PTC547295	2R1011LPC	LP	18.5 x 20.25"	10 x 11	Unlimited
PTC547500	2R109C	HP	18.5 x 16.75"	10 x 9	Unlimited
PTC547505	2R109LPC	LP	18.5 x 16.75"	10 x 9	Unlimited
PTC548000	2R910C	HP	16.75 x 18.5"	9 x 10	Unlimited
PTC548100	2R910LPC	LP	16.75 x 18.5"	9 x 10	Unlimited
PTC548400	2R99C	HP	16.25 x 16.75"	9 x 9	Unlimited
PTC548500	2R99LPC	LP	16.25 x 16.25"	9 x 9	Unlimited

ROHO® QUADTRO SELECT CUSHION

- ROHO Dry Floatation promotes circulation, and pressure relief whilst flexible neoprene cells reduce shear friction
- Quadtro Select cushions employ ISOFLO Memory Control which locks air into compartments to help maintain better postural integrity in addition to pressure reduction
- Quick-and-easy air compartment adjustment maximizes function and utility

LOW PROFILE (LP)

• Greater stability

HIGH PROFILE (HP)

- 4"(10cm) High Cells 2" (5cm) High Cells
- Higher pressure disbursement





ISOFLO - Open

ISOFLO - Closed

All models include

- 2 way stretch cover
- Hand Pump & Repair Kit
- Operating Instructions
- 36 Month Warranty

	<u>0</u>			DI	EPTH (Cells / Size	e")		
WIDTH (Cells / Size")	Profile	12 / 22″w	11 / 20″	10 / 18.5″	9 / 16.75″	8 / 15″	7 / 13″	6 / 11″
12 / 22″	ΗP	PTC589900		PTC589600	PTC589500			
12,22	LP			PTC589700				
11 / 20″	HP	PTC589200	PTC589000	PTC588800	PTC589400	PTC589300		
11 / 20	LP		PTC589100	PTC588900	PTC589500			
10 / 18.5″	ΗP	PTC588200	PTC588000	PTC587700	PTC588500	PTC588300		
10 / 10.5	LP		PTC588100	PTC587900	PTC588700	PTC588400		
9 / 16.75″	HP	PTC592600	PTC592400	PTC592200	PTC593000	PTC592800	PTC592700	
7 10.75	LP		PTC592500	PTC592300	PTC593200	PTC592900		
8 / 15″	HP		PTC591500	PTC591300	PTC592000	PTC591800	PTC591600	
0713	LP			PTC591400	PTC592100	PTC591900	PTC591700	
7 / 13"	HP				PTC591200	PTC591000	PTC590700	
,,15	LP					PTC591100	PTC590800	
6 / 11"	HP						PTC590400	PTC590200
0711	LP					PTC590600	PTC590500	PTC590300

(Available in a Mid-profile version contact Aidacare C.A.S.E)







Please consult with Aidacare's Complex Rehab Specialists for this item.

ROHO® - HYBRID ELITE CUSHION

- Two trusted brands, one effective solution. The Hybrid Elite combines the skin protection and shape matching ability of ROHO Dry Floatation with a contoured foam base designed for ROHO by JAY
- The design allows for improved pelvic alignment and stability while maintaining optimal immersion and dispersion of pressure away from the pelvic area

All models include

- 2-way stretch cover
- Hand Pump & Repair Kit

- The stability provided by the JAY foam base also assists with achieving safer transfers
- The ROHO Dry Floatation system and 2-way stretch cover help reduce friction shear whilst maximising blood flow
- Appropriate for clients with moderate to high risk of skin breakdown who require positioning support



Hyrbid Elite Cushion Assembly

- Operating Instructions
- 24 Month Warranty

Aidacare Code	Supplier Code	Cushion Size (width x depth)
PTC594845	HE-1RHY1414C	14 x 14
PTC594850	HE-1RHY1415C	14 x 15
PTC594860	HE-1RHY1416C	14 x 16
PTC594899	HE-1RHY1515C	15 x 15
PTC594900	HE-1RHY1516C	15 x 16
PTC594901	HE-1RHY1517C	15 x 17
PTC594909	HE-1RHY1615C	16 x 15
PTC594910	HE-1RHY1616C	16 × 16
PTC594911	HE-1RHY1617C	16 x 17
PTC594920	HE-1RHY1618C	16 x 18
PTC594921	HE-1RHY1620C	16 × 20
PTC594929	HE-1RHY1716C	17x 16
PTC594930	HE-1RHY1717C	17 x 17
PTC594931	HE-1RHY1718C	17 x 18
PTC594935	HE-1RHY1816C	18 x 16
PTC594939	HE-1RHY1817C	18 x 17
PTC594940	HE-1RHY1818C	18 x 18
PTC594945	HE-1RHY1820C	18 × 20
PTC594950	HE-1RHY2018C	20 × 18
PTC594951	HE-1RHY2020C	20 × 20
PTC594946	HE-1RHY2218C	22 x 18
PTC594947	HE-1RHY2220C	22 × 20
PTC594948	HE-1RHY2418C	24 x 18
PTC594949	HE-1RHY2420C	24 × 20







Closed-cell ultra-lightweight base
'Jay Flow' fluid Tripad

3 JAY AirExchange or Incontinence cover (optional)

JAY - J2 CUSHION

- Designed for the client who has moderate to more aggressive postural needs that may change over time
- The J2 pre-contoured, closed-cell, lightweight, foam base provides excellent lateral and forward/rearward stability creating a fountation for postural stability
- Designed for the client with moderate to aggressive postural needs that may change over time
- The J2 pre-contoured, closed-cell, lightweight, foam base provides excellent lateral and forward/rearward stability creating a fountation for postural stability
- The base fluid pad can be modified/customised on site to accommodate fixed postures or to correct flexible postures. Modifications can support changing needs over time

- Optional solid seat inserts provide extra postural stability discouraging the cushion from slinging with the chair upholstery
- Postural accessories such as medial and lateral thigh supports, lateral pelvic supports and pelvic obliquity supports can be added to the cushion
- The 'Jay Flow' fluid Tripad conforms to bony prominences, distributes pressure more evenly and reduces peak pressures. The Tripad design assists in maintaining fluid under the at-risk bony areas
- Most suitable for the client who is at higher risk of skin breakdown and who has poor skin integrity
- Optional Balistic Stretch; AirExchange or Incontinence cover.

Please consult with Aidacare's

Complex Rehab Specialists for

this item about the J2 Recline option for reclining wheelchairs.

• 24mths Warranty (6mths on cover)

Cushion Comparison

Lateral Stability

Forward Stability

Pressure Management

Aidacare **Chair Size** Max User Supplier Height Code Code (width x depth) (inches) Weight 16 x 16" PTC562500 J2105 3″ 114kg PTC562700 J2106 16 x 18" 3″ 114kg 3″ PTC563800 J2117 17 x 17" 114kg J2100 3″ PTC561100 18 x 16" 114kg PTC561300 J2101 18 x 18" 3″ 114kg 18 x 20" 3″ PTC562100 J2103 114kg PTC563700 J2116 20 x 16" 3″ 114kg 3″ PTC564100 J2118 20 x 18" 114kg



JAY - J2 DEEP CONTOUR (DC) CUSHION

- The J2-DC takes the J2 to a higher level of pressure and postural management
- \bullet Cushion depth is increased in the DC from 3" to 4"
- Pelvic contour depth increases from 2.5" to 3"
- An optional solid drop-seat (shown bottom) can be permanently attached to promote a level platform

Aidacare Code	Supplier Code	Chair Size (width x depth)	Height (inches)	Max User Weight
PTC566000	J2466	16 x 16″	4″	114 kg
PTC566100	J2468	16 x 18"	4″	114 kg
PTC566400	J2477	17 x 17"	4″	114 kg
PTC566600	J2486	18 x 16″	4″	114 kg
PTC566700	J2488	18 x 18″	4″	114 kg
PTC566500	J2480	18 x 20"	4″	114 kg
PTC565150	J2406	20 x 16"	4″	114 kg
PTC565000	J2408	20 x 18"	4″	114 kg

- '3D Jay Flow' fluid Tripad provides 3" deepimmersion pressure and postural management
- Most suitable for the client who is at the highest level of risk for skin breakdown or who has very poor skin integrity
- Optional Balistic Stretch; AirExchange or Incontinence cover. 24mths Warranty on Cushion (6mths on cover)

	Solid seat insert (optional)
2	Closed-cell ultra-lightweight base
3	Leg Trough modification for hip ROM.
4	3D Jay Flow' fluid Tripad

5 JAY AirExchange or Incontinence covers

Cushion Comparison

Lateral Stability	1	2	3	4	5
Forward Stability	1	2	3	4	5
Pressure Management	1	2	3	4	5

5

5

4 5

MOBILITY SCOOTERS



CELEBRITY XL

- High performance mobility scooter
- Large knobby tyres for excellent outdoor traction
- Deluxe front-mounted high back swivel seat
- Lumber support and manual recline for increased comfort
- State of the art rear suspension
- Unique easily adjusted steering mount with modern LED lighting
- Rear view mirror and full directional lighting package including front and rear turn signal

Code	Description	Colour	Length	Width	Max.Speed	Warranty	Product Weight	Max. UserWeight
SCT670540	Celebrity XL	Red	1300 mm	640 mm	9.5 km/h	24 Months (12 months batteries)	91 kg (without batteries)	159 kg



CELEBRITY DX

- High performance mobility scooter
- Large knobby tyres for excellent outdoor traction
- Deluxe front-mounted high back swivel seat
- Lumber support and manual recline for increased comfort
- State of the art rear suspension
- Unique easily adjusted steering mount with modern LED lighting

Code	Description	Colour	Length	Width	Max. Speed	Warranty	Product Weight	Max. User Weight
SCT670530	Celebrity DX	Red	1190 mm	610 mm	Up to 9.5 km/h	24 months	77 kg (without batteries)	159 kg



DX SPORT

- State-of-the-art rear suspension
- Larger rear wheels
- Top speed of up to 12km/h
- Maximum weight capacity of 159k
- Ultra simple, auto-latching lockup for quick disassembly
- Easily accessible tie-down points
- Front basket and rear vision mirror is standard

- Tight turning circle radius of 145cm
- Includes standard angle-adjustable dual bulb kerb light
- Front frame-mounted seat post for optimal stability
- High visibility reflectors
- Tiller rotation lock permits tiller to be used as a handle during transport

Code	Description	Colour	Seat Width	Seat Depth	Length	Max Speed	Warranty	Product Weight	Max User Weight
SCT670535	DX Sport	Yellow	610 mm	250 - 510 mm	1190 mm	12 km/h	24 months	68 kg (without batteries)	159 kg



DES Enquiries:

1300 133 120 E des@aidacare.com.au

W www.aidacare.com.au





PATHRIDER 130XL

- Top of the range sleek and sporty powered mobility scooters
- Large 13" pneumatic tyres to handle rugged outdoor terrain
- Luxurious suspension seat with arm rest, head rest and lumbar support for stability

Code	Description	Length	Width	Max. Speed	Warranty	Product Weight	Max. User Weight
SCT670572	Pathrider 130XL	1352 mm	640 mm	10 knm/h	24 months (12 months batteries)	101 kg (without batteries)	181 kg



PATHRIDER 140XL

- Top of the range sleek and sporty powered mobility scooters
- Large 13" pneumatic tyres to handle rugged outdoor terrain
- Luxurious suspension seat with arm rest, head rest and lumbar support for stability

Code	Description	Colour	Length	Width	Max. Speed	Warranty	Product Weight	Max. User Weight
SCT670578	Pathrider 140XL	Yellow	1428 mm	692 mm	10 km/h	24 months	117.70 kg (without batteries)	181 kg
SCT670576	Pathrider 140XL	Red	1428 mm	692 mm	10 km/h	24 months	117.70 kg (without batteries)	181 kg
SCT670577	Pathrider 140XL	Silver	1428 mm	692 mm	10 km/h	24 months	117.70 kg (without batteries)	181 kg



AIDACARE SOLUTIONS / June 2014 aidacare solut Complex Rehab

Subscribe to our monthly solutions e-bulletin by visiting www.aidacare.com.au to recieve the latest equipment and industry news.

- New product releases and reviews
- Industry news & events
- Aidacare professional training (APT) sessions
- Aidacare news and views

You can also like our new facebook page!





Are you keen to learn about the latest products and technology whilst accruing CPD points at the same time?

- Learn from our highly skilled APT educators/presenters
- Participate in hands-on sessions
- Take home your Certificate of Attendance crucial for your CPD activities and professional registration
- Get to know our Specialists who can assist with advice and/or problem solving and network with fellow therapists.

All our sessions are free of charge! Don't miss out, subscribe to our mailing list now to ensure you receive emails on the upcoming events in your area!



About this Catalogue

Due care has been taken to ensure all the information contained in this catalogue is correct at the time of printing. All measurements, pictures, colours and weight capacities are to be used as a guide only. Our suppliers reserve the right to modify the design or appearance of any product displayed in this catalogue without prior notice. © Copyright Aidacare Pty Ltd 2014 Version 1.0

No part of this catalogue may be copied, distributed, reproduced, scanned or stored on any electronic database, whether in whole or in part in any form or by any means.