

Pressure Care & Circulation



Contents

P 436	Protection - General
P 440	Protection - Hip
P 445	Protection - Skin + Limb
P 449	Compression - Legwear
	P 451 - Dr. Comfort P 451 - GloriaMed P 455 - Ibici P 457 - Jobst P 459 - Mediven P 465 - Sigvaris P 467 - Solaris P 472 - TED Anti-Embolism P 473 - Venosan P 480 - Stocking Accessories
P 482	Compression - Armwear

- P 483 Gloves (Norco/Isotoner/Handeze) P 484 - Jobst P 487 - Mediven P 490 - Sigvaris P 491 - Solaris
- P 493 Compression Lymphoedema





Protection - General









Sheepskin Natural Gold **Overlay**

- Genuine Australian lambs wool overlay provides additional comfort to an existing mattress
- Soft pile allows air circulation to keep skin dry



Sheepskin Hi Temp Classic **Mattress Overlay**

- A full covered medical sheepskin bed overlay, designed to cover head to toe
- Suitable for long term patients where shear and friction are a concern

Code	Description	Length	Width	DVA:	Code	Description	Length	Width	DV
PTA502500	Sheepskin Overlay	1000 mm	750 mm	AB11	PTA521400	Sheepskin Mattress Overlay	1250 mm	750 mm	AB



Sheepskin Hi Temp **Medical Overlav**

- Medical grade sheepskin is a versatile product that can be used in the bed, on a chair or on a car seat
- Assists in pressure reduction, will absorb water, transfer heat away from the body and reduce shear and friction

Code	Description	Length	Width	DVA:
PTA521500	Medical Sheepskin Overlay	1000 mm	650 mm	AB11



Sheepskin Leg Rest Cover

- High quality medical grade sheepskin covers provide natural comfort for any wheelchair
- Provide excellent skin protection and comfort

Code	Description	DVA:
PTA522100	Wheelchair Cover	AB11
PTA522101	Wheelchair Calf Protector Cover	AB11



Code

Gel & Foam Underarm **Crutch Covers**

• Fit easily over the existing foam of the underarm/axilla crutches providing a comfortable support surface for the user



Sheepskin Cushion/Back Cover

- The Cushion It Sheepskin Cover can be used directly on a seat or back canvas of a wheelchair
- Can be placed over standard foam cushions to increase pressure care and absorb perspiration

Code	Description	Length	Width	DVA:
PTA521300	Sheepskin Hi Temp - Cushion/Back Cover	500 mm	500 mm	AB11

Description

PTA509020 Gel and Foam Underarm Crutch Axilla Covers PTA509021 Gel Underarm Crutch Handle Covers

Gel Armrest Pads (Pair)

Gel Armrest Pads (Pair)

PTA508450 Gel Armrest Pads (Pair)





PTA508400

PTA508425

Gel Armrest Pads

- Ultimate gel comfort and protection for the elbows and wrists
- Distributes pressure efficiently to help avoid harmful pressure-related problems
- Durable neoprene nylon surface for easv cleaning

88.9 mm 304.8 mm

88.9 mm 355.6 mm

88.9 mm 406.4 mm

Lenath

DVA

AP14

AP14

AP14



Gel Armrest Trough Liner

- Fits onto the Ottobock trough channel forearm pad armrest
- Provides a comfortable, supportive gel support pad for the user

Code	Description	Length
PTA508470	Gel Armrest Trough Liner - Medium	330.2 mm
PTA508480	Gel Armrest Trough Liner - Large	482.6 mm



Gel Calf Support Panel

- Provide comfortable support to the lower legs in the seating position
- Single piece model features large gel pad, with 6 velcro secure straps
- Split Panel model allows use of swing away legrests without removal of the calf panel

Code	Description	Fits Wheelchair	DVA:
PTA508500	Gel Calf Support Panel (Single)	40.64 - 50.8 cm	AP14
PTA508550	Gel Calf Support Panel (Pair)	40.64 - 50.8 cm	-



Gel & Foam Footplate **Protector Covers**

- Designed to fit over existing footplates on a wheelchair
- Highly efficient at distributing pressure thereby avoiding harmful pressurerelated problems

Code	Description	Size	DVA:
PTA508600	Gel Covers (Pair)	203.2 x 152.4 mm	AP14
PTA508610	Gel Covers w/ Adhesive Back (Pair)	203.2 x 152.4 mm	-
PTA508615	Gel Cover (Single)	279.4 x 228.6 mm	-



Heelift Suspension Boot

- Designed by orthopaedic surgeons specifically for prevention and treatment of heel pressure ulcers
- Effectively eliminates heel pressure by suspending the heel in the air
- Two options available: Smooth or Convoluted

Code	Description	DVA:
PTA509200	Heelift Suspension Boot Smooth	AB11
PTA509100	Heelift Suspension Boot Convoluted	AB11



DVA Contact Details

AB11



Pressure Care Assist

- Can be used anywhere in the bed to assist with pressure care
- Made from medical grade sheepskin
- Helps reduce the risk of pressure area breakouts in the heels, sacrum, shoulders or head

Code	Description	DVA:
PTA521600	Pressure Care Assist	AB11



Gel Knee Protector Buttons

- Provide a soft surface to protect the users knees from bumping against learests
- Attach to the inside of the legrests on a wheelchair
- Includes mounting hardware

Code	Description	DVA:
PTA508800	Gel Knee Protector Buttons - (3.18) Swing Away Legrests	AP14
PTA508820	Gel Knee Protector Buttons - (3.18) Elevating Legrests	AP14
PTA508810	Gel Knee Protector Buttons - (1.98) Swing Away Legrest	AP14
PTA508830	Gel Knee Protector Buttons - (1.98) Elevating Legrest	AP14



I

ndis

Sheepskin Short Slipper Boot

- Assists in the prevention of pressure ulcers on the feet with an existing injury or condition
- Designed for use primarily in bed, wheelchair use or recliner chair use

Code	Description	Shoe Size US	DVA:
PTA521897	Short Slipper Boot - XX-Small	3 - 4	AB11
PTA521898	Short Slipper Boot - X-Small	5 - 6	AB11
PTA521899	Short Slipper Boot - Small	7 - 8	AB11
PTA521900	Short Slipper Boot - Medium	9 - 10	AB11
PTA521901	Short Slipper Boot - Large	11 - 12	AB11
PTA521902	Short Slipper Boot - X-Large	13 - 14	AB11
PTA521903	Short Slipper Boot - XX-Large	14 - 15	AB11













Sheepskin Wrap Around Boot

- Wrap Around Sheepskin Boots designed for those at risk of foot ulcers and other types of tissue damage
- Constructed with low profile EVA sole, these boots can be worn outside



Sheepskin Pressure Care Boot

- Sheepskin Boots designed for those at risk of foot ulcers and other types of tissue damage
- · Constructed with low profile EVA sole, these boots can be worn outside

Shoe Size US

3 - 4

5 - 6

7 - 8

9 - 10

11 - 12

13 - 14

14 - 15

AB11

AB11

AB11

AB11

AB11

AB11

AB11

Code	Description	Shoe Size US	DVA:	Code	Description
PTA521704	Wrap Around Boot - XX-Small	3 - 4	AB11	PTA521697	Pressure Care Boot - XX-Small
PTA521705	Wrap Around Boot - X-Small	5 - 6	AB11	PTA521698	Pressure Care Boot - X-Small
PTA521706	Wrap Around Boot - Small	7 - 8	AB11	PTA521699	Pressure Care Boot - Small
PTA521707	Wrap Around Boot - Medium	9 - 10	AB11	PTA521700	Pressure Care Boot - Medium
PTA521708	Wrap Around Boot - Large	11 - 12	AB11	PTA521701	Pressure Care Boot - Large
PTA521709	Wrap Around Boot - X-Large	13 - 14	AB11	PTA521702	Pressure Care Boot - X-Large
PTA521710	Wrap Around Boot - XX-Large	14 - 15	AB11	PTA521703	Pressure Care Boot - XX-Large



- This product has been designed to protect the entire lower limb
- The Tall Slipper Sheepskin Boot can be worn in bed, in a chair or wheelchair and allows the use of compression bandages or support hosiery

Æ		
	1	
	0	

Sheepskin Sovereign Snug

- Beautifully crafted Sheepskin Snugs will accommodate mild foot oedema whilst offering support and comfort
- Made especially from carefully selected lambskin material

Code	Description	DVA:	
PTA520500	Tall Slipper Boot - One Size Fits All	AB11	

Code	Description	Shoe Size US	DVA:
PTA521905	Sovereign Snug - Small	7 - 8	AB11
PTA521906	Sovereign Snug - Medium	9 - 10	AB11
PTA521907	Sovereign Snug - Large	11 - 12	AB11



Heel Protectors Silicone® Fibre

- Provides comfort and assists in the prevention of pressure build up in the heel area
- Silicone hollow core fibres provide soft comfortable support, while still keeping the rest of the foot cool in warm weather
- Sold as a pair

Code	Description	Length	Width	Height	DVA:
PTA503350	Heel Protectors Silicone® Fibre	140 mm	100 mm	30 mm	AB11



General Enquiries

ROHO Heelpad

- ROHO Heelpad offers an adjustable therapeutic environment which aids in the healing of heel and elbow ischemic ulcers
- Made from soft, latex free neoprene the heelpad is extremely durable

Code	Description	DVA:
PTC541600	ROHO Heelpad	AB17



Bed Heel Elevator

- Wedge-shaped contoured foam provides excellent support for the lower lea
- Rigid foam aids blood flow and reduces pressure

Code	Description	DVA:
PTA509050	Bed Heel Elevator	AB17



Sheepskin Heel Supports

- Protect the heel from ulcerations caused by bumps and knocks
- Made from medical grade sheepskin
- Wraps around the heel and ankle leaving the rest of the foot unencumbered

Code	Description	Show Size US	DVA:
PTA521199	Heel Supports - Small (Pair)	7 - 8	AB11
PTA521200	Heel Supports - Medium (Pair)	9 - 10	AB11
PTA521201	Heel Supports - Large (Pair)	11 - 12	AB11







Code

PTA530200

Description

Spenco

Heel Protectors

Heel Protectors Spenco

- Made from medically treated, washable natural cotton / polyester
- One size fits all, with a velcro band closer

Width

100 mm

Height

30 mm

AB11



Heel/Elbow Protector (With Gel Pad)

- Gel Elbow/Heel Protectors are made from a comfortable, stretchable acrylic/spandex
- Designed to prevent trauma or wearing of the heel or elbow
- Sold as a pair

Code	Description	Size	DVA:
PTA506400	Gel Heel/Elbow Protector - Small	230 mm	AR29
PTA506500	Gel Heel/Elbow Protector - Large	280 mm	AR29
PTA506600	Gel Heel/Elbow Protector - X-Large	330 mm	AR29



Elbow Lift Suspension Pad

- Cushions and elevates the elbow to min
- imise surface contact
- Helps protect the olecranon and olecranon bursa from friction burns and contusions
- Latex free

Length

140 mm

and the	
60	2
107	

Sheepskin Elbow Protector

- Provides elbow pressure point control
- Assists in avoiding skin tears, bruising and other tissue damage
- Easy and convenient velcro fastening
- Sold as a pair

Code	Description	Code	Description	DVA:
PTA508200	Elbow Lift Suspension Pad	PTA521100	Elbow Protector (One Size Fits All)	AB11



Sheepskin Palm Protector

- Provides comfort for patients with hand contractures
- Designed to offer moisture absorption, reducing the likelihood of maceration

• One size fits all

Code	Description
PTA521104	Palm Protector (Ambidextrous)

More from the Recovery Range



Aspire Vogue Carbon Fibre Seat Walker WAF705350 (Refer to Page 247)



DVA Contact Details

Foam Mattresses (Refer to Page 71)



Aspire Shower Commode with Platform Footrest BTC066000 (Refer to Page 20)



Pillows & Supports (Refer to Page 89)



Ottobock Cushions (Refer to Page 149)

DVA

AB11



Aspire Evoke 2 Wheelchairs - Range (Refer to Page 266)



Brochures also available.

Please contact us for a copy

Aspire Bedside Commode BEC045005 (Refer to Page 78)



Aspire Signature Lift Recline Chair - Space Saver CHP197700 (Refer to Page 111)



pressure

CORE 1300 133 120

Inida core

Aspire High Back Day Chair CHP208050 (Refer to Page 106)



Pressure Care & Circulation



(Refer to Page 86)

439

I For NDIS Item Code, refer to page: 500



HipSaver

Overview



- One of the most serious injuries, due to a fall, is a fractured hip especially amongst those 65 and older
- While medications, often increase the risk of falls and fall-related injuries, poor vision and difficulty in moving put the elderly at risk of a fall that results in a fatal or debilitating hip injury or hip fracture
- Seniors often develop a fear of falling, causing them to limit activities and leading to reduced mobility. HipSaver products restore confidence of motion to the elderly, both preserving physical health and increasing senior mobility
- Designed with the user in mind, HipSaver products are incredibly soft and comfortable to wear, whilst offering the most effective and superior protection
- HipSaver is available in a variety of designs to suit all types of wearers needs, ranging from the elderly in nursing home settings to those for seniors who daily participate in exercise programs and even robust activities such as jogging
- Whether you need hip protection for your personal use, or for your loved ones, or you are seeking hip protection for your Group Home, Nursing Home, or other Institution, HipSaver has a product to suit your needs

To Order

Use the Size Guide to select your size, then choose your preferred model and the appropriate code

Size Guide

	Size	Hips (cm)	Size	Hips (cm)
Sizing of all HipSaver products is from this guide	XS	70 - 79	L	100 - 112
Measure around widest part of hip (over incontinence products if applicable	S	80 - 89	XL	113 - 127
	Μ	90 - 99	XXL	128 - 146

Protection - Hip Pressure Care & Circulation



DVA: AR28

DVA: AR28

HipSaver - Nursing Home (Sewn-In Pads)

Designed for the mobility and comfort of residents plus convenience and peace of mind for carers

- Full, soft fabric crotch perfectly accommodates incontinence pads or diapers
- Elastic waistband and large leg openings for easy dressing and toileting
- Wearable 24 hours a day comfortable to sleep in •
- A popular all-round choice for most wearers
- Choose from High Compliance or High Compliance With Tailbone Pad for • those at risk of falling backward
- The Veterans Kit consists of 6 garments



Size	Nursing Home - High Compliance		Nursing Home - High Compliance With Tailbone Pad		Nursing Home - Veterans Kit (6 Garments)	
		Code - Female	Code - Male	Code - Female	Code - Male	Code - Female
XS	PTA510810	PTA510300	PTA492155	PTA492149	PTA513906	PTA513900
S	PTA510900	PTA509900	PTA492156	PTA492150	PTA513907	PTA513901
Μ	PTA510700	PTA509700	PTA492157	PTA492151	PTA513908	PTA513902
L	PTA510500	PTA509500	PTA492158	PTA492152	PTA513909	PTA513903
XL	PTA511100	PTA510100	PTA492159	PTA492153	PTA513910	PTA513904
XXL	PTA511225	PTA510220	PTA492160	PTA492154	PTA513911	PTA513905

Nursing Home - Pads & Pants - Separate

Size	Nursing Hom	Nursing Home - Pants Only		Nursing Home - Starter Kit (3 Pants & 1 Pair Pads)		Nursing Home - Veterans Kits (6 Pants & 2 Pairs Pads)	
		Code - Female	Code - Male	Code - Female		Code - Female	
XS	PTA512094	PTA512045	PTA511010	PTA510400	PTA512406	PTA512412	
S	PTA512095	PTA512050	PTA511000	PTA510000	PTA512407	PTA512413	
Μ	PTA512090	PTA512040	PTA510800	PTA509800	PTA512408	PTA512414	
L	PTA512080	PTA512030	PTA510600	PTA509600	PTA512409	PTA512415	
XL	PTA512070	PTA512020	PTA511200	PTA510200	PTA512410	PTA512416	
XXL	PTA512060	PTA512010	PTA511210	PTA510250	PTA512411	PTA512417	

HipSaver - Slim Fit (Sewn-In Pads)

Ideal for fashion conscious and active elders - those who are at risk of hip fracture, but reject uncomfortable and obvious hip protection products Looks good, feels good and pulls on just like regular underwear

- Soft, thin (13 mm) pads, tapered for protection with style
- Made from ultra-soft breathable polycotton with form fitting Spandex®
- SlimFit can accommodate an incontinence pad
- Fully machine washable and dryable, even at high temperature
- Choose from High Compliance or High Compliance with Tailbone Pad . forthose at risk of falling backward
- The Veterans Kit consists of 6 garments



1

Size	SlimFit - Higl	SlimFit - High Compliance		SlimFit - High Compliance With Tailbone Pad		SlimFit - Veterans High Compliance Kit (6 Garments)	
		Code - Female	Code - Male	Code - Female		Code - Female	
XS	PTA493705	PTA493701	PTA513803	PTA513801	PTA513880	PTA513820	
S	PTA493710	PTA493700	PTA513804	PTA513798	PTA513881	PTA513821	
Μ	PTA493720	PTA493600	PTA513805	PTA513799	PTA513882	PTA513822	
L	PTA493730	PTA493500	PTA513806	PTA513800	PTA513883	PTA513823	
XL	PTA493740	PTA493450	PTA513807	PTA513797	PTA513884	PTA513824	
XXL	PTA493745	PTA493400	PTA513808	PTA513802	PTA513885	PTA513825	





DVA: AR28

441



Slim Fit - Pads & Pants - Separate

Unit i						JVA. ANZO
Size	Slim Fit - Pants Only		Slim Fit - Pants Only (3 Pants & 1 Pair Pads)		Slim Fit - Veterans Kits (6 Pants & 2 Pairs Pads)	
	Code - Male	Code - Female		Code - Female	Code - Male	Code - Female
XS	PTA494090	PTA494053	PTA513650	PTA513210	PTA512418	PTA512424
S	PTA494080	PTA494055	PTA513600	PTA513200	PTA512419	PTA512425
Μ	PTA494075	PTA494050	PTA513500	PTA513100	PTA512420	PTA512426
L	PTA494070	PTA494045	PTA513400	PTA513250	PTA512421	PTA512427
XL	PTA494065	PTA494040	PTA513700	PTA513300	PTA512422	PTA512428
XXL	PTA494060	PTA494035	PTA513750	PTA513350	PTA512423	PTA512429

HipSaver - Quick Change (Sewn-In Pads)

Provides the same features and benefits of the Nursing Home model, but with the convenience of a snap down crotch panel.

- Incontinence pads and diapers can be changed easily and quickly either while the wearer is standing or lying down
- Fabric covered low profile snaps for ease and comfort
- Not suitable for wearers toileting themselves
- Laundry safe for washing and drying up to 95°C
- The Veterans Kit consists of 6 garments





Size	Quick Change - High Compliance		ge - High Compliance Quick Change - High Compliance With Tailbone Pad		Quick Change - Veterans High Compliance Kit (6 Garments)	
		Code - Female		Code - Female		Code - Female
XS	PTA512975	PTA512968	PTA513000	PTA512994	PTA512988	PTA512982
S	PTA512976	PTA512969	PTA513001	PTA512995	PTA512989	PTA512983
Μ	PTA512978	PTA512970	PTA513002	PTA512996	PTA512990	PTA512984
L	PTA512979	PTA512971	PTA513003	PTA512997	PTA512991	PTA512985
XL	PTA512980	PTA512972	PTA513004	PTA512998	PTA512992	PTA512986
XXL	PTA512981	PTA512973	PTA513005	PTA512999	PTA512993	PTA512987

Quick Change - Pads & Pants - Separate

Quick Change	ge - Pants Only Quick Change - Starter Kit (3 Pants & 1 Pair Pads)		Quick Change - Veterans Kits (6 Pants & 2 Pairs Pads)		
Code - Male	Code - Female	Code - Male	Code - Female	Code - Male	Code - Female
PTA494054	PTA494151	PTA513912	PTA512840	PTA512430	PTA512436
PTA494056	PTA494152	PTA513913	PTA512900	PTA512431	PTA512437
PTA494057	PTA494153	PTA513914	PTA512800	PTA512432	PTA512438
PTA494058	PTA494154	PTA513915	PTA512810	PTA512433	PTA512439
PTA494059	PTA494155	PTA513916	PTA512820	PTA512434	PTA512440
PTA494061	PTA494156	PTA513917	PTA512830	PTA512435	PTA512441
	Code - Male PTA494054 PTA494056 PTA494057 PTA494058 PTA494059	PTA494054 PTA494151 PTA494056 PTA494152 PTA494057 PTA494153 PTA494058 PTA494154 PTA494059 PTA494155	Code - Male Code - Female Code - Male PTA494054 PTA494151 PTA513912 PTA494056 PTA494152 PTA513913 PTA494057 PTA494153 PTA513914 PTA494058 PTA494154 PTA513915 PTA494059 PTA494155 PTA513916	Code - Male Code - Female Code - Male Code - Female PTA494054 PTA494151 PTA513912 PTA512840 PTA494056 PTA494152 PTA513913 PTA512900 PTA494057 PTA494153 PTA513914 PTA512800 PTA494058 PTA494154 PTA513915 PTA512810 PTA494059 PTA494155 PTA513916 PTA512820	Code Change - Pants Only (3 Pants & 1 Pair Pads) (6 Pants & 2 Code - Male Code - Female Code - Male Code - Male Code - Male PTA494054 PTA494151 PTA513912 PTA512840 PTA512430 PTA512430 PTA494056 PTA494152 PTA513913 PTA512900 PTA512431 PTA494057 PTA494153 PTA513914 PTA512800 PTA512432 PTA494058 PTA494154 PTA513915 PTA512810 PTA512433 PTA494059 PTA494155 PTA513916 PTA512820 PTA512434

HipSaver - Open Bottom (Sewn-In Pads)

Requested by nursing staff for residents needing a garment that permits easy self-toileting.

- The open crotch design is perfect for patients with urge incontinence or challenges in getting clothing on or off
- This garment is worn against the skin underneath any incontinence products, and allows a resident's normal underwear to be worn over the top of the protector
- Reduces accidents due to the 'fumble factor' ideal for Parkinson's sufferers

Open Bottom - High Compliance

PTA512398

PTA512399

PTA512401

PTA512402

PTA512403

PTA512404

- Laundry safe for washing and drying up to 95°C
- Suitable for both male and female users
- The Veterans Kit consists of 6 garments



XXL

442

Pressure Care & Circula

General Enquiries

DVA: AR28 ance Kit - Unisex (6 Garments

	-	
- Veterans High	Compliance	Kit - Unis
Р	TA492701	

PTA492702

PTA492703

PTA492704

PTA492705

PTA492706

Open Bottor

Protection - Hip Pressure Care & Circulation



DVA: AR28

Open Bottom - Pads & Pants - Separate

Open Bottom - Pants Only - Unisex	Open Bottom - Starter Kit - Unisex (3 Pants & 1 Pair Pads)	Open Bottom - Veterans Kits - Unisex (6 Pants & 2 Pairs Pads)
PTA512550	PTA512405	PTA512442
PTA512540	PTA512400	PTA512443
PTA512530	PTA512100	PTA512444
PTA512520	PTA512150	PTA512445
PTA512510	PTA512200	PTA512446
PTA512500	PTA512300	PTA512447
	PTA512550 PTA512540 PTA512530 PTA512520 PTA512510	Open Bottom - Pants Only - Onisex (3 Pants & 1 Pair Pads) PTA512550 PTA512405 PTA512540 PTA512400 PTA512530 PTA512100 PTA512520 PTA512150 PTA512510 PTA512200

HipSaver - Interim Overpant (Sewn-In Pads)

Immediate hip protection.

- 3/4 length unisex pants with sewn-in hip protection pad •
- Immediate hip protection for new admissions, acute care and short stay patients
- Residents with temporary changes in risk factors
- As a 'stop gap' component of your Falls Safety Programme while you order regular HipSaver for ongoing hip protection
- Suitable for both male and female users
- Colour grey



Size	Interim Overpant - Unisex	Interim Overpant - With Tailbone Protection - Unisex
XS	PTA534107	PTA534108
S	PTA534107S	PTA534108S
Μ	PTA534107M	PTA534108M
L	PTA534107L	PTA534108L
XL	PTA534107XL	PTA534108XL
XXL	PTA534107XXL	PTA534108XXL

HipSaver - Track Pants (Sewn-In Pads)

Stylish hip protection for active seniors (Unisex)

- HipSaver Track Pants have the same all-soft hip protectors as other HipSaver garments
- Smartly designed in classic navy blue & pre-shrunk to fit properly (washable with regular laundry)
- The loose ankle hem prevents legs from binding and provides added comfort
- 2 lengths regular (height under 170 cm), long (over 170 cm)
- Suitable for both male and female users

Options:

- Hip Protection with Tailbone Padding (protects against backward falls) •
- Hip Protection with Knee Pads (proper positioning and protects against forward falls)

	0.75.0		\checkmark	DVA: AR28			
Size	Track Pants	Track Pants with Knee Pad	Track Pants with Tailbone Pad	Track Pants with Tailbone & Knee Pad			
	Regular (User Height Under 170cm)						
XS	PTA534109	PTA534111	PTA534110	PTA534112			
S	PTA534109RS	PTA534111S	PTA534110RS	PTA534112S			
Μ	PTA534109RM	PTA534111M	PTA534110RM	PTA534112M			
L	PTA534109RL	PTA534111L	PTA534110RL	PTA534112L			
XL	PTA534109RXL	PTA534111XL	PTA534110RXL	PTA534112XL			
XXL	PTA534109RXXL	PTA534111XXL	PTA534110RXXL	PTA534112XXL			
		Long (User Height C	Over 170cm)				
XS	PTA534109LXS	PTA534111LXS	PTA534110LXS	PTA534112LXS			
S	PTA534109LS	PTA534111LS	PTA534110LS	PTA534112LS			
Μ	PTA534109LM	PTA534111LM	PTA534110LM	PTA534112LM			
L	PTA534109LL	PTA534111LL	PTA534110LL	PTA534112LL			
XL	PTA534109LXL	PTA534111LXL	PTA534110LXL	PTA534112LXL			
XXL	PTA534109LXXL	PTA534111LXXL	PTA534110LXXL	PTA534112LXXL			



DVA Contact Details



HipSaver - Shorts (Sewn-In Pads) - Unisex

Immediate hip protection.

- HipSaver Shorts have the same all-soft hip protectors as other HipSaver garments
- Smartly designed in classic navy blue
- Pre-shrunk to fit properly, and washable with regular laundry
- Six hip sizes offer a perfect fit for any active body
- Also available with tailbone protection which offers additional security for those who tend to fall backwards

		DVA. ARZo
Size	Shorts - Unisex	Shorts with Tailbone Protection - Unisex
XS	PTA534113	PTA534114
S	PTA534113S	PTA534114S
Μ	PTA534113M	PTA534114M
L	PTA534113L	PTA534114L
XL	PTA534113XL	PTA534114XL
XXL	PTA534113XXL	PTA534114XXL

Hornsby Comfy Hips

Overview

• Provides protection from injuries caused by falls. Australian designed, owned and manufactured

To Order

Use the Size Guide to select your size, then choose your preferred model and the appropriate code

Size Guide

Size	Hips (cm)
XXS	50 - 60
XS	60 - 70
S	70 - 80
Μ	80 - 90
L	90 - 100
XL	100 - 110
XXL	110 - 120

 $\ensuremath{\textbf{Measure}}$ around the waist when standing (cm)

- Available with either Sewn-In Shields, or as pants with removable shields (shields sold separately)
- Available in either unisex styles, or male style (with fly opening)
- A Real
 - Comfortable fabric 97% cotton
 - Range of 7 sizes with durable woven elastic waist band

DVA: AR28

Size	Comfy Hip (Sewn-In Shields)		Comfy Hip (Pants Only)	
	Code - Unisex		Code - Unisex	Code - Male
XXS	PTA512448	PTA512449	PTA512450	PTA512451
XS	PTA512448XS	PTA512449XS	PTA512450XS	PTA512451XS
S	PTA512448S	PTA512449S	PTA512450S	PTA512451S
Μ	PTA512448M	PTA512449M	PTA512450M	PTA512451M
L	PTA512448L	PTA512449L	PTA512450L	PTA512451L
XL	PTA512448XL	PTA512449XL	PTA512450XL	PTA512451XL
XXL	PTA512448XXL	PTA512449XXL	PTA512450XXL	PTA512451XXL
One Size	Protector Shields - Comfy Hip (Pair)	PTA5	12452



PRODUCT AND SERVICES

Protection -Skin & Limb

ndis

 \checkmark

<



Overview

Protects from pressure, friction and minor traumas while providing ample airflow •

To Order

Choose your preferred model, select your size and order using the appropriate code

DermaSaver - Hand & Finger Garments

				F	
Size	Knuckle Protector	Dermamitt	Palm Pillow	Finger Separator (One Size Fits All)	
S/M	PTA533301SM	PTA533302SM	PTA533303SM	DTAF22200	
L/XL	PTA533301LXL	PTA533302LXL	PTA533303LXL	PTA533300	

DermaSaver - Forearm & Elbow Garments



Measurement: All measurements are circumferences







DVA: AR29

For NDIS Item Code, refer to page: 500



(]

DVA: AR29

DermaSaver - Arm Garments

masaver - Ann Gannenis			(DVA: AR29
				i la
	Forearm (cm) (widest part)	Style	Arm Tube	Arm Tube (with Knuckle Protector)
	13 - 18	with Single Elbow	PTA533308	PTA533311
		with Double Elbow	PTA533010	FTA033311
	40.05	with Single Elbow	PTA533308S	PTA533311S
	18 - 25	with Double Elbow	PTA533020	F 1A3333115
	25 - 30	with Single Elbow	PTA533308M	PTA533311M
	25 - 30	with Double Elbow	PTA533030	F TA9533 T TW
	30 - 38	with Single Elbow	PTA533308L	PTA533311L
	30 - 38	with Double Elbow	PTA533040	FIA003311L
	38 - 46	with Single Elbow	PTA533308XL	PTA533311XL
38 - 46		with Double Elbow	PTA533050	FIA5533TIAL
rer	nent: All measurements are circ	umferences		

Measurement: All measurements are circumferences

XS

S

Μ

L

XL

DermaSaver - Leg Garments

						5	
Size	Calf (cm)	Style (if applicable)	Leg Tube	Knee Tube	Shin Tube	Shin to Knee Tube	
XS	XS 18 - 23	with Single Knee	PTA533317	PTA533320	PTA533322	PTA533323	
A3	10 - 23	with Double Knee	PTA533318	PTA533321	T 1A555522	F 1A555525	
S	23 - 30	with Single Knee	PTA533317S	PTA533320S	PTA533322S	PTA533323S	
3	23 - 30	with Double Knee	PTA533318S	PTA533321S	F 1A5555225		
M	20. 20	with Single Knee	PTA533317M	PTA533320M	PTA533322M	DTAF22222M	
IVI	30 - 38	with Double Knee	PTA533318M	PTA533321M	PTA533322IVI	PTA533323M	
L		with Single Knee	PTA533317L	PTA533320L	PTA533322L	PTA533323L	
L	38 - 48	with Double Knee	PTA533318L	PTA533321L	F1A033322L	F1A533323L	
XL	48 - 58	with Single Knee	PTA533317XL	PTA533320XL	PTA533322XL	PTA533323XL	
AL	48 - 98	with Double Knee	PTA533318XL	PTA533321XL	FIA5333ZZAL	P 1A533323AL	

Measurement: All measurements are circumferences

DermaSaver - Foot & Heel Garments

 \checkmark Stay-Put Heel Protector XS < 28 PTA533315 PTA533313 S 28 - 30 PTA533315S PTA533313S PTA533316 Μ 30 - 35 PTA533315M PTA533313M 35 - 38 PTA533315L PTA533313L PTA533316L L XL > 39 PTA533315XL PTA533313XL

Measurement: All measurements are circumferences



DermaSaver - Ankle Garments



\checkmark	DVA: AR29				
Size	Ankle (cm)	Ankle Bumper			
S	< 23	PTA533314			
Μ	23 - 28 PTA533314M				
L	28 - 33	PTA533314L			
XL	> 33	PTA533314XL			

Measurement: All measurements are circumferences

Protect-A-Limb - Limb & Skin Protectors



DermaSaver - Amputee Stump Cover



Measurement: Around the leg, 5cm above amputation

Overview & Features:

Protect-A-Limb® protects the limbs from skin tears and skin damage.

- Protection against friction and impact
- Durable and comfortable to wear
- Constructed of Coolmax® fabric on the inner surface, and slightly retentive cotton on the outer surface
- Moisture wicking and prevents maceration
- Sealed ends prevent fraying
- Supplied as pairs

Size	Code	Size Guide
S	PTA491600	Suits arm and small leg
Μ	PTA491700	Suits arm and average lower leg
L	PTA491800	Full leg length

Rolyan

Overview

- Shock-absorbing gel padding: Effective following sudden injury or cumulative trauma injury to the elbow or heel
- Helps protect the joint by diffusing and distributing pressure
- Includes 6.4 mm shock-absorbing gel pad, which won't bottom out
- Sold in pairs

To Order

Choose your preferred model, select your size and order using the appropriate code. For borderline measurements, choose the larger size for extra comfort.

Heel/Elbow Protector - Gel



		DVA: AR29	
Size	Code	Size Guide	
S	PTA506400	Fits most small and medium women	
M/L	PTA506500	Fits larger women and average-sized men	
XL	PTA506600	Fits larger men	

lolyaya n





GelBodies

Overview

- Gel provides pressure-relieving skin protection and protection from pressure injury
- Uses a pressure-relieving gel found in high-end wheelchair cushions
- Cotton-lycra "coolband" stretches for ideal fit and feels great against the skin
- Machine washable
- Supplied as pairs

To Order

Choose your preferred model, select your size and order using the appropriate code. For borderline measurements, choose the larger size for extra comfort.

GelBodies - Heel & E	DVA: AR29		
The second s	Size	Heel or Elbow (cm)	Code
	S	20 - 25	PTA533324
	Μ	25 - 30	PTA533324M
	L	30 - 36	PTA533324L
	XL	36 - 41	PTA533324XL

Measurement: All measurements are circumferences

GelBodies - Knee			DVA: AR29
	Size	Knee (cm)	Code
	S	23 - 28	PTA533325
	Μ	28 - 33	PTA533325M
	L	33 - 41	PTA533325L

Measurement: All measurements are circumferences

GelBodies - Calf	DVA: AR29		
	Size	Calf (cm)	Code
	Μ	25 - 33	PTA533326
	L	30 - 41	PTA533326L

Measurement: All measurements are circumferences

GelBodies - Ankle			DVA: AR29
(Annual States)	Size	Ankle (cm)	Code
A A	S	20 - 25	PTA533327
	Μ	25 - 30	PTA533327M
	L	30 - 36	PTA533327L

Measurement: All measurements are circumferences







Compression - Legwear

ndis

For NDIS Item Code refer to page: 500



- Improve blood flow in the lower limbs and assist in the prevention and management of venous disease •
- Widely used to treat fluid build-up associated with lymphoedema
- Lymphoedema occurs if the lymphatic system is damaged or blocked. Excess fluid, which would normally be drained • away by the lymphatics, builds up in the tissues (most commonly the arms or legs)
- Compression plays a major role in the treatment of lymphoedema by encouraging lymph flow

Compression Values

- Medical stockings feature "graduated compression" - strongest at the ankle & progressively reducing up the leg
- Compression is measured in mmHg (millimetres of mercury compression) at the ankle
- Stockings < 20mmHg (preventative): suitable for general leg relief, mild varicose veins, safeguard against DVT
- Stockings > 20mmHg (therapeutic): more-serious or pre-existing conditions
- Stockings > 30 mmHg: medical professional prescription required

PREVENTATIVE		THERAPEUTIC	
15 - 20mmHg	20 - 30mmHg (class 1)	30 - 40mmHg (class 2)	40-50 mmHg (class 3)
 Moderate Tired, aching legs Mild varicose veins Mild varicose veins during pregnancy Mild fluid retention Post sclerotherapy Travel - for those with otherwise healthy legs 	 Firm Mild/moderate varicose veins Varicosities tending towards mild oedema Varicose veins during and after pregnancy Prevention of thrombosis (ambulatory) Post sclerotherapy 	 Extra Firm Mild chronic venous insufficiency Moderate/severe varicosities during pregnancy Phlebitis with mild oedema After deep vein thrombosis Diffuse varicosities Varicose ulcers Mild lymphoedema 	 Strong Pronounced varicosities with oedema Pronounced chronic venous insufficiency After deep thrombophlebitis Valvular insufficiency of the deep veins Varicose ulcers & lymphoedema
Contraindications	Arteriopathy, hypodermatitis in acute thrombophlebitis.	state, wet dermatosis, cutaneous infecti	ons, unreduced oedema, acute

449

DVA Contact Details



I ndis

For NDIS Item Code, refer to page: 500 to page



Compression - Legwear Pressure Care & Circulation



	DVA #	DVA Approval Needed?	DVA Approval Rules	Eligible Prescriber*
	AR22	No (unless Aidacare advises)	Limit of 3 pairs every 6 months	S, RN, Physio, OT, LMO
D\	/A Comments	*For lymphoedema garments, health providers shoul management in order to access, measure, fit and rev		ema/lymphoedema

Styles Available

- Closed-toe (full-foot, like a sock) Most common. Appropriate with swelling of the foot/toes
- Open-toe (toes exposed) Suitable where large foot size relative to leg measurement



Fitting & Selection Tips

- Correct measurement and sizing ARE VERY important
- Measurements should be taken first thing in the morning
- Cotton is best for sensitive skin
- Knee High Ensure the top-band is flat against the leg, about 2 3cm below the knee crease. Do not allow the stocking to roll over at the top and never double the stocking over at the top if it seems too long
- Sizing Selection Select size according to the largest measurement taken (eg. if ankle is 'small' and calf is 'medium', select the medium). If more than one size difference, contact Aidacare.

Refer To Pages 480 For Donning And Doffing Aids



Sigvaris Doff N' Donner PTS642660



Jobst Donning Gloves FBS303490



Easy Slide 'Magnide' On /Off Aid (Closed Toe) PTS642512



PRODUCT AND SERVICES CATALOGUE

Dr. Comfort Diabetic Support Socks - Unisex (15 - 20mmHg)

- Diabetic support socks for men and women
- Helps to alleviate early lower limb conditions
- Prevents swelling in the legs, ankles and feet
- Compression: 15 20mmHa

To Order

Choose your preferred model, then use the Size Guide to select the appropriate size and code

	Size Guide									
Size	Ankle (cm)	Calf (cm)	Men's Shoe Size	Women's Shoe Size						
Small	18 - 21	28 - 38	6 - 7.5	5 - 7						
Medium	21 - 25	30 - 42	8 - 10	7.5 - 9.5						
Large	25 - 29	32 - 46	10.5 - 12.5	10 - 12						
X Large	29 - 33	34 - 50	13+	12.5+						

Measurement: All measurements are circumferences. Shoe Size - USA sizes

		DVA: AR22
Size	Colour	Code
Small	Black	PTS642502
Sman	White	PTS642506
Medium	Black	PTS642501
wearum	White	PTS642505
Lorgo	Black	PTS642500
Large	White	PTS642504
Vlarga	Black	PTS642503
X Large	White	PTS642507

GloriaMed - Medical Legwear

Overview

Made in Italy, GloriaMed is one of the oldest and most respected manufacturers of compression legwear

To Order

Choose your preferred model, then use the Size Guide to select the appropriate size and code

GloriaMed - Stockings - Unisex

GloriaMed stockings offer a choice of a regular fit or wide fit (good for when there is increased / excessive swelling)



DVA Contact Details

			C C	DVA: A	AR22						
Size Guide											
Size	Ankle (cm)	Below	/ Knee	Upper	Upper Thigh						
5120	Ankie (cm)	Regular Fit (cm)	Wide Fit (cm)	Regular Fit (cm)	Wide Fit (cm)						
X Small	18 -19	28 - 30	31 - 34	42 - 46	47 - 51						
Small	20 - 22	31 - 34	35 - 38	47 - 51	52 - 57						
Medium	23 - 25	35 - 38	39 - 42	52 - 57	58 - 63						
Large	26 - 28	39 - 42	43 - 47	58 - 63	64 - 70						
X Large	29 - 31	43 - 47	48 - 52	64 - 70	71 - 78						

Measurement: All measurements are circumferences. Ankle - Narrowest part (just above ankle bone). Below knee - 2cm below knee. Upper Thigh - 5cm below groin Note: Thigh measurement relevant to thigh style only





P 1300 888 052 | F 1300 787 052 | E dva@aidacare.com.au







451



DVA: AR22

DVA: AR22

GloriaMed - Stockings - Unisex - Knee (Regular Fit)

- Available in soft, comfortable nylon OR hypoallergenic 100% cotton
- A choice of short (up to 38cm) or long (39cm+) lengths for a better fit
- Available in **beige only**

 Availab 	ie in beige only	¥				JVA. ANZZ
		20 - 30)mmHg	30 - 40)mmHg	40 - 50mmHg
		Open Toe	Closed Toe	Open Toe	Closed Toe	Open Toe
Size	Material			Short (up to 38cm)		
X Small	Nylon	PTS641860XSRS	PTS641861XSRS	PTS641865XSRS	PTS641866XSRS	PTS641870XSRS
∧ SIIIdii	Cotton	PTS641873XSRS	-	PTS641875XSRS	-	-
Small	Nylon	PTS641860SRS	PTS641861BE1SRS	PTS641865SRS	PTS641866BE1SRS	PTS641870SRS
Silidii	Cotton	PTS641873SRS	-	PTS641875SRS	-	-
Medium	Nylon	PTS641860MRS	PTS641861BE2MRS	PTS641865MRS	PTS641866BE2MRS	PTS641870MRS
Wedium	Cotton	PTS641873MRS	-	PTS641875MRS	-	-
Lorgo	Nylon	PTS641860LRS	PTS641861BE3LRS	PTS641865LRS	PTS641866BE3LRS	PTS641870LRS
Large	Cotton	PTS641873LRS	-	PTS641875LRS	-	-
X Large	Nylon	PTS641860XLRS	PTS641861BE4XLRS	PTS641865XLRS	PTS641866BE4XLRS	PTS641870XLRS
A Large	Cotton	PTS641873XLRS	-	PTS641875XLRS	-	-
Size	Material			Long (39cm+)		
X Small	Nylon	PTS641860XSRL	PTS641861XSRL	PTS641865XSRL	PTS641866XSRL	PTS641870XSRL
X Smail	Cotton	PTS641873XSRL	-	PTS641875XSRL	-	-
Small	Nylon	PTS641860SRL	PTS641861BE1SRL	PTS641865SRL	PTS641866BE1SRL	PTS641870SRL
SIIIdii	Cotton	PTS641873SRL	-	PTS641875SRL	-	-
Medium	Nylon	PTS641860MRL	PTS641861BE2MRL	PTS641865MRL	PTS641866BE2MRL	PTS641870MRL
Medium	Cotton	PTS641873MRL	-	PTS641875MRL	-	-
Lorgo	Nylon	PTS641860LRL	PTS641861BE3LRL	PTS641865LRL	PTS641866BE3LRL	PTS641870LRL
Large	Cotton	PTS641873LRL	-	PTS641875LRL	-	-
VLarga	Nylon	PTS641860XLRL	PTS641861BE4XLRL	PTS641865XLRL	PTS641866BE4XLRL	PTS641870XLRL
X Large	Cotton	PTS641873XLRL	-	PTS641875XLRL	-	-

GloriaMed - Stockings - Unisex - Knee (Wide Fit)

- Available in soft, comfortable nylon OR hypoallergenic 100% cotton
- A choice of short (up to 38cm) or long (39cm+) lengths for a better fit
- Wide fit will better accommodate increased / excessive swelling
- Available in **beige only**

Availab	ie in beige on	y				
		20 - 30	lmmHg	30 - 40	lmmHg	40 - 50mmHg
		Open Toe	Closed Toe	Open Toe	Closed Toe	Open Toe
Size	Material			Short (up to 38cm)		
V Creatil	Nylon	PTS641860XSWS	PTS641861XSWS	PTS641865XSWS	PTS641866XSWS	PTS641870XSWS
X Small	Cotton	PTS641873XSWS	-	PTS641875XSWS	-	-
Small	Nylon	PTS641860SWS	PTS641861BE1SWS	PTS641865SWS	PTS641866BE1SWS	PTS641870SWS
Sman	Cotton	PTS641873SWS	-	PTS641875SWS	-	-
Medium	Nylon	PTS641860MWS	PTS641861BE2MWS	PTS641865MWS	PTS641866BE2MWS	PTS641870MWS
weatum	Cotton	PTS641873MWS	-	PTS641875MWS	-	-
Lavaa	Nylon	PTS641860LWS	PTS641861BE3LWS	PTS641865LWS	PTS641866BE3LWS	PTS641870LWS
Large	Cotton	PTS641873LWS	-	PTS641875LWS	-	-
X Lawrence	Nylon	PTS641860XLWS	PTS641861BE4XLWS	PTS641865XLWS	PTS641866BE4XLWS	PTS641870XLWS
X Large	Cotton	PTS641873XLWS	-	PTS641875XLWS	-	-
Size	Material			Long (39cm+)		
V. Cara all	Nylon	PTS641860XSWL	PTS641861XSWL	PTS641865XSWL	PTS641866XSWL	PTS641870XSWL
X Small	Cotton	PTS641873XSWL	-	PTS641875XSWL	-	-
Small	Nylon	PTS641860SWL	PTS641861BE1SWL	PTS641865SWL	PTS641866BE1SWL	PTS641870SWL
Sman	Cotton	PTS641873SWL	-	PTS641875SWL	-	-
Medium	Nylon	PTS641860MWL	PTS641861BE2MWL	PTS641865MWL	PTS641866BE2MWL	PTS641870MWL
weatum	Cotton	PTS641873MWL	-	PTS641875MWL	-	-
Lorgo	Nylon	PTS641860LWL	PTS641861BE3LWL	PTS641865LWL	PTS641866BE3LWL	PTS641870LWL
Large	Cotton	PTS641873LWL	-	PTS641875LWL	-	-
VLorgo	Nylon	PTS641860XLWL	PTS641861BE4XLWL	PTS641865XLWL	PTS641866BE4XLWL	PTS641870XLWL
X Large	Cotton	PTS641873XLWL	-	PTS641875XLWL	-	-



DVA: AR22

DVA: AR22

GloriaMed - Stockings - Unisex - Thigh (Regular Fit)

- Available in soft comfortable nylon
- A choice of short (up to 71cm) or long (72cm+) lengths for a better fit
- Available in **beige only** •

				OPEN TO	DE ONLY			
		20 - 30mmHg			30 - 40mmHg		40 - 50	mmHg
	Thigh high (Grip top)	Thigh Belt Top (Left)	Thigh Belt Top (Right)	Thigh high (Grip top)	Thigh Belt Top (Left)	Thigh Belt Top (Right)	Thigh Belt Top (Left)	Thigh Belt Top (Right)
Size				Short (up	to 71cm)			
XS	PTS641863XSRS	PTS641864LXSRS	PTS641864RXSRS	PTS641868XSRS	PTS641869LXSRS	PTS641869RXSRS	PTS641871LXRS	PTS641871RXSRS
S	PTS641863SRS	PTS641864LSRS	PTS641864RSRS	PTS641868SRS	PTS641869LSRS	PTS641869RSRS	PTS641871LSRS	PTS641871RSRS
Μ	PTS641863MRS	PTS641864LMRS	PTS641864RMRS	PTS641868MRS	PTS641869LMRS	PTS641869RMRS	PTS641871LMRS	PTS641871RMRS
L	PTS641863LRS	PTS641864LLRS	PTS641864RLRS	PTS641868LRS	PTS641869LLRS	PTS641869RLRS	PTS641871LLRS	PTS641871RLRS
XL	PTS641863XLRS	PTS641864LXLRS	PTS641864RXLRS	PTS641868XLRS	PTS641869LXLRS	PTS641869RXLRS	PTS641871LXLRS	PTS641871RXLRS
Size				Long (72cm+)			
XS	PTS641863XSRL	PTS641864LXSRL	PTS641864RXSRL	PTS641868XSRL	PTS641869LXSRL	PTS641869RXSRL	PTS641871LXSRL	PTS641871RXSRL
S	PTS641863SRL	PTS641864LSRL	PTS641864RSRL	PTS641868SRL	PTS641869LSRL	PTS641869RSRL	PTS641871LSRL	PTS641871RSRL
Μ	PTS641863MRL	PTS641864LMRL	PTS641864RMRL	PTS641868MRL	PTS641869LMRL	PTS641869RMRL	PTS641871LMRL	PTS641871RMRL
L	PTS641863LRL	PTS641864LLRL	PTS641864RLRL	PTS641868LRL	PTS641869LLRL	PTS641869RLRL	PTS641871LLRL	PTS641871RLRL
XL	PTS641863XLRL	PTS641864LXLRL	PTS641864RXLRL	PTS641868XLRL	PTS641869LXLRL	PTS641869RXLRL	PTS641871LXLRL	PTS641871RXLRL

GloriaMed - Stockings - Unisex - Thigh (Wide Fit)

- Available in soft comfortable nylon
- A choice of short (up to 71cm) or long (72cm+) lengths for a better fit •
- Wide fit will better accommodate increased / excessive swelling
- Available in beige only •

	OPEN TOE ONLY								
	20 - 30mmHg				30 - 40mmHg		40 - 50mmHg		
	Thigh high (Grip top)	Thigh Belt Top (Left)	Thigh Belt Top (Right)	Thigh high (Grip top)	Thigh Belt Top (Left)	Thigh Belt Top (Right)	Thigh Belt Top (Left)	Thigh Belt Top (Right)	
Size				Short (up	to 71cm)				
XS	PTS641863XSWS	PTS641864LXSWS	PTS641864RXSWS	PTS641868XSWS	PTS641869LXSWS	PTS641869RXSWS	PTS641871LXWS	PTS641871RXSWS	
S	PTS641863SWS	PTS641864LSWS	PTS641864RSWS	PTS641868SWS	PTS641869LSWS	PTS641869RSWS	PTS641871LSWS	PTS641871RSWS	
Μ	PTS641863MWS	PTS641864LMWS	PTS641864RMWS	PTS641868MWS	PTS641869LMWS	PTS641869RMWS	PTS641871LMWS	PTS641871RMWS	
L	PTS641863LWS	PTS641864LLWS	PTS641864RLWS	PTS641868LWS	PTS641869LLWS	PTS641869RLWS	PTS641871LLWS	PTS641871RLWS	
XL	PTS641863XLWS	PTS641864LXLWS	PTS641864RXLWS	PTS641868XLWS	PTS641869LXLWS	PTS641869RXLWS	PTS641871LXLWS	PTS641871RXLWS	
Size				Long (72cm+)				
XS	PTS641863XSWL	PTS641864LXSWL	PTS641864RXSWL	PTS641868XSWL	PTS641869LXSWL	PTS641869RXSWL	PTS641871LXSWL	PTS641871RXSWL	
S	PTS641863SWL	PTS641864LSWL	PTS641864RSWL	PTS641868SWL	PTS641869LSWL	PTS641869RSWL	PTS641871LSWL	PTS641871RSWL	
Μ	PTS641863MWL	PTS641864LMWL	PTS641864RMWL	PTS641868MWL	PTS641869LMWL	PTS641869RMWL	PTS641871LMWL	PTS641871RMWL	
L	PTS641863LWL	PTS641864LLWL	PTS641864RLWL	PTS641868LWL	PTS641869LLWL	PTS641869RLWL	PTS641871LLWL	PTS641871RLWL	
XL	PTS641863XLWL	PTS641864LXLWL	PTS641864RXLWL	PTS641868XLWL	PTS641869LXLWL	PTS641869RXLWL	PTS641871LXLWL	PTS641871RXLWL	

DVA Contact Details



For NDIS Item Code, refer to page: 500



Compression - Legwear Pressure Care & Circulation



GloriaMed - Micro Mesh - Pantyhose - Women

- 入
- Fashionable and sheer
- Suits all standard body shapes
- Closed toe only

Available in beige or black

 Two compression classes available (20 - 30mmHg or 30 - 40mmHg)

270 C						DVA: AR22
		Size Guide		20 - 30mmHg	30 - 40mmHg	
Size	Ankle (cm)	Below Knee (cm)	Upper Thigh (cm)	Colour	Code	
Small	20 - 21	31 - 34	52 - 57	Beige	PTS641877	PTS641878
Sman	20 - 21	31 - 34	52 - 57	Black	PTS641877BL0	PTS641878BL0
Medium	22 - 23	33 - 36	55 - 60	Beige	PTS641877BE1	PTS641878BE1
Medium	22 - 23	33 - 30	55 - 60	Black	PTS641877BL1	PTS641878BL1
Louno	04 05	25 20	58 - 63	Beige	PTS641877BE2	PTS641878BE2
Large	24 - 25	35 - 38	58 - 63	Black	PTS641877BL2	PTS641878BL2
X Lorgo	26 - 27	37 - 40	61 - 66	Beige	PTS641877BE3	PTS641878BE3
X Large	20 - 27	37 - 40	01-00	Black	PTS641877BL3	PTS641878BL3

Measurement: All measurements are circumferences. Ankle - Narrowest part (just above ankle bone). Below Knee - 2cm below knee. Upper Thigh - 5cm below groin

GloriaMed - Lady Gloria - Thigh - Women

- Attractive sheer and comfortable
- Suitable for young and dynamic women
- Closed toe only

- Available in **beige or black**
- Compression class available (15 20mmHg)

 \checkmark

15 - 20mmHg Code

PTS641881

PTS641881BL1

PTS641881BE2

PTS641881BL2

PTS641881BE3

PTS641881BL3

PTS641881BE4

PTS641881BL4

PTS641881BE5

PTS641881BL5





DVA: AR22

Measurement: Determine size from the height / weight grid

Beige

Black

Beige

Black

Beige

Black

Beige

Black

Beige

Black

1

2

3

4

5

General Enquiries

IBICI - Medical Legwear

Overview

• Designed and manufactured in Italy, IBICI's compression legwear range is ideal for mild vascular conditions, or as a preventative measure. All IBICI garments are **closed toe (full foot) only**

To Order

Choose your preferred model, then use the Size Guide to select the appropriate size and code

Segreta 140 & Segreta PLUS - Knee High - Women

- Segreta 140: Fashionably sheer, in a range of colours and styles. Closed toe only. Compression: 15 - 20mmHg
- **Segreta PLUS:** Comparable to 20 30, these give additional compression at the ankle. **Closed toe only.** Compression: 23 27mmHg

					Segreta 140 (Knee high)	Segreta PLUS (Knee high)
		Size Guide		15 - 20mmHg	23 - 27mmHg	
Size	Shoe Size	Ankle (cm)	Calf (cm)	Colour	Co	de
1	34 - 35	17 - 20	28 - 38	Beige	PTS641879	PTS641885
1	34 - 35	17 - 20	20 - 30	Black	PTS641879BL1	-
2	35 - 36	19 - 22	20. 20	Beige	PTS641879BE2	PTS6418852
2	35 - 30	19 - 22	30 - 39	Black	PTS641879BL2	-
3	27 20	20 - 24	31 - 40	Beige	PTS641879BE3	PTS6418853
3	37 - 38	20 - 24	31 - 40	Black	PTS641879BL3	-
4	39 - 40	21 - 26	32 - 41	Beige	PTS641879BE4	PTS6418854
4	39 - 40	21-20	32 - 41	Black	PTS641879BL4	-
-	44 40	00 00	22 42	Beige	PTS641879BE5	PTS6418855
5	41 - 42	22 - 28	33 - 42	Black	PTS641879BL5	-
C	40 40	22 20	24 42	Beige	PTS641879BE6	-
6	42- 43	23 - 29	34 - 43	Black	PTS641879BL6	

Measurement: Shoe Size - Metric European. All other measurements are circumferences. Ankle - Narrowest part (just above ankle bone). Calf - Widest part

Segreta 140 & Segreta PLUS - Pantyhose - Women

- Segreta 140: Fashionably sheer, in a range of colours and styles. Closed toe only. Compression: 15 20mmHg. Available in 5 regular sizes and 3 wider sizes (2W, 3W, 4W)
- Segreta PLUS: Comparable to 20 30, these give additional compression at the ankle. The Segreta PLUS are semi-sheer. Closed toe only. Compression: 23 27mmHg. Available in 5 regular sizes only
- Maternity pantyhose are not on the DVA/RAP schedule





Measurement:

- First, determine size from the height/weight grid
- Then check that size against the ankle measurements on the
- Size Guide on the previous page

- If one measurement indicates a larger size than the other, select the larger size
- If more than 1 size difference, contact Aidacare for assistance

455



DVA: AR22

PRODUCT AND SERVICES







Segreta 140 & Segreta PLUS - Pantyhose - Women (contd)

				DVA: AR22
		Segreta 140 (Pantyhose - Regular Fit)	Segreta 140 (Pantyhose - Wider Fit)	Segreta PLUS (Pantyhose)
		15 - 20mmHg	15 - 20mmHg	23 - 27mmHg
Size	Colour		Code	
1	Beige	PTS641882	-	PTS641887
1	Black	PTS641882BL1	-	-
0	Beige	PTS641882BE2	PTS641883BE2	PTS6418872
2	Black	PTS641882BL2	-	-
0	Beige	PTS641882BE3	PTS641883BE3	PTS6418873
3	Black	PTS641882BL3	-	-
4	Beige	PTS641882BE4	PTS641883BE4	PTS6418874
4	Black	PTS641882BL4	-	-
F	Beige	PTS641882BE5	-	PTS6418875
5	Black	PTS641882BL5	-	-

SKYWALK - Knee High - Unisex (15 - 20mmHg)

- Designed as a travel sock, but also suitable for day-to-day wear for the prevention and relief of mild vascular conditions
- Similar to REPOMEN in appearance and compression, but featuring a high cotton content (33%) for added softness and comfort

				DVA: AR22
	Size Gui	de		15 - 20mmHg
Size	Ankle (cm)	Calf (cm)	Colour	Code
X Small	18 - 21	27 - 37	Black	PTS641889
Small	20 - 23	30 - 40	Black	PTS641889S
Medium	22 - 25	33 - 43	Black	PTS641889M
Large	24 - 27	36 - 46	Black	PTS641889L
X Large	26 - 29	39 - 49	Black	PTS641889XL

Measurement: All measurements are circumferences. Ankle - Narrowest part (just above ankle bone). Calf - Widest part

REPOMEN - Knee High - Men (15 - 20mmHg)

- Has the stylish appearance of dress socks
- Two colours to choose from (beige or black)
- Compression: 15 20mmHg

General Enquiries



					DVA: AR22
		15 - 20mmHg			
Size	Shoe Size	Ankle (cm)	Calf (cm)	Colour	Code
X Small	37 - 38	18 - 21	30 - 40	Beige	PTS641888
X Sinan	Silidii 37 - 30	10 - 21	50 - 40	Black	PTS641888BLXS
Small	mall 39 - 40	20 - 23	32 - 42	Beige	PTS641888BES
Sman		20 - 23		Black	PTS641888BLS
Medium	41 - 42	21 - 25	33 - 43	Beige	PTS641888BEM
Wearum	41-42	21-25		Black	PTS641888BLM
Large	43 - 44	22 - 27	34 - 46	Beige	PTS641888BEL
Large	45 - 44	22 - 21	54 - 40	Black	PTS641888BLL
X Large	45 - 46	23 - 28	35 - 48	Beige	PTS641888BEXL
A Large	45 - 40		55 - 40	Black	PTS641888BLXL

Measurement: Shoe Size - Metric European. All other measurements are circumferences. Ankle - Narrowest part (just above ankle bone). Calf - Widest part

Jobst - Medical Legwear

Overview

Jobst stockings focus on improving the guality of life for those suffering from venous and lymphatic diseases

To Order

Choose your preferred model, then use the Size Guide to select the appropriate size and code

Size Guide

KNEE HIGH:

Sizing is the same for all Jobst knee high stockings, but there are an additional 2 larger sizes in the 'Relief' knee high range.

Take the ankle and calf measurements to determine size:

- Ankle the narrowest point of the ankle, immediately above the ankle bone
- Calf the widest part of the lower leg

	Small	Medium	Large	X Large	'Relief' Only	'Relief' Only
Ankle (cm)	18 - 21	21 - 25	25 - 29	29 - 33	25 - 30	30 - 36
Calf (cm)	28 - 38	30 - 42	32 - 46	34 - 50	46 - 61	46 - 61
Groin (cm)	40 - 62	46 - 70	54 - 78	60 - 81	-	-
Hip (cm)	71 - 117	76 - 127	81 - 137	102 - 166	-	-

THIGH & WAIST HIGH:

the leg (groin)

hip measurement

In addition to ankle and calf measurements.

Take the measurement around the top of

The waist high (pantyhose), also require a

Measurement: All measurements are circumferences.

Jobst UltraSheer for Women

- Sheer, comfortable and shaped for women. These are **closed toe only** stockings
- The UltraSheer range is available in knee high, thigh high and waist high styles

WELET AND	AND IN COLUMN	CLOSED TOE ONLY									
201	and the second	15 - 20mmHg				20 - 30mmHg			30 - 40mmHg		
		Knee high	Thigh high (Grip top)	Waist high (Panty- hose)	Knee high	Thigh high (Grip top)	Waist high (Pantyhose)	Knee high	Thigh high (Grip top)	Waist high (Pantyhose)	
Size	Colour					Code					
Small	Black	PTS640140	PTS640330B	PTS640380	PTS640610	PTS640710	PTS640779	PTS640870	PTS640945	PTS640983	
SIIIdii	Natural	PTS640220	PTS640330	PTS640450	PTS640650	PTS640740	PTS640789	PTS640910	PTS640949	PTS640969	
Medium	Black	PTS640150	PTS640300	PTS640390	PTS640620	PTS640720	PTS640780	PTS640880	PTS640946	PTS640984	
weulum	Natural	PTS640230	PTS640335	PTS640460	PTS640660	PTS640750	PTS640790	PTS640920	PTS640950	PTS640970	
Lorgo	Black	PTS640160	PTS640310	PTS640400	PTS640630	PTS640730	PTS640781	PTS640890	PTS640947	PTS640985	
Large	Natural	PTS640240	PTS640350	PTS640470	PTS640670	PTS640760	PTS640800	PTS640930	PTS640955	PTS640980	
Viargo	Black	PTS640170	PTS640315	PTS640410	PTS640640	PTS640735	PTS640782	PTS640900	PTS640948	PTS640986	
X Large	Natural	PTS640250	PTS640360	PTS640480	PTS640680	PTS640761	PTS640810	PTS640940	PTS640960	PTS640981	

Jobst for Men

- Comfortable, lightweight and stylish ribbed design. These are **closed-toe only** stockings of nylon fabric
- The Jobst for Men range is available in knee high or thigh high styles

19225			CLOSED TOE ONLY				
1		15 - 20mmHg	20 - 30	mmHg	30 - 40	mmHg	
0.		Knee high	Knee high	Thigh high (Grip top)	Knee high	Thigh high (Grip top)	
Size	Colour			Code			
Small	Black	PTS640040	PTS640530	PTS640601	PTS640830	PTS640605	
Medium	Black	PTS640050	PTS640540	PTS640602	PTS640840	PTS640606	
Large	Black	PTS640060	PTS640550	PTS640603	PTS640850	PTS640607	
X Large	Black	PTS640070	PTS640560	PTS640604	PTS640860	PTS640608	





DVA: AR22



457





PRODUCT AND SERVICES



Jobst Relief - Unisex - Knee High

- A more durable range, with compression values from 15 20mmHg to 30 40mmHg
- Available in open or closed toe

						U DVA	: AR22
	15 - 20mmH		mmHg	mHg 20 - 30mmHg			mmHg
Size	Colour	Open Toe	Closed Toe	Open Toe	Closed Toe	Open Toe	Closed Toe
Small	Beige	PTS642516	PTS642522	PTS641150	PTS641110	PTS641190	-
Small	Black	-	-	-	-	-	-
Medium	Beige	PTS642517	PTS642523	PTS641160	PTS641120	PTS641200	PTS641240
weatum	Black	-	-	-	-	-	-
Lorgo	Beige	PTS642518	PTS642524	PTS641170	PTS641130	PTS641210	PTS641250
Large	Black	-	-	-	-	-	-
VIana	Beige	PTS642519	PTS642525	PTS641180	PTS641140	PTS641220	PTS641260
X Large	Black	-	-	-	-	-	-
Large	Beige	PTS642520	PTS642526	PTS641181	PTS641141	PTS641221	PTS641261
(full calf)	Black	-	-	-	-	-	-
X Large	Beige	PTS642521	PTS642527	PTS641182	PTS641142	PTS641222	PTS641262
(full calf)	Black	-	-	-	-	-	-

Jobst Relief - Unisex - Thigh High (Grip Top)

- A more durable range, with compression values from 20 30mmHg to 30 40mmHg
- Available in open or closed toe

			DVA:	AR22	
		20 - 30	20 - 30mmHg		lmmHg
Size	Colour	Open Toe	Closed Toe	Open Toe	Closed Toe
Small	Beige	PTS641300	PTS641265	PTS641940	PTS641944
Medium	Beige	PTS641270	PTS641310	PTS641941	PTS641945
Large	Beige	PTS641280	PTS641319	PTS641942	PTS641946
X Large	Beige	PTS641290	PTS641320	PTS641943	PTS642264

Jobst Relief - Unisex - Chaps

• A more durable range, with compression values from 20 - 30mmHg to 30 - 40mmHg

				DVA:	AR22		
		OPEN TOE ONLY					
		20 - 30	mmHg	30 - 40	mmHg		
Size	Colour	Left	Right	Left	Right		
Small	Beige	PTS641349	PTS641359	PTS641371	PTS641374		
Medium	Beige	PTS641351	PTS641360	PTS641380	PTS641375		
Large	Beige	PTS641352	PTS641365	PTS641372	PTS641376		
X Large	Beige	PTS641350	PTS641370	PTS641373	PTS641377		



Jobst Relief - Unisex - Pantyhose

- A more durable range, with compression values from 20 30mmHg to 30 40mmHg
- Available in open or closed toe

/ Wanac		DVA: AR22			
		20 - 30	mmHg	30 - 40	mmHg
Size	Colour	Open Toe	Closed Toe	Open Toe	Closed Toe
Small	Beige	PTS641339	PTS641938	PTS641947	-
Medium	Beige	PTS641340	PTS641330	PTS641948	-
Large	Beige	PTS641341	PTS641939	PTS641949	-
X Large	Beige	PTS641342	PTS641335	PTS641950	-





Jobst 'UlcerCare' (40mmHg)

- Designed to help in the healing of venous leg ulcers, the Jobst UlcerCare two-stocking system provides the gradient • compression required for a venous ulcer to heal
- Starter Kit includes: stocking and 2 liners. Replacement Liners come in a 3 pack
- Important: Ensure the zipper is on the opposite side to the wound

			1.1		
	Size Guide		Starter Kit (Left Zip)	Starter Kit (Right Zip)	Replacement Liners
Size	Ankle (cm)	Calf (cm)		Code	
Small	18 - 21	29 - 35	PTS641951	PTS640990	PTS641025
Medium	21 - 24	34 - 41	PTS641952	PTS641000	PTS641030
Large	24 - 27	40 - 47	PTS641953	PTS641010	PTS641040
X Large	27 - 30	46 - 54	PTS641954	PTS641020	PTS641031
XX Large	30 - 33	48 - 56	PTS641955	PTS641021	PTS641032

Measurement: All measurements are circumferences.

DVA: AR22

Jobst Travel Socks - Unisex

- These socks help reduce the health risk associated with long distance travel
- They prevent swollen ankles and feet and reduce the risk of traveller's thrombosis
- Circulation is improved and legs feel great during and after travelling

Size	Colour	Code
1	Beige	PTS640026
I	Black	PTS640021
2	Beige	PTS640027
2	Black	PTS640022
3	Beige	PTS640028
3	Black	PTS640023
4	Beige	PTS640029
4	Black	PTS640024
5	Beige	PTS640030
5	Black	PTS640025



1

	CIRCUMFERENCE (cm)				
Size	l I	II	III	IV	V
Calf	27 - 37	30 - 40	33 - 43	36 - 46	39 - 49
Ankle	18 - 21	20 - 23	22 - 25	24 - 27	26 - 29

Mediven - Medical Legwear

Overview

Mediven offer a variety of compression stockings in a choice of colours and styles to suit a wide range of users

To Order

Choose your preferred model, then use the Size Guide to select the appropriate size and code

Size Guide

Take circumferential measurements to determine size (thigh measurement only needed for thigh and waist high stockings)

	MEDIVEN ELEGANCE - MEDIVEN COMFORT - MEDIVEN PLUS (UNISEX)						
		CIRCUMFERENCE (cm)					
Size	l I	II	III	IV	V	VI	VII
Ankle	18 - 20	20 - 22	22 - 24	24 - 26	26 - 28	28 - 30	30 - 32
Calf	28 - 34	30 - 37	33 - 40	35 - 43	37 - 46	39 - 49	41 - 51
Below knee	27 - 33	29 - 36	32 - 39	34 - 42	36 - 45	38 - 48	40 - 50
Thigh - plain top	43 - 57	45 - 62	49 - 67	53 - 72	56 - 77	60 - 81	64 - 85
Thigh - grip top	43 - 48	45 - 52	49 - 56	53 - 60	56 - 64	60 - 68	64 - 72

Measurement: All measurements are circumferences. Ankle - Narrowest part (just above ankle bone). Calf - Widest part. Below Knee - 2cm below knee. Thigh - 5cm below groin Note: Thigh measurement relevant for thigh, pantyhose and chaps styles only

P 1300 888 052



Ι



med



Mediven - Elegance Series - Knee High

- Sheer, elegant, fashionable and versatile
- Suitable for the female wearer

- Available in short or regular sizes
- Available in **beige only**

		DVA: AR22			
	CLOSED TOE ONLY				
	20 - 30mmHg	30 - 40mmHg			
Size	Short (34	4 - 38cm)			
1	PTS642001	PTS642003			
Ш	PTS642001BEII	PTS642003BEII			
Ш	PTS642001BEIII	PTS642003BEIII			
IV	PTS642001BEIV	PTS642003BEIV			
V	PTS642001BEV	PTS642003BEV			
VI	PTS642001BEVI	PTS642003BEVI			
VII	PTS642001BEVII	PTS642003BEVII			
Size	Regular (3	39 - 44cm)			
1	PTS642000	PTS642002			
Ш	PTS642000BEII	PTS642002BEII			
Ш	PTS642000BEIII	PTS642002BEIII			
IV	PTS642000BEIV	PTS642002BEIV			
V	PTS642000BEV	PTS642002BEV			
VI	PTS642000BEVI	PTS642002BEVI			
VII	PTS642000BEVII	PTS642002BEVII			



Measurement: To determine Short or Regular, measure from the floor to 2cm below the knee

Mediven - Elegance Series - Thigh High

- Sheer, elegant, fashionable and versatile
- Suitable for the female wearer

- Available in short or regular sizes
- Available in **beige only**

DVA: AR22

		TOE ONLY					
	20 - 30	mmHg	30 - 40	30 - 40mmHg			
Size		Short (62 - 71cm)					
5120	Plain Top	Grip Top	Plain Top	Grip Top			
1	PTS642005	PTS642009	PTS642007	PTS642011			
Ш	PTS642005BEII	PTS642009BEII	PTS642007BEII	PTS642011BEII			
Ш	PTS642005BEIII	PTS642009BEIII	PTS642007BEIII	PTS642011BEIII			
IV	PTS642005BEIV	PTS642009BEIV	PTS642007BEIV	PTS642011BEIV			
V	PTS642005BEV	PTS642009BEV	PTS642007BEV	PTS642011BEV			
VI	PTS642005BEVI	PTS642009BEVI	PTS642007BEVI	PTS642011BEVI			
VII	PTS642005BEVII	PTS642009BEVII	PTS642007BEVII	PTS642011BEVII			
Size	Regular (72 - 83cm)						
5120	Plain Top	Grip Top	Plain Top	Grip Top			
1	PTS642004	PTS642008	PTS642006	PTS642010			
Ш	PTS642004BEII	PTS642008BEII	PTS642006BEII	PTS642010BEII			
Ш	PTS642004BEIII	PTS642008BEIII	PTS642006BEIII	PTS642010BEIII			
IV	PTS642004BEIV	PTS642008BEIV	PTS642006BEIV	PTS642010BEIV			
V	PTS642004BEV	PTS642008BEV	PTS642006BEV	PTS642010BEV			
VI	PTS642004BEVI	PTS642008BEVI	PTS642006BEVI	PTS642010BEVI			
VII	PTS642004BEVII	PTS642008BEVII	PTS642006BEVII	PTS642010BEVII			



Measurement: To determine Short or Regular, measure from the floor to 5cm below the groin

General Enquiries



Mediven - Elegance Series - Pantyhose

- Sheer, elegant, fashionable and versatile
- Suitable for the female wearer

- Available in short or regular sizes
- Available in beige only

		DVA: AR22				
	CLOSED TOE ONLY					
	20 - 30mmHg	30 - 40mmHg				
Size	Short (62	2 - 71cm)				
1	PTS642013	PTS642015				
П	PTS642013II	PTS642015II				
Ш	PTS642013III	PTS642015III				
IV	PTS642013IV	PTS642015IV				
V	PTS642013V	PTS642015V				
VI	PTS642013VI	PTS642015VI				
VII	PTS642013VII	PTS642015VII				
Size	Regular (72 - 83cm)				
1	PTS642012	PTS642014				
П	PTS642012BEII	PTS642014II				
Ш	PTS642012BEIII	PTS642014III				
IV	PTS642012BEIV	PTS642014IV				
V	PTS642012BEV	PTS642014V				
VI	PTS642012BEVI	PTS642014VI				
VII	PTS642012BEVII	PTS642014VII				



Measurement: To determine Short or Regular, measure from the floor to 5cm below the groin

Mediven - Comfort Series - Knee High

- Highly elasticised opaque knit comfortable to wear
- Available in **beige only**

- Suitable for both males and females
- Two compression options

			\checkmark	DVA: AR22			
	20 - 30	mmHg	30 - 40mmHg				
Size	Short (34 - 38cm)						
Size	Open Toe	Closed Toe	Open Toe	Closed Toe			
1	PTS642019	PTS642017	PTS642035	PTS642033			
Ш	PTS642019BEII	PTS642017BEII	PTS642035BEII	PTS642033BEII			
Ш	PTS642019BEIII	PTS642017BEIII	PTS642035BEIII	PTS642033BEIII			
IV	PTS642019BEIV	PTS642017BEIV	PTS642035BEIV	PTS642033BEIV			
V	PTS642019BEV	PTS642017BEV	PTS642035BEV	PTS642033BEV			
VI	PTS642019BEVI	PTS642017BEVI	PTS642035BEVI	PTS642033BEVI			
VII	PTS642019BEVII	PTS642017BEVII	PTS642035BEVII	PTS642033BEVII			
Size	Regular (39 - 44cm)						
5126	Open Toe	Closed Toe	Open Toe	Closed Toe			
I	PTS642018	PTS642016	PTS642034	PTS642032			
Ш	PTS642018BEII	PTS642016BEII	PTS642034BEII	PTS642032BEII			
Ш	PTS642018BEIII	PTS642016BEIII	PTS642034BEIII	PTS642032BEIII			
IV	PTS642018BEIV	PTS642016BEIV	PTS642034BEIV	PTS642032BEIV			
V	PTS642018BEV	PTS642016BEV	PTS642034BEV	PTS642032BEV			
VI	PTS642018BEVI	PTS642016BEVI	PTS642034BEVI	PTS642032BEVI			
VII	PTS642018BEVII	PTS642016BEVII	PTS642034BEVII	PTS642032BEVII			

Measurement: To determine Short or Regular, measure from the floor to 2cm below the knee





461



Mediven - Comfort Series - Thigh High (Grip Top)

- Highly elasticised opaque knit comfortable to wear •
- Available in **beige only** •

- Suitable for both males and females
- Two compression options

			\checkmark	DVA: AR22		
	20 - 30	mmHg	30 - 40	mmHg		
Size	Short (62 - 71cm)					
5126	Open Toe	Closed Toe	Open Toe	Closed Toe		
1	PTS642023	PTS642021	PTS642039	PTS642037		
Ш	PTS642023BEII	PTS642021BEII	PTS642039BEII	PTS642037BEII		
III	PTS642023BEIII	PTS642021BEIII	PTS642039BEIII	PTS642037BEIII		
IV	PTS642023BEIV	PTS642021BEIV	PTS642039BEIV	PTS642037BEIV		
V	PTS642023BEV	PTS642021BEV	PTS642039BEV	PTS642037BEV		
VI	PTS642023BEVI	PTS642021BEVI	PTS642039BEVI	PTS642037BEVI		
VII	PTS642023BEVII	PTS642021BEVII	PTS642039BEVII	PTS642037BEVII		
Size	Regular (72 - 83cm)					
5120	Open Toe	Closed Toe	Open Toe	Closed Toe		
1	PTS642022	PTS642020	PTS642038	PTS642036		
II	PTS642022BEII	PTS642020BEII	PTS642038BEII	PTS642036BEII		
Ш	PTS642022BEIII	PTS642020BEIII	PTS642038BEIII	PTS642036BEIII		
IV	PTS642022BEIV	PTS642020BEIV	PTS642038BEIV	PTS642036BEIV		
V	PTS642022BEV	PTS642020BEV	PTS642038BEV	PTS642036BEV		
VI	PTS642022BEVI	PTS642020BEVI	PTS642038BEVI	PTS642036BEVI		
VII	PTS642022BEVII	PTS642020BEVII	PTS642038BEVII	PTS642036BEVII		



Measurement: To determine Short or Regular, measure from the floor to 5cm below the groin

Mediven - Comfort Series - Pantyhose

- Highly elasticised opaque knit comfortable to wear •
- Available in beige only

General Enquiries

- Suitable for both males and females
- Two compression options

			\checkmark	DVA: AR22		
	20 - 30	mmHg	30 - 40	30 - 40mmHg		
Size	Short (62 - 71cm)					
	Open Toe	Closed Toe	Open Toe	Closed Toe		
1	PTS642031	PTS642029	PTS642047	PTS642045		
Ш	PTS642031BEII	PTS642029BEII	PTS642047BEII	PTS642045BEII		
Ш	PTS642031BEIII	PTS642029BEIII	PTS642047BEIII	PTS642045BEIII		
IV	PTS642031BEIV	PTS642029BEIV	PTS642047BEIV	PTS642045BEIV		
V	PTS642031BEV	PTS642029BEV	PTS642047BEV	PTS642045BEV		
VI	PTS642031BEVI	PTS642029BEVI	PTS642047BEVI	PTS642045BEVI		
VII	PTS642031BEVII	PTS642029BEVII	PTS642047BEVII	PTS642045BEVII		
Size	Regular (72 - 83cm)					
5126	Open Toe	Closed Toe	Open Toe	Closed Toe		
1	PTS642030	PTS642028	PTS642046	PTS642044		
Ш	PTS642030BEII	PTS642028BEII	PTS642046BEII	PTS642044BEII		
Ш	PTS642030BEIII	PTS642028BEIII	PTS642046BEIII	PTS642044BEIII		
IV	PTS642030BEIV	PTS642028BEIV	PTS642046BEIV	PTS642044BEIV		
V	PTS642030BEV	PTS642028BEV	PTS642046BEV	PTS642044BEV		
VI	PTS642030BEVI	PTS642028BEVI	PTS642046BEVI	PTS642044BEVI		
VII	PTS642030BEVII	PTS642028BEVII	PTS642046BEVII	PTS642044BEVII		



Measurement: To determine Short or Regular, measure from the floor to 5cm below the groin



Mediven - Plus Series - Knee High

- Soft but durable with antibacterial action for fresher feet
- Available in **beige only**

- Suitable for both males and females
- Three compression options available

			JVA: ARZZ
		OPEN TOE ONLY	
	20 - 30mmHg	30 - 40mmHg	40 - 50mmHg
Size		Short (34 - 38cm)	
1	PTS642051	PTS642052	PTS642053
II	PTS642051BEII	PTS642052BEII	PTS642053BEII
III	PTS642051BEIII	PTS642052BEIII	PTS642053BEIII
IV	PTS642051BEIV	PTS642052BEIV	PTS642053BEIV
V	PTS642051BEV	PTS642052BEV	PTS642053BEV
VI	PTS642051BEVI	PTS642052BEVI	PTS642053BEVI
VII	PTS642051BEVII	PTS642052BEVII	PTS642053BEVII
Size		Regular (39 - 44cm)	
1	PTS642048	PTS642049	PTS642050
II	PTS642048BEII	PTS642049BEII	PTS642050BEII
III	PTS642048BEIII	PTS642049BEIII	PTS642050BEIII
IV	PTS642048BEIV	PTS642049BEIV	PTS642050BEIV
V	PTS642048BEV	PTS642049BEV	PTS642050BEV
VI	PTS642048BEVI	PTS642049BEVI	PTS642050BEVI
VII	PTS642048BEVII	PTS642049BEVII	PTS642050BEVII



Measurement: To determine Short or Regular, measure from the floor to 2cm below the knee

Mediven - Plus Series - Thigh High

- Soft but durable with antibacterial action for fresher feet •
- Available in beige only •
- Choice of Plain Top or Grip Top
- Three compression options available •



	OPEN TOE ONLY						
	20 - 30	mmHg	30 - 40	ImmHg	40 - 50	40 - 50mmHg	
Size	Short (62 - 71cm)						
0126	Plain Top	Grip Top	Plain Top	Grip Top	Plain Top	Grip Top	
1	PTS642059	PTS642065	PTS642060	PTS642066	PTS642061	PTS642067	
11	PTS642059BEII	PTS642065BEII	PTS642060BEII	PTS642066BEII	PTS642061BEII	PTS642067BEII	
111	PTS642059BEIII	PTS642065BEIII	PTS642060BEIII	PTS642066BEIII	PTS642061BEIII	PTS642067BEIII	
IV	PTS642059BEIV	PTS642065BEIV	PTS642060BEIV	PTS642066BEIV	PTS642061BEIV	PTS642067BEIV	
V	PTS642059BEV	PTS642065BEV	PTS642060BEV	PTS642066BEV	PTS642061BEV	PTS642067BEV	
VI	PTS642059BEVI	PTS642065BEVI	PTS642060BEVI	PTS642066BEVI	PTS642061BEVI	PTS642067BEVI	
VII	PTS642059BEVII	PTS642065BEVII	PTS642060BEVII	PTS642066BEVII	PTS642061BEVII	PTS642067BEVII	
Size	Regular (72 - 83cm)						
Size	Plain Top	Grip Top	Plain Top	Grip Top	Plain Top	Grip Top	
1	PTS642056	PTS642062	PTS642057	PTS642063	PTS642058	PTS642064	
11	PTS642056BEII	PTS642062BEII	PTS642057BEII	PTS642063BEII	PTS642058BEII	PTS642064BEII	
111	PTS642056BEIII	PTS642062BEIII	PTS642057BEIII	PTS642063BEIII	PTS642058BEIII	PTS642064BEIII	
IV	PTS642056BEIV	PTS642062BEIV	PTS642057BEIV	PTS642063BEIV	PTS642058BEIV	PTS642064BEIV	
V	PTS642056BEV	PTS642062BEV	PTS642057BEV	PTS642063BEV	PTS642058BEV	PTS642064BEV	
VI	PTS642056BEVI	PTS642062BEVI	PTS642057BEVI	PTS642063BEVI	PTS642058BEVI	PTS642064BEVI	
VII	PTS642056BEVII	PTS642062BEVII	PTS642057BEVII	PTS642063BEVII	PTS642058BEVII	PTS642064BEVII	

Measurement: To determine Short or Regular, measure from the floor to 5cm below the groin

DVA Contact Details

I For NDIS Item Code, refer to page: 500 ndis



463



Mediven - Plus Series - Chaps

- Soft but durable with antibacterial action for fresher feet •
- Available in beige only .

- Choose between left or right for adequate fit
- Three compression options available •

						OVA: AR22	
			OPEN T	OE ONLY			
	20 - 30mmHg 30 - 40mmHg			40 - 50	ImmHg		
Size			Short (6	2 - 71cm)	2 - 71cm)		
0126	Left	Right	Left	Right	Left	Right	
1	PTS642071	PTS642077	PTS642072	PTS642078	PTS642073	PTS642079	
П	PTS642071BEII	PTS642077BEII	PTS642072BEII	PTS642078BEII	PTS642073BEII	PTS642079BEII	
Ш	PTS642071BEIII	PTS642077BEIII	PTS642072BEIII	PTS642078BEIII	PTS642073BEIII	PTS642079BEIII	
IV	PTS642071BEIV	PTS642077BEIV	PTS642072BEIV	PTS642078BEIV	PTS642073BEIV	PTS642079BEIV	
V	PTS642071BEV	PTS642077BEV	PTS642072BEV	PTS642078BEV	PTS642073BEV	PTS642079BEV	
VI	PTS642071BEVI	PTS642077BEVI	PTS642072BEVI	PTS642078BEVI	PTS642073BEVI	PTS642079BEVI	
VII	PTS642071BEVII	PTS642077BEVII	PTS642072BEVII	PTS642078BEVII	PTS642073BEVII	PTS642079BEVII	
Size			Regular (lar (72 - 83cm)			
3126	Left	Right	Left	Right	Left	Right	
1	PTS642068	PTS642074	PTS642069	PTS642075	PTS642070	PTS642076	
П	PTS642068BEII	PTS642074BEII	PTS642069BEII	PTS642075BEII	PTS642070BEII	PTS642076BEII	
Ш	PTS642068BEIII	PTS642074BEIII	PTS642069BEIII	PTS642075BEIII	PTS642070BEIII	PTS642076BEIII	
IV	PTS642068BEIV	PTS642074BEIV	PTS642069BEIV	PTS642075BEIV	PTS642070BEIV	PTS642076BEIV	
V	PTS642068BEV	PTS642074BEV	PTS642069BEV	PTS642075BEV	PTS642070BEV	PTS642076BEV	
VI	PTS642068BEVI	PTS642074BEVI	PTS642069BEVI	PTS642075BEVI	PTS642070BEVI	PTS642076BEVI	
VII	PTS642068BEVII	PTS642074BEVII	PTS642069BEVII	PTS642075BEVII	PTS642070BEVII	PTS642076BEVII	

Measurement: To determine Short or Regular, measure from the floor to 5cm below the groin

Mediven - Plus Series - Pantyhose

- Soft but durable with antibacterial action for fresher feet •
- Available in **beige only**
- Three compression options available

Thre	e compression options avai	lable 🕜	DVA: AR22
		OPEN TOE ONLY	
	20 - 30mmHg	30 - 40mmHg	40 - 50mmHg
Size		Short (62 - 71cm)	
I	PTS642083	PTS642084	PTS642085
II	PTS642083BEII	PTS642084BEII	PTS642085BEII
III	PTS642083BEIII	PTS642084BEIII	PTS642085BEIII
IV	PTS642083BEIV	PTS642084BEIV	PTS642085BEIV
V	PTS642083BEV	PTS642084BEV	PTS642085BEV
VI	PTS642083BEVI	PTS642084BEVI	PTS642085BEV
VII	PTS642083BEVII	PTS642084BEVII	PTS642085BEVI
Size		Regular (72 - 83cm)	
I	PTS642080	PTS642081	PTS642082
11	PTS642080BEII	PTS642081BEII	PTS642082BEII
111	PTS642080BEIII	PTS642081BEIII	PTS642082BEIII
IV	PTS642080BEIV	PTS642081BEIV	PTS642082BEIV
V	PTS642080BEV	PTS642081BEV	PTS642082BEV
VI	PTS642080BEVI	PTS642081BEVI	PTS642082BEVI
VII	PTS642080BEVII	PTS642081BEVII	PTS642082BEVI
Measurem	ent: To determine Short or Regular, me	asure from the floor to 5cm below the	groin



Sigvaris

To Order

Choose your preferred model. Then, order using: Code, Description, Styles, Compression (mmHg), Size

Sigvaris - Cotton - Unisex

- Particularly suitable for users with sensitive skin
- Colour available: Natural

The mix of fibres also ensures that the stocking slides on and off easily

Size Guide							
Size	Ankle (cm)	Below Knee (cm) (Calf)	Upper Thigh (cm) (Thigh Grip Tops)	Upper Thigh (cm) (Pantyhose & Chaps)			
XS	20 - 22	27 - 34	47 - 56	47 - 56			
S	22 - 24	31 - 38	52 - 61	52 - 64			
Μ	24 - 26	34 - 41	56 - 65	56 - 68			
L	26 - 29	37 - 44	61 - 71	61 - 74			



DVA: AR22

Measurement: All measurements are circumferences. Ankle - Narrowest part (just above ankle bone). Below Knee - 2cm below knee. Upper Thigh - 5cm below groin

		Compression (mmHg)				
Code	Description	Style	18 - 21 (Class 1)	23 - 32 (Class 2)	34 - 46 (Class 3)	Size
PTS642310	Cotton AD Calf	Open / Closed Toe	\checkmark	\checkmark	\checkmark	XS / S / M / L
PTS642315	Cotton Extra AD Calf	Open / Closed Toe	\checkmark	\checkmark	\checkmark	XS / S / M / L
PTS642320	Cotton Calf with Grip Top	Open / Closed Toe	-	\checkmark	-	XS / S / M / L
PTS642325	Cotton Thigh with Grip Top	Open / Closed Toe	\checkmark	\checkmark	-	XS / S / M / L
PTS642330	Cotton Pantyhose	Open / Closed Toe	-	\checkmark	-	XS / S / M / L
PTS642335	Cotton Chaps	Left or Right	\checkmark	\checkmark	\checkmark	XS / S / M / L

Sigvaris - Traditional - Unisex

- Reliable safe and effective option
- Colour available: Beige

The elasticity and breathability mean that they are very comfortable to wear

Size Guide							
Size	Ankle (cm)	Below Knee (cm) (Calf)	Mid Thigh (cm) (Thigh)	Upper Thigh (cm) (Pantyhose & Chaps)			
XS	17 - 19	26 - 32	39 - 48	40 - 50			
S	19 - 23	28 - 34	41 - 50	42 - 52			
Μ	23 - 26	31 - 37	48 - 57	50 - 60			
L	26 - 29	34 - 40	55 - 64	58 - 68			

Measurement: All measurements are circumferences. Ankle - Narrowest part (just above ankle bone). Below Knee - 2cm below knee. Mid Thigh - Around mid thigh. Upper Thigh - 5cm below groin



DVA Contact Details



I



Pressure Care & Circulation

465



PRODUCT AND SERVICES CATALOGUE



Sigvaris - Magic - Women

- The elegant look and subtle sheer quality of Sigvaris Magic appeals to the senses •
- This versatile compression stocking is exceptionally easy to handle and offers a perfect fit •
- Colour available: Skin

Size Guide						
Size	Ankle (cm)	Below Knee (cm) (Calf)	Upper Thigh (cm) (Thigh Grip Tops)	Upper Thigh (cm) (Pantyhose & Chaps)		
XS	20 - 22	27 - 34	47 - 56	47 - 56		
S	22 - 24	31 - 38	52 - 61	52 - 64		
Μ	24 - 26	34 - 41	56 - 65	56 - 68		
L	26 - 29	37 - 44	61 - 71	61 - 74		



Measurement: All measurements are circumferences. Ankle - Narrowest part (just above ankle bone). Below Knee - 2cm below knee. Upper Thigh - 5cm below groin

					DVA: AR22
			Compression (mmHg)		
Code	Description	Style	18 - 21 (Class 1)	23 - 32 (Class 2)	Size
PTS642380	Magic Calf	Open / Closed Toe	\checkmark	\checkmark	XS / S / M / L
PTS642385	Magic Thigh with Grip Top	Open / Closed Toe	\checkmark	\checkmark	XS / S / M / L
PTS642390	Magic Panty	-	\checkmark	\checkmark	XS / S / M / L

Sigvaris - Ulcer X - Unisex

- Ulcer X is a compression stocking that is easy to don, comfortable and discrete
- Offers quicker healing rates compared to conventional bandaging
- Colour available: Beige •

Size Guide						
Size	Ankle (cm)	Below Knee (cm) (Calf)				
S	19 - 23	28 - 34				
Μ	23 - 26	31 - 37				
L	26 - 29	34 - 40				

Measurement: All measurements are circumferences. Ankle - Narrowest part (just above ankle bone) Below Knee - 2cm below knee



S/M/L

One Size

Compression (mmHg)

(Class 2)

1

Ulcer X Kit - Calf Open Toe PTS642395 (1 stocking; 2 underliners) PTS642400 Ulcer X Underliners

Enquiries



General



Sigvaris - Post Surgery - Unisex

- This stocking is suitable for patients who have undergone leg surgery
- Available in traditional quality for maximum pressure stability throughout the day
- Colour available: **Beige**

Size Guide							
Size	Ankle (cm)	Below Knee (cm) (Calf)	Mid Thigh (cm) (Thigh)	Upper Thigh (cm) (Pantyhose & Chaps)			
S	19 - 23	28 - 34	41 - 50	42 - 52			
Μ	23 - 26	31 - 37	48 - 57	50 - 60			
L	26 - 29	34 - 40	55 - 64	58 - 68			

Measurement: All measurements are circumferences. Ankle - Narrowest part (just above ankle bone). Below Knee - 2cm below knee Mid Thigh - Around mid thigh. Upper Thigh - 5cm below groin

		Style	Compressi	Ci-r	
Code Description	23 - 32 (Class 2)		34 - 46 (Class 3)	Size	
PTS642405	Post Surgery Kit (1 stocking; 2 underliners)	Open Toe	\checkmark	\checkmark	S / M / L
PTS642410	Post Surgery Underliners	-	-	-	One Size

Solaris

SOLARIS

DVA: AR22

Overview

- The Solaris brand offers a wide range of compression garments to suit users with different needs
- The range includes ReadyWrap, ExoStrong, ExoFusion, Caresia

To Order

Choose your preferred model, then use the Size Guide to select the appropriate size and code

ReadyWrap (Inelastic Compression)

- ReadyWrap garments are designed to aid people living with venous disease or insufficiencies
- Created to replace daytime compression bandaging or serve as an alternative or supplement to elastic compression stockings
- Silver microfibre under-garment liner is available for added comfort
- Contraindications: arterial insufficiency or degeneration, deep vein thrombosis, untreated congestive heart failure, localised or systemic infection or cancer, absent or severely impaired sensation

Solaris - ReadyWrap - Toe

- Unique design allows for adjustable toe girths
- Provides full foot coverage when paired with a ReadyWrap foot unit

				A: AR22	
	Size Guide		Left	Right	
Size	MTP Girth (cm)	Colour	Code		
Medium	20.5 - 23.5	Beige	PTS642645	PTS642648	
Large	23.6 - 26	Beige	PTS642644	PTS642647	
X Large	26.1 - 28.5	Beige	PTS642646	PTS642649	







For NDIS Item Code, refer to page: 500



Compression - Legwear Pressure Care & Circulation



Solaris - ReadyWrap - Foot

- The simple two-strap design of the ReadyWrap foot unit fits into most shoes and allows for easier donning
- It is the perfect alternative for those having difficulty donning elastic stockings or whose skin may be compromised by shear forces
- Left and right specific design provides enhanced foot contouring



						DVA: AR22
Size Guide					Left	Right
Size	e Foot Length (cm) Y Girth (cm) Ankle Girth (cm) Colour			Co	de	
o "	un to 01	28 - 33	18 - 23	Beige	FBS303597	FBS303605
Small	up to 21	28 - 33		Black	FBS303598	FBS303606
N d a all'a sura		20. 20	20 - 28	Beige	FBS303595	FBS303603
Medium	up to 23	30 - 38		Black	FBS303596	FBS303604
1		05 40	05 00	Beige	FBS303593	FBS303601
Large	up to 26	35 - 43	25 - 33	Black	FBS303594	FBS303602
V Lawrence		40 40	20, 20,	Beige	FBS303599	FBS303607
X Large up to 29 40	40 - 48	30 - 38	Black	FBS303600	FBS303608	

Solaris - ReadyWrap - Calf

- Maximum therapeutic benefits, excellent contouring and complete comfort
- The strapping system mimics the 50% overlap of short-stretch bandaging and reduces shear force pressure when being applied to the limb
- Velcro strap placement is designed to minimise strap shifting and enhance ease of donning
- Low-stretch materials provide simultaneous low resting/high working pressure, maximising calf muscle function, especially when exercising or walking
- Can be worn alone or on top of compression hosiery



					DVA. AII22
	Size Guide			Average (length up to 30cm)	Tall (length up to 35cm)
Size	Ankle Girth (cm)	Calf Girth (cm)	Colour	Co	de
Small	18 - 23	28 - 38	Beige	FBS303617	FBS303619
SIIIdii	10 - 23	28 - 38	Black	FBS303618	FBS303620
N.4 a alfavora	00 00	28 33 - 43	Beige	FBS303613	FBS303615
Medium	20 - 28		Black	FBS303614	FBS303616
Lours	05 00	38 - 48	Beige	FBS303609	FBS303611
Large	25 - 33	38 - 48	Black	FBS303610	FBS303612
X I away	20. 20	42 50	Beige	FBS303621	FBS303623
X Large	30 - 38	43 - 58	Black	FBS303622	FBS303624
XX Lawrence		52 62	Beige	FBS303625	FBS303627
XX Large	35 - 43	53 - 68	Black	FBS303626	FBS303628

General Enquiries



Solaris - ReadyWrap - Knee

- The ReadyWrap knee unit provides support for those with localised patella swelling or acts as a bridge connecting the thigh and calf units to provide full leg coverage
- The simple 3-strap design makes that unit easy to don and the 50% overlap mimics short-stretch bandaging

				DVA: AR22	
	Size	Guide			
Size	Calf Girth (cm)	Above Knee Girth (cm)	Colour	Code	
	20 20	38 - 48	Beige	FBS303633	
Small	28 - 38		Black	FBS303634	
Medium	33 - 43	45 - 58	Beige	FBS303631	
weatum	33 - 43		Black	FBS303632	
Large	38 - 48 52 -	52 69	Beige	FBS303629	
		52 - 68	Black	FBS303630	
X Large	43 - 58	40 50 50 70	50 70	Beige	FBS303635
		59 - 78	Black	FBS303636	
	53 - 68		<u> </u>	Beige	FBS303637
XX Large		68 66 -88	Black	FBS303638	



Solaris - ReadyWrap - Thigh

- The perfect addition to the ReadyWrap knee, calf and foot units for those with full leg oedema
- When used on its own, it is also a great option to prevent localised swelling from the superior patella border to the groin
- The strapping system mimics the 50% overlap of short-stretch bandaging and reduces shear force pressure when being applied to the limb

					DVA: AR22
	Size	Guide		Average (length up to 25cm)	Tall (length up to 30cm)
Size	Above Knee Girth (cm)	Groin Girth (cm)	Colour	Code	
Small	38 - 48	48 - 63	Beige	FBS303647	FBS303649
Small 38 - 48	30 - 40	40 - 03	Black	FBS303648	FBS303650
Medium	45 - 58	58 - 73	Beige	FBS303643	FBS303645
Wealum 45 - 58	58 - 73	Black	FBS303644	FBS303646	
Large	52 - 68	52 - 68 68 - 83	Beige	FBS303639	FBS303641
_arge 52 - 68	5 00-03	Black	FBS303640	FBS303642	
(Largo	50 70	78 - 93	Beige	FBS303651	FBS303653
X Large 59 - 78	78 - 93	Black	FBS303652	FBS303654	
	00 100	Beige	FBS303655	FBS303657	
(X Large 66 - 88		88 88 - 103	Black	FBS303656	FBS303658

Solaris - Extender Strap

- Provides an additional 10cm to a standard ReadyWrap strap
- The unique forked design ensures a secure fit •
- Available in black or beige

	DVA: AR22
Description	Code
Extender Strap - Beige	FBS303591
Extender Strap - Black	FBS303592

Solaris - Silver Liner

- These non-compressive Solaris Silver Liners are specifically designed to be used under the Solaris ReadyWrap for additional comfort
- Manufactured with unique silver, antimicrobial properties for maximum comfort and protection •
- Ordered per pair fits left or right
- Colour available: Black

		D W W WILL		
Description	Girth Range (cm)	Calf Girth (cm)	Code	
Below Knee	Up to 68	28 - 38	PTS642641	
Thigh High - Universal	48 - 74	33 - 43	PTS642643	
Thigh High - Large	74 - 103	38 - 48	PTS642642	



469









Solaris - ExoStrong

- ExoStrong garments are designed to move with you
- With a functional design at the ankle and knee, you can enjoy the freedom to move
- ExoStrong garments use flat seams. This helps reduce skin irritation and keeps you comfortable
- ExoStrong provides 20 30mmHg of compression
- An open toe design
- Available in **beige**

Solaris - ExoStrong - Knee (20 - 30mmHg)



Size Guide				Short Foot Length: up to 17cm Calf Length: up to 38cm	Average Foot Length: up to 19cm Calf Length: up to 41cm	Tall Foot Length: up to 22cm Calf Length: up to 44cm
Size	Ankle Girth (cm)	Calf Girth (cm)	Style	Code		
Small	18 - 20	27 - 36	Knit Top	-	PTS642595	PTS642605
Small	18 - 20		Silicone Top	PTS642587	PTS642596	PTS642606
Madium	ledium 20 - 24	29 - 29	Knit Top	PTS642588	PTS642597	PTS642607
weatum			Silicone Top	-	PTS642598	PTS642608
Laura	24 - 28 33	- 28 32 - 43	Knit Top	PTS642589	PTS642599	PTS642609
Large			Silicone Top	PTS642590	PTS642600	PTS642610
VIerre	00 01	33 - 45	Knit Top	PTS642591	PTS642601	PTS642611
X Large 28 - 31	33 - 45	Silicone Top	PTS642592	PTS642602	PTS642612	
VV Lorgo	22 24	07 50	Knit Top	PTS642593	PTS642603	PTS642613
XX Large 32 - 34		37 - 50	Silicone Top	PTS642594	PTS642604	PTS642614

Solaris - ExoStrong - Thigh (20 - 30mmHg)

Size Guide				Short Foot Length: up to 17cm Thigh Length: up to 71cm	Average Foot Length: up to 19cm Thigh Length: up to 74cm	Tall Foot Length: up to 22cm Thigh Length: up to 77cm
Size	Ankle Girth (cm)	Calf Girth (cm)	Thigh Girth (cm)	Code		
Small	18 - 20	27 - 36	41 - 55	PTS642615S	PTS642616S	PTS642617S
Medium	20 - 24	29 - 29	45 - 61	PTS642615M	PTS642616M	PTS642617M
Large	24 - 28	32 - 43	50 - 66	PTS642615L	PTS642616L	PTS642617L
X Large	28 - 31	33 - 45	56 - 71	PTS642615XL	PTS642616XL	PTS642617XL

Solaris - ExoFusion - Knee (20 - 30mmHg)

- Focused compression over the foot makes ExoFusion easier to put on than a regular compression stocking
- ExoFusion is made with material that helps regulate temperature and moisture
- ExoFusion provides a focused 20 30mmHg compression over the foot
- Available in **black**

	IdCK	DVA: AR22		
Size	Guide	Short (Foot Length: up to 21cm)	Tall (Foot Length: up to 28cm)	
Size	MTP Girth (cm)	Code		
Medium	20 - 24	PTS642581	PTS642584	
Large	24 - 28	PTS642582	PTS642585	
X Large	28 - 32	PTS642583	-	




 \bigcirc

Solaris - Caresia - Lymphoedema Compression Bandage Liners

- Simplify complex multi-layered bandaging protocols by replacing and improving upon foam liners, stockinette and gauze
- The unique foam chips used in every Caresia garment help resolve tissue induration, stimulate interstitial fluid movement, and diffuse irregular wrapping pressures. The foam also helps provide padding protection to limbs

Foot	DV/	DVA: AR22		
	Size Guide			
Size	Foot Length (cm)	Ankle Girth (cm)	Code	
Small	20 - 25	22 - 28	PTS6426560	
Medium	25 - 30	28 - 34	PTS6426561	
Large	30 - 35	34 - 40	PTS6426562	



DVA: AR22

Below Knee

Size Guide			Short (Calf Length: 35 - 41cm)	Average (Calf Length: 41 - 47cm)	Tall (Calf Length: 47 - 53cm)
Size	Ankle Girth (cm)	Calf Girth (cm)		Code	
Small	22 - 28	30 - 40	PTS642230	PTS642230AS	PTS642230TS
Medium	28 - 34	40 - 52	PTS642230M	PTS642230AM	PTS642230TM
Large	34 - 40	52 - 62	PTS642230L	PTS642230AL	PTS642230TL

Thigh

			C C	
Size Guide		Short Average (Thigh Length: 23 - 28cm) (Thigh Length: 28 - 33cm)		Tall (Thigh Length: 33 - 38cm)
Size	Groin Girth (cm)		Code	
Small	45 - 65	PTS642625	PTS642622	PTS642628
Medium	65 - 85	PTS642624	PTS642621	PTS642627
Large	85 - 99.9	PTS642623	PTS642620	PTS642626









TED Anti-Embolism Legwear

Overview

- TED anti-embolism stockings with graduated compression are designed to ensure optimum blood flow velocity when the wearer is immobile in bed for long periods of time
- They are widely used in hospitals for non-ambulatory or convalescing persons to avoid the danger of thromboembolic disease (blood clotting)
- TED stockings are full-foot, but feature an inspection opening under the toes
- Unisex and available in **white**

To Order

Choose your preferred model, then use the Size Guide to select the appropriate size and code

TED - Knee High

- TED knee high stockings are designed to reduce venous stasis below the knee when thigh length styles are medically contraindicated
- Knee length with soft top band to minimise slippage without constricting circulation
- Available in 6 sizes, with choice of either regular or long fit



	S					
Size	Calf (cm)	Leg Length (cm) (Regular)	Leg Length (cm) (Long)	Product Code (Regular)	Product Code (Long)	
Small	to 30.5	to 40.6	40.6+	PTS641570	PTS641600	
Medium	30.5 - 38.1	to 43.2	43.2+	PTS641580	PTS641610	
Large	38.1 - 44.5	to 45.7	45.7+	PTS641590	PTS641620	
X Large	44.5 - 50.8	to 45.7	45.7+	PTS641591	PTS641621	
XX Large	50.8 - 58.4	to 45.7	45.7+	PTS641592	PTS641622	
XXX Large	58.4 - 66.0	to 45.7	45.7+	PTS641593	PTS641623	

Measurement: All measurements are circumferences. **Calf** - Widest part. **Leg Length** - From heel to bend in back of knee. Note: if one measurement indicates a larger size than the other, select the larger size.

TED - Thigh High - Grip Top

- Pressure break at knee, protecting the critical popliteal area and promoting free leg movement in the recumbent patient
- Available in 5 sizes, with choice of either, short, regular or long fit

							0	
Size Guide								
Size Call	Calf (cm)	Thigh (are)		Leg Length (cm)		Code (Short)	Code (Regular)	Code (Long)
		Thigh (cm)	Short	Regular	Long	(0.1014)	(insgenar)	
Small	to 30.5	to 63.5	to 73.7	73.7 - 83.8	83.8+	PTS641471	PTS641476	PTS641482
Medium	30.5 - 38.1	to 63.5	to 73.7	73.7 - 83.8	83.8+	PTS641475	PTS641477	PTS641483
Large	38.1 - 44.5	to 63.5	to 73.7	73.7 - 83.8	83.8+	PTS641472	PTS641478	PTS641484
X Large	44.5 - 54.6	63.5 - 81.3	to 73.7	73.7 - 83.8	83.8+	PTS641473	PTS641479	PTS641485
XX Large	54.6 - 66.0	81.3 - 91.4	to 73.7	73.7 - 83.8	83.8+	PTS641474	PTS641481	PTS641486

Measurement: All measurements are circumferences. Calf - Widest part. Thigh - 5cm below groin. Leg Length - From heel to gluteal furrow (just below groin) Note: In thigh length stockings, when applied, the stitch change (change in fabric sheerness) should fall between 2.5 to 5cm below the bend in the knee. If thigh exceeds 91.4cm, select knee length style





TED - Thigh High - Waist Belt

- The TED thigh high with waist belt is very suitable for thin or obese patients, where stocking slippage is often a problem
- One-piece, snap-on belt is easily adjustable
- Pressure break at knee, protecting the critical popliteal area and promoting free leg movement in the recumbent patient
- Available in 5 sizes, with choice of either regular or long fit
- Supplied per pair



DVA: AR22

		Size Guide	Code	Code			
Size	Calf (cm)	Thigh (cm)	Leg Length (cm) (Regular)	Leg Length (cm) (Long)	(Regular)	(Long)	
X Small	to 25.4	to 63.5	to 71.1	71.1+	PTS641425	PTS641501	
Small	25.4 - 30.5	to 63.5	to 72.4	72.4+	PTS641430	PTS641502	
Medium	30.5 - 38.1	to 63.5	to 72.4	72.4+	PTS641460	PTS641503	
Large	38.1 - 44.5	to 63.5	to 73.7	73.7+	PTS641470	PTS641504	
X Large	38.1 - 44.5	63.5 - 81.3	to 72.4	72.4+	PTS641500	PTS641505	

Measurement: All measurements are circumferences. Calf - Widest part. Thigh - 5cm below groin. Leg Length - From heel to gluteal furrow (just below groin) Note: In thigh length stockings, when applied, the stitch change (change in fabric sheerness) should fall between 2.5 to 5cm below the bend in the knee.

TED - Continuing Care

- Designed for the continuing care of the recuperating patient after leaving hospital
- To be worn only on physician's advice and during the time period recommended
- Unisex, knee length, closed toe without inspection opening
- Available in 4 sizes, with choice of either regular or long fit
- Supplied per pair
- Colour: White

/	<u> </u>					
		Size Guide	Ocale	0-1-		
Size	Calf (cm)	Leg Length (cm) (Regular)	Leg Length (cm) (Long)	Code (Regular)	Code (Long)	
Small	to 30.5	to 40.6	40.6+	PTS641630	PTS641634	
Medium	30.5 - 38.1	to 43.2	43.2+	PTS641631	PTS641635	
Large	38.1 - 44.5	to 45.7	45.7+	PTS641632	PTS641636	
X Large	44.5 - 50.8	to 45.7	45.7+	PTS641633	PTS641637	

Measurement: All measurements are circumferences. Calf - Widest part. Leg Length - From heel to bend in back of knee.

Venosan - Medical Legwear

Overview

 Venosan compression garments have been manufactured since 1942 and are distributed in over 70 countries

To Order

Choose your preferred model, then use the Size Guide to select the appropriate size and code

Size Guide

• Take circumferential measurements to determine size (thigh measurement only needed for thigh & waist high stockings)

Size Guide									
Circumference	Small	Medium	Large	X Large					
Ankle (cm)	19 - 23	23 - 26	26 - 29	29 - 31					
Calf (cm)	29 - 39	33 - 42	36 - 45	39 - 47					
Below Knee (cm)	28 - 38	32 - 41	35 - 44	38 - 46					
Thigh (cm)	48 - 63	54 - 67	58 - 72	63 - 74					

Measurement: All measurements are circumferences. Ankle - Narrowest part (just above ankle bone). Calf - Widest part. Below Knee - 2cm below knee. Upper Thigh - 5cm below groin







Venosan 'Perlflor' Medical Legwear

- VENOSAN® Perlflor is a soft, strong and robust garment with a high cotton content
- Suitable for: lymphoedema, severe swelling, lipoedema & ulcer-care
- Available in beige and open toe only. Unisex sizes
- (√ Size Guide Code Ankle (cm) Calf (cm) Below Knee (cm) 20 - 22 31 - 36 PTS641999 29 - 34 22 - 24 34 - 39 32 - 37 PTS641999M 24 - 26 37 - 42 35 - 40 PTS641999L

38 - 43

Measurement: All measurements are circumferences. Ankle - Narrowest part (just above ankle bone). Calf - Widest part. Below Knee - 2cm below knee

Venosan 4000 Series - Nylon - Unisex

26 - 28

- Venosan 4000 is manufactured in a refined sandwich construction and is the first compression stocking in the world to be awarded the TACTEL climate effect label
- Suitable for: varicose veins, ulceration, thrombosis, sclerotherapy, pregnancy, work & travel

40 - 45

The Venosan 4000 Series is available in: knee high, thigh high, chaps, pantyhose styles

Venosan 4000 Series - Knee High

Size

Small

Large

Medium

X Large

		18 - 22	mmHg	23 - 32mmHg		
Size	Colour	Open Toe	Closed Toe	Open Toe	Closed Toe	
Small	Black	PTS641956OBLS	PTS641956CBLS	PTS641957OBLS	PTS641957CBLS	
Small	Marokko	PTS641956	PTS641956CMAS	PTS641957	PTS641957CMAS	
Medium	Black	PTS641956OBLM	PTS641956CBLM	PTS641957OBLM	PTS641957CBLM	
weatum	Marokko	PTS641956OMAM	PTS641956CMAM	PTS641957OMAM	PTS641957CMAM	
Lavaa	Black	PTS641956OBLL	PTS641956CBLL	PTS641957OBLL	PTS641957CBLL	
Large	Marokko	PTS641956OMAL	PTS641956CMAL	PTS641957OMAL	PTS641957CMAL	
V Lawrea	Black	PTS641956OBXL	PTS641956CBXL	PTS641957OBXL	PTS641957CBXL	
X Large	Marokko	PTS641956OMAXL	PTS641956CMAXL	PTS641957OMAXL	PTS641957CMAXL	

Venosan 4000 Series - Thigh High - Grip Top

		18 - 22	mmHg	23 - 32mmHg		
Size	Colour	Open Toe	Closed Toe	Open Toe	Closed Toe	
Creatil	Black	PTS641959OBLS	PTS641959PTBLS	PTS641960OBLS	PTS641960PTBLS	
Small	Marokko	PTS641959	PTS641959PTMAS	PTS641960	PTS641960PTMAS	
Medium	Black	PTS641959OBLM	PTS641959PTBLM	PTS641960OBLM	PTS641960PTBLM	
weatum	Marokko	PTS641959OMAM	PTS641959PTMAM	PTS641960OMAM	PTS641960PTMAM	
Lawre	Black	PTS641959OBLL	PTS641959PTBLL	PTS641960OBLL	PTS641960PTBLL	
Large	Marokko	PTS641959OMAL	PTS641959PTMAL	PTS641960OMAL	PTS641960PTMAL	
VIerre	Black	PTS641959OBXL	PTS641959PTBLXL	PTS641960OBXL	PTS641960PTBLXL	
X Large	Marokko	PTS641959OMAXL	PTS641959PTMAXL	PTS641960OMAXL	PTS641960PTMAXL	

Venosan 1000 Series - Chans - Open Toe

nes - Chaps - Open loe				
1	23 - 32mmHg			
Colour	Code (Left)	Code (Right)		
Mexico	PTS641958LS	PTS641958		
Mexico	PTS641958LM	PTS641958RM		
Mexico	PTS641958LL	PTS641958RL		
Mexico	PTS641958LXL	PTS641958RXL		
ries - Pantyhose - Closed Toe	, 🧭	DVA: AR22		
	18 - 22mmHg	23 - 32mmHg		
Colour	Code	Code		
Black	PTS641963BLS	PTS641964BLS		
Marokko	PTS641963	PTS641964		
Black	PTS641963BLM	PTS641964BLM		
Marokko	PTS641963MAM	PTS641964MAM		
Black	PTS641963BLL	PTS641964BLL		
Marokko	PTS641963MAL	PTS641964MAL		
Black	PTS641963BLXL	PTS641964BLXL		
Manual Jun	DTCC410C2MA VI	PTS641964MAXL		
	Mexico Mexico Mexico ries - Pantyhose - Closed Toe Black Marokko Black Marokko Black Marokko Black Marokko Black Marokko Black Marokko Black	Colour Code (Left) Mexico PTS641958LS Mexico PTS641958LM Mexico PTS641958LL Mexico PTS641958LL Mexico PTS641958LXL		



PTS641999XL



DVA: AR22

DVA: AR22

PRODUCT AND SERVICES CATALOGUE

Venosan 5000 Series - Cotton - Unisex

- The combination of the thermal properties of X-static, cotton and TACTEL provide the wearer with optimal comfort
- These compression garments are suitable for active men and women
- Styles available: knee high or chaps, in a choice of two colours •

Venosan 5000 Series - Knee High

		18 - 22	18 - 22mmHg		23 - 32mmHg		36 - 48mmHg	
Size	Colour	Open Toe	Closed Toe	Open Toe	Closed Toe	Open Toe	Closed Toe	
Care all	Black	PTS641965OBLS	PTS641965BLS	PTS641966OBLS	PTS641966BLS	PTS641967OBLS	PTS641967BLS	
Small	Marokko	PTS641965	PTS641965MAS	PTS641966	PTS641966MAS	PTS641967	PTS641967MAS	
Medium	Black	PTS641965OBLM	PTS641965BLM	PTS641966OBLM	PTS641966BLM	PTS641967OBLM	PTS641967BLM	
weatum	Marokko	PTS641965OMAM	PTS641965MAM	PTS641966OMAM	PTS641966MAM	PTS641967OMAM	PTS641967MAM	
Large	Black	PTS641965OBLL	PTS641965BLL	PTS641966OBLL	PTS641966BLL	PTS641967OBLL	PTS641967BLL	
Large	Marokko	PTS641965OMAL	PTS641965MAL	PTS641966OMAL	PTS641966MAL	PTS641967OMAL	PTS641967MAL	
VLarga	Black	PTS641965OBXL	PTS641965BLXL	PTS641966OBXL	PTS641966BLXL	PTS641967OBXL	PTS641967BLXL	
X Large	Marokko	PTS641965OMAXL	PTS641965MAXL	PTS641966OMAXL	PTS641966MAXL	PTS641967OMAXL	PTS641967MAXL	

Venosan 5000 Series - Chaps - Open Toe

		10 22	i i i i i g
Size	Colour	Code (Left)	Code (Right)
Small	Black	PTS641970LBLS	PTS641970RBLS
Siliali	Marokko	PTS641970LMAS	PTS641970
Medium	Black	PTS641970LBLM	PTS641970RBLM
Wedrum	Marokko	PTS641970LMAM	PTS641970RMAM
Lourse	Black	PTS641970LBLL	PTS641970RBLL
Large	Marokko	PTS641970LMAL	PTS641970RMAL
X Larga	Black	PTS641970LBLXL	PTS641970RBLXL
X Large	Marokko	PTS641970LMAXL	PTS641970RMAXL

Venosan 6000 Series - Nylon / Lycra - Unisex

- Fulfils the strict requirements of RAL certification
- Cost effective without compromise, optimal fit, easy care
- The Venosan 6000 Series is available in: knee high, thigh high, pantyhose styles
- Available in **black or beige**, with a choice of three compression ratings

Venosan 6000 Series - Knee High



		18 - 22	mmHg	23 - 32	mmHg	36 - 48mmHg
Size	Colour	Open Toe	Closed Toe	Open Toe	Closed Toe	Open Toe
Size	Colour			Short (34 - 39cm)		
Small	Black	PTS641974	PTS641974CSBLS	PTS641975	PTS641975CSBLS	-
Small	Beige	PTS641974OSBES	PTS641974CSBES	PTS641975OSBES	PTS641975CSBES	PTS641976OSBES
Medium	Black	PTS641974OSBLM	PTS641974CSBLM	PTS641975OSBLM	PTS641975CSBLM	-
weatum	Beige	PTS641974OSBEM	PTS641974CSBEM	PTS641975OSBEM	PTS641975CSBEM	PTS641976OSBEM
Lange	Black	PTS641974OSBLL	PTS641974CSBLL	PTS641975OSBLL	PTS641975CSBLL	-
Large	Beige	PTS641974OSBEL	PTS641974CSBEL	PTS641975OSBEL	PTS641975CSBEL	PTS641976OSBEL
N I	Black	PTS641974OSBLXL	PTS641974CSBLXL	PTS641975OSBLXL	PTS641975CSBLXL	-
X Large	Beige	PTS641974OSBEXL	PTS641974CSBEXL	PTS641975OSBEXL	PTS641975CSBEXL	PTS641976OSBEXL
				Long (39 - 44cm)		
Small	Black	PTS641974OLBLS	PTS641974CLBLS	PTS641975OLBLS	PTS641975CLBLS	-
SIIIdii	Beige	PTS641974OLBES	PTS641974CLBES	PTS641975OLBES	PTS641975CLBES	PTS641976OLBES
Medium	Black	PTS641974OLBLM	PTS641974CLBLM	PTS641975OLBLM	PTS641975CLBLM	-
Wealum	Beige	PTS641974OLBEM	PTS641974CLBEM	PTS641975OLBEM	PTS641975CLBEM	PTS641976OLBEM
Large	Black	PTS641974OLBLL	PTS641974CLBLL	PTS641975OLBLL	PTS641975CLBLL	-
Large	Beige	PTS641974OLBEL	PTS641974CLBEL	PTS641975OLBEL	PTS641975CLBEL	PTS641976OLBEL
X Large	Black	PTS641974OLBLXL	PTS641974CLBLXL	PTS641975OLBLXL	PTS641975CLBLXL	-
A Large	Beige	PTS641974OLBEXL	PTS641974CLBEXL	PTS641975OLBEXL	PTS641975CLBEXL	PTS641976OLBEXL

Measurement: To determine Short or Long, measure from the floor to 2cm below the knee

DVA Contact Details





PRODUCT AND SERVICES CATALOGUE





Venosan 6000 Series - Thigh High - Grip Top

						DVA: AR22
		18 - 22	mmHg	23 - 32	23 - 32mmHg	
0:	Calaur	Open Toe	Closed Toe	Open Toe	Closed Toe	Open Toe
Size	Colour			Short (59 - 67cm)		
Small	Black	PTS641977	PTS641977CSBLS	PTS641978	PTS641978CSBLS	-
Sman	Beige	PTS641977OSBES	PTS641977CSBES	PTS641978OSBES	PTS641978CSBES	PTS641979OSBES
Medium	Black	PTS641977OSBLM	PTS641977CSBLM	PTS641978OSBLM	PTS641978CSBLM	-
Medium	Beige	PTS641977OSBEM	PTS641977CSBEM	PTS641978OSBEM	PTS641978CSBEM	PTS641979OSBEM
1	Black	PTS641977OSBLL	PTS641977CSBLL	PTS641978OSBLL	PTS641978CSBLL	-
Large	Beige	PTS641977OSBEL	PTS641977CSBEL	PTS641978OSBEL	PTS641978CSBEL	PTS641979OSBEL
VIerre	Black	PTS641977OSBLXL	PTS641977CSBLXL	PTS641978OSBLXL	PTS641978CSBLXL	-
X Large	Beige	PTS641977OSBEXL	PTS641977CSBEXL	PTS641978OSBEXL	PTS641978CSBEXL	PTS641979OSBEXL
				Long (68 - 78cm)		
Small	Black	PTS641977OLBLS	PTS641977CLBLS	PTS641978OLBLS	PTS641978CLBLS	-
Sman	Beige	PTS641977OLBES	PTS641977CLBES	PTS641978OLBES	PTS641978CLBES	PTS641979
Medium	Black	PTS641977OLBLM	PTS641977CLBLM	PTS641978OLBLM	PTS641978CLBLM	-
Medium	Beige	PTS641977OLBEM	PTS641977CLBEM	PTS641978OLBEM	PTS641978CLBEM	PTS641979OLBEM
Lawaa	Black	PTS641977OLBLL	PTS641977CLBLL	PTS641978OLBLL	PTS641978CLBLL	-
Large	Beige	PTS641977OLBEL	PTS641977CLBEL	PTS641978OLBEL	PTS641978CLBEL	PTS641979OLBEL
Viana	Black	PTS641977OLBLXL	PTS641977CLBLXL	PTS641978OLBLXL	PTS641978CLBLXL	-
X Large	Beige	PTS641977OLBEXL	PTS641977CLBEXL	PTS641978OLBEXL	PTS641978CLBEXL	PTS641979OLBEXL

Measurement: To determine Short or Long, measure from the floor to 5cm below the groin

Venosan 6000 Series - Pantyhose

		CLOSED TOE ONLY				
Size	Colour	18 - 22mmHg	23 - 32mmHg	36 - 48mmHg		
Size	Colour		Short (67 - 74cm)			
Small	Black	PTS641980	PTS641981	-		
Sman	Beige	PTS641980CSBES	PTS641981CSBES	PTS641982CSBES		
Medium	Black	PTS641980CSBLM	PTS641981CSBLM	-		
weatum	Beige	PTS641980CSBEM	PTS641981CSBEM	PTS641982CSBEM		
1	Black	PTS641980CSBLL	PTS641981CSBLL	-		
Large	Beige	PTS641980CSBEL	PTS641981CSBEL	PTS641982CSBEL		
V Lawrence	Black	PTS641980CSBLXL	PTS641981CSBLXL	-		
X Large	Beige	PTS641980CSBEXL	PTS641981CSBEXL	PTS641982CSBEXL		
			Long (74 - 85cm)			
Small	Black	PTS641980CLBLS	PTS641981CLBLS	-		
Small	Beige	PTS641980CLBES	PTS641981CLBES	PTS641982		
Medium	Black	PTS641980CLBLM	PTS641981CLBLM	-		
weatum	Beige	PTS641980CLBEM	PTS641981CLBEM	PTS641982CLBEM		
1	Black	PTS641980CLBLL	PTS641981CLBLL	-		
Large	Beige	PTS641980CLBEL	PTS641981CLBEL	PTS641982CLBEL		
VLorgo	Black	PTS641980CLBLXL	PTS641981CLBLXL	-		
X Large	Beige	PTS641980CLBEXL	PTS641981CLBEXL	PTS641982CLBEXL		



×

Pressure Care & Circulation

Measurement: To determine Short or Long, measure from the floor to the groin

General Enquiries



Venosan 7000 Series - Opaque - Unisex

- Compression stockings with high stiffness and increased massage effect. Suitable for: lymphoedema/heavy swelling & ulcerated legs with compromised skin condition
- Antibacterial and fungicidal action inhibits the growth of bacteria and fungi



Venosan 7000 Series - Knee High

	23 - 32mmHg		36 - 48mmHg				
Colour	Open Toe	Closed Toe	Open Toe	Closed Toe			
Colour		Short (34	4 - 39cm)				
Black	PTS641983	PTS641983CSBLS	PTS641984	PTS641984CSBLS			
Beige	PTS641983OSBES	PTS641983CSBES	PTS641984OSBES	PTS641984CSBES			
Black	PTS641983OSBLM	PTS641983CSBLM	PTS641984OSBLM	PTS641984CSBLM			
Beige	PTS641983OSBEM	PTS641983CSBEM	PTS641984OSBEM	PTS641984CSBEM			
Black	PTS641983OSBLL	PTS641983CSBLL	PTS641984OSBLL	PTS641984CSBLL			
Beige	PTS641983OSBEL	PTS641983CSBEL	PTS641984OSBEL	PTS641984CSBEL			
Black	PTS641983OSBLXL	PTS641983CSBLXL	PTS641984OSBLXL	PTS641984CSBLXL			
Beige	PTS641983OSBEXL	PTS641983CSBEXL	PTS641984OSBEXL	PTS641984CSBEXL			
		Long (39) - 44cm)				
Black	PTS641983OLBLS	PTS641983CLBLS	PTS641984OLBLS	PTS641984CLBLS			
Beige	PTS641983OLBES	PTS641983CLBES	PTS641984OLBES	PTS641984CLBES			
Black	PTS641983OLBLM	PTS641983CLBLM	PTS641984OLBLM	PTS641984CLBLM			
Beige	PTS641983OLBEM	PTS641983CLBEM	PTS641984OLBEM	PTS641984CLBEM			
Black	PTS641983OLBLL	PTS641983CLBLL	PTS641984OLBLL	PTS641984CLBLL			
Beige	PTS641983OLBEL	PTS641983CLBEL	PTS641984OLBEL	PTS641984CLBEL			
Black	PTS641983OLBLXL	PTS641983CLBLXL	PTS641984OLBLXL	PTS641984CLBLXL			
Beige	PTS641983OLBEXL	PTS641983CLBEXL	PTS641984OLBEXL	PTS641984CLBEXL			
	Beige Black Beige Black Beige Black Beige Black Beige Black Beige Black Beige Black	Open ToeBlackPTS641983BeigePTS641983OSBESBlackPTS641983OSBEMBeigePTS641983OSBLMBeigePTS641983OSBLBlackPTS641983OSBLLBeigePTS641983OSBLLBiackPTS641983OSBLLBlackPTS641983OSBLSBlackPTS641983OLBLSBiackPTS641983OLBLSBiackPTS641983OLBLSBiackPTS641983OLBLSBiackPTS641983OLBLMBeigePTS641983OLBLMBiackPTS641983OLBLMBiackPTS641983OLBLLBiackPTS641983OLBLLBiackPTS641983OLBLL	Open ToeClosed ToeShort (3)BlackPTS641983BeigePTS641983OSBESBlackPTS641983OSBESBlackPTS641983OSBLMBeigePTS641983OSBLMBeigePTS641983OSBLMBiackPTS641983OSBLMBiackPTS641983OSBLBrigePTS641983OSBLBrigePTS641983OSBLBiackPTS641983OSBLBiackPTS641983OSBLBiackPTS641983OSBLBiackPTS641983OSBLBiackPTS641983OSBLXBrigePTS641983OSBLXBiackPTS641983OLBLSBiackPTS641983OLBLSBiackPTS641983OLBLSBiackPTS641983OLBLMPTS641983OLBLMPTS641983CLBLMBeigePTS641983OLBLMPTS641983OLBLPTS641983CLBLMBiackPTS641983OLBLPTS641983OLBLPTS641983CLBLLBiackPTS641983OLBLPTS641983OLBLPTS641983CLBLLBiackPTS641983OLBLPTS641983OLBLPTS641983CLBLLBiackPTS641983OLBLPTS641983OLBLXPTS641983CLBLL	Open ToeClosed ToeOpen ToeShort (34 - 39cm)BlackPTS641983PTS641983CSBLSPTS641984BeigePTS641983OSBESPTS641983CSBESPTS641984OSBESBlackPTS641983OSBLMPTS641983CSBLMPTS641984OSBLMBeigePTS641983OSBLMPTS641983CSBLMPTS641984OSBLMBeigePTS641983OSBLPTS641983CSBLPTS641984OSBLBlackPTS641983OSBLPTS641983CSBLPTS641984OSBLBiackPTS641983OSBLPTS641983CSBLPTS641984OSBLBiackPTS641983OSBLXPTS641983CSBLXPTS641984OSBLXLBiackPTS641983OSBEXLPTS641983CSBEXLPTS641984OSBEXLBiackPTS641983OSBEXLPTS641983CSBEXLPTS641984OSBEXLBiackPTS641983OLBLSPTS641983CLBLSPTS641984OLBLSBiackPTS641983OLBLSPTS641983CLBLSPTS641984OLBLSBiackPTS641983OLBLMPTS641983CLBLMPTS641984OLBLSBiackPTS641983OLBLMPTS641983CLBLMPTS641984OLBLMBiackPTS641983OLBLMPTS641983CLBLMPTS641984OLBLMBiackPTS641983OLBLMPTS641983CLBLMPTS641984OLBLMBiackPTS641983OLBLMPTS641983CLBLMPTS641984OLBLLBiackPTS641983OLBLPTS641983CLBLPTS641984OLBLLBiackPTS641983OLBLPTS641983OLBLPTS641984OLBLLBiackPTS641983OLBLPTS641983OLBLPTS641984OLBLLBiackPTS641983OLBLPTS641983OLBLPTS641984OLBLLBiack			

Measurement: To determine Short or Long, measure from the floor to 2cm below the knee

Venosan 7000 Series - Thigh High - Grip Top

	23 - 32r	nmHg	36 - 48	ßmmHg
Colour	Open Toe	Closed Toe	Open Toe	Closed Toe
		Short (59 - 67cm)		
Black	PTS641986	PTS641986CSBLS	PTS641987	PTS641987CSBLS
Beige	PTS641986OSBES	PTS641986CSBES	PTS641987OSBES	PTS641987CSBES
Black	PTS641986OSBLM	PTS641986CSBLM	PTS641987OSBLM	PTS641987CSBLM
Beige	PTS641986OSBEM	PTS641986CSBEM	PTS641987OSBEM	PTS641987CSBEM
Black	PTS641986OSBLL	PTS641986CSBLL	PTS641987OSBLL	PTS641987CSBLL
Beige	PTS641986OSBEL	PTS641986CSBEL	PTS641987OSBEL	PTS641987CSBEL
Black	PTS641986OSBLXL	PTS641986CSBLXL	PTS641987OSBLXL	PTS641987CSBLXL
Beige	PTS641986OSBEXL	PTS641986CSBEXL	PTS641987OSBEXL	PTS641987CSBEXL
		Long (68 - 78cm)		
Black	PTS641986OLBLS	PTS641986CLBLS	PTS641987OLBLS	PTS641987CLBLS
Beige	PTS641986OLBES	PTS641986CLBES	PTS641987OLBES	PTS641987CLBES
Black	PTS641986OLBLM	PTS641986CLBLM	PTS641987OLBLM	PTS641987CLBLM
Beige	PTS641986OLBEM	PTS641986CLBEM	PTS641987OLBEM	PTS641987CLBEM
Black	PTS641986OLBLL	PTS641986CLBLL	PTS641987OLBLL	PTS641987CLBLL
Beige	PTS641986OLBEL	PTS641986CLBEL	PTS641987OLBEL	PTS641987CLBEL
Black	PTS641986OLBLXL	PTS641986CLBLXL	PTS641987OLBLXL	PTS641987CLBLXL
Beige	PTS641986OLBEXL	PTS641986CLBEXL	PTS641987OLBEXL	PTS641987CLBEXL
	Black Beige Black Beige Black Beige Black Beige Black Beige Black Beige Black Beige Black Beige	ColourOpen ToeBlackPTS641986BeigePTS641986OSBESBlackPTS641986OSBEMBiackPTS641986OSBLMBiackPTS641986OSBLLBeigePTS641986OSBLLBiackPTS641986OSBELBlackPTS641986OSBELBlackPTS641986OSBELSBiackPTS641986OLBLSBeigePTS641986OLBLSBlackPTS641986OLBLSBlackPTS641986OLBLSBlackPTS641986OLBLSBlackPTS641986OLBLSBlackPTS641986OLBLMBeigePTS641986OLBLLBlackPTS641986OLBLLBiackPTS641986OLBLLBiackPTS641986OLBLLBiackPTS641986OLBLL	Short (59 - 67cm)BlackPTS641986PTS641986CSBLSBeigePTS641986OSBESPTS641986CSBESBlackPTS641986OSBLMPTS641986CSBLMBeigePTS641986OSBLMPTS641986CSBLMBlackPTS641986OSBLLPTS641986CSBLLBeigePTS641986OSBLLPTS641986CSBLLBeigePTS641986OSBLLPTS641986CSBLLBiackPTS641986OSBLXPTS641986CSBLXLBeigePTS641986OSBEXLPTS641986CSBLXLBeigePTS641986OLBLSPTS641986CLBLSBiackPTS641986OLBLSPTS641986CLBLSBiackPTS641986OLBLMPTS641986CLBLMBeigePTS641986OLBLMPTS641986CLBLMBiackPTS641986OLBLMPTS641986CLBLMBiackPTS641986OLBLPTS641986CLBLLBiackPTS641986OLBLPTS641986CLBLLBiackPTS641986OLBLPTS641986CLBLLBiackPTS641986OLBLLPTS641986CLBLLBiackPTS641986OLBLLPTS641986CLBLLBiackPTS641986OLBLLPTS641986CLBLLBiackPTS641986OLBLLPTS641986CLBLLBiackPTS641986OLBLLPTS641986CLBLLBiackPTS641986OLBLLPTS641986CLBLLBiackPTS641986OLBLLPTS641986CLBLLBiackPTS641986OLBLLPTS641986CLBLLBiackPTS641986OLBLLPTS641986CLBLLBiackPTS641986OLBLXPTS641986CLBLXL	ColourOpen ToeClosed ToeOpen ToeShort (59 - 67cm)BlackPTS641986PTS641986CSBLSPTS641987BeigePTS641986OSBESPTS641986CSBESPTS641987OSBESBlackPTS641986OSBLMPTS641986CSBLMPTS641987OSBLMBeigePTS641986OSBLMPTS641986CSBLMPTS641987OSBLMBeigePTS641986OSBLMPTS641986CSBLMPTS641987OSBLMBiackPTS641986OSBLLPTS641986CSBLLPTS641987OSBLLBiackPTS641986OSBLLPTS641986CSBLLPTS641987OSBLLBiackPTS641986OSBLXLPTS641986CSBLXLPTS641987OSBLXLBiackPTS641986OSBLXLPTS641986CSBEXLPTS641987OSBLXLBiackPTS641986OLBLSPTS641986CLBLSPTS641987OLBLSBiackPTS641986OLBLSPTS641986CLBLSPTS641987OLBLSBiackPTS641986OLBLSPTS641986CLBLSPTS641987OLBLSBiackPTS641986OLBLMPTS641986CLBLSPTS641987OLBLMBiackPTS641986OLBLMPTS641986CLBLMPTS641987OLBLMBiackPTS641986OLBLMPTS641986CLBLMPTS641987OLBLMBiackPTS641986OLBLMPTS641987OLBLMPTS641987OLBLMBiackPTS641986OLBLLPTS641986CLBLPTS641987OLBLLBiackPTS641986OLBLLPTS641986CLBLPTS641987OLBLLBiackPTS641986OLBLLPTS641986OLBLLPTS641987OLBLLBiackPTS641986OLBLLPTS641986CLBLLPTS641987OLBLLBiackPTS641986OLBLLPTS641986CLBLLPTS641987OLB

Measurement: To determine Short or Long, measure from the floor to 5cm below the groin





Compression - Legwear Pressure Care & Circulation

Venosan 7000 Series - Chaps - Open Toe

			DVA: AR22
		23 - 32	mmHg
Size	Colour	Code (Left)	Code (Right)
	Short (6	7 - 74cm)	
Small	Beige	PTS641985LSBES	PTS641985
Medium	Beige	PTS641985LSBEM	PTS641985RSBEM
Large	Beige	PTS641985LSBEL	PTS641985RSBEL
X Large	Beige	PTS641985LSBEXL	PTS641985RSBEXL
	Long (74	l - 85cm)	
Small	Beige	PTS641985LLBES	PTS641985RLBES
Medium	Beige	PTS641985LLBEM	PTS641985RLBEM
Large	Beige	PTS641985LLBEL	PTS641985RLBEL
X Large	Beige	PTS641985LLBEXL	PTS641985RLBEXL

Measurement: To determine Short or Long, measure from the floor to the groin

Venosan 'Legline' Series

- Venosan 'Legline 20' Series: suitable for everyday wear, work & travel. An effective compression stocking with an extremely fine mesh. Closed toe only
- Venosan 'Legline 30' Series: suitable for work, travel & post-sclerotherapy. An effective compression stocking with an extremely fine mesh. Closed toe only



PRODUCT AND SERVICES

Size Guide						
Circumference	Small	Medium	Large	X Large		
Ankle (cm)	22 - 23	24 - 25	26 - 27	28 - 29		
Calf (cm)	34 - 36	37 - 39	40 - 42	43 - 45		
Below Knee (cm)	32 - 34	35 - 37	38 - 40	41 - 43		
Thigh (cm)	45 - 47	48 - 50	51 - 53	54 - 60		

Measurement:

All measurements are circumferences. Ankle - Narrowest part (just above ankle bone). Calf - Widest part. Below Knee - 2cm below knee. Thigh - 5cm below groin

			CLOSED TOE ONLY					
		Venosan	'Legline 20' Series (2	20mmHg)	Venosan	'Legline 30' Series (3	0mmHg)	
		Knee high	Thigh high (Plain top)	Pantyhose	Knee high	Thigh high (Plain top)	Pantyhose	
Size	Colour			Cc	de			
Creatil	Black	PTS641988	PTS641990	PTS641992	PTS641989	PTS641991	PTS641993	
Small Nude	Nude	PTS641988NUS	PTS641990NUS	PTS641992NUS	PTS641989NUS	PTS641991NUS	PTS641993NUS	
NA - diama	Black	PTS641988BLM	PTS641990BLM	PTS641992BLM	PTS641989BLM	PTS641991BLM	PTS641993BLM	
Medium	Nude	PTS641988NUM	PTS641990NUM	PTS641992NUM	PTS641989NUM	PTS641991NUM	PTS641993NUM	
	Black	PTS641988BLL	PTS641990BLL	PTS641992BLL	PTS641989BLL	PTS641991BLL	PTS641993BLL	
Large	Nude	PTS641988NUL	PTS641990NUL	PTS641992NUL	PTS641989NUL	PTS641991NUL	PTS641993NUL	
X I	Black	PTS641988BLXL	PTS641990BLXL	PTS641992BLXL	PTS641989BLXL	PTS641991BLXL	PTS641993BLXL	
X Large	Nude	PTS641988NUXL	PTS641990NUXL	PTS641992NUXL	PTS641989NUXL	PTS641991NUXL	PTS641993NUXL	



Venosan 'Microfibre' Medical Legwear

- Cotton on the inside with a reinforced heel and toe
- Microfibre softness with a comfortable non-restricting top band
- **Closed toe only**

Size Guide							
Circumference		Small	Medium	Large	X Large		
Ankle (cm)	Female	16.5 - 20.5	20.5 - 24.5	24.5 - 28	28 - 32		
	Male	19.5 - 21.5	21.5 - 24	24 - 26.5	26.5 - 32		
Colf (am)	Female	26.5 - 37	29 - 39.5	32 - 43	34.5 - 47		
Calf (cm)	Male	28 - 38	33 - 43.5	35.5 - 46	38.5 - 51.5		
Shoe Size	Female	34 - 36	37 - 39	40 - 42	43 - 45		
	Male	37 - 39	40 - 43	44 - 47	48 - 51		



DVA: AR22

Measurement: All measurements are circumferences. Ankle - Narrowest part (just above ankle bone). Calf - Widest part. Shoe Size - European sizes

		Female	Male		
		15 -20mmHg	15 -20mmHg	18 -22mmHg	
Size	Colour		Code		
Small	Black	PTS641997BLS	PTS641996	PTS641998	
Silidii	Beige	PTS641997	-	-	
Medium	Black	PTS641997BLM	PTS641996M	PTS641998M	
weatum	Beige	PTS641997BEM	-	-	
Longo	Black	PTS641997BLL	PTS641996L	PTS641998L	
Large	Beige	PTS641997BEL	-	-	
Vierre	Black	PTS641997BLXL	PTS641996XL	PTS641998XL	
X Large	Beige	PTS641997BEXL	-	-	

Venosan 'Silverline' Medical Legwear (20 - 30mmHg)

- X-static® Silver yarn helps eliminate bacteria and fungi in the socks associated with foot odour, athlete's foot and fungal nail disease. Anti-microbial effect will not wash or wear out
- Superior moisture wicking for dry, comfortable feet and legs. Cool in summer, warm in winter

	Size Guide							
Circumference		Small	Medium	Large	X Large			
	Female	16.5 - 20.5	20.5 - 24.5	24.5 - 28	28 - 32			
Ankle (cm)		19.5 - 21.5	21.5 - 24	24 - 26.5	26.5 - 32			
Calf (ana)	Female	26.5 - 37	29 - 39.5	32 - 43	34.5 - 47			
Calf (cm)	Male	28 - 38	33 - 43.5	35.5 - 46	38.5 - 51.5			
Shoe Size	Female	34 - 36	37 - 39	40 - 42	43 - 45			
	Male	37 - 39	40 - 43	44 - 47	48 - 51			

Measurement: All measurements are circumferences. Ankle - Narrowest part (just above ankle bone). Calf - Widest part. Shoe Size - European sizes

			DVA: AR22
		CLOSED TOE ONLY	
		Female	Male
Size	Colour	Co	de
Small	Black	PTS641995BLS	PTS641994
Smail	Beige	PTS641995	-
Medium	Black	PTS641995BLM	PTS641994M
weatum	Beige	PTS641995BEM	-
Laura	Black	PTS641995BLL	PTS641994L
Large	Beige	PTS641995BEL	-
X Large	Black	PTS641995BLXL	PTS641994XL
A Large	Beige	PTS641995BEXL	-



DVA: AR22

DVA Contact Details







PRODUCT AND SERVICES CATALOGUE

Stocking Accessories

Easy 'Sim - Slide' On / Off Aid (Open Toe)

-	
	- 2
	-
	-
ARION	

- The ultimate donning and doffing aid for open toe compression stockings and tights
- Reduces the physical strain for both the caregiver and user

Available in 4 sizes

	DVA: AU13	
Size	Shoe Size	Code
Small	34 - 37	PTS642508
Medium	38 - 41	PTS642509
Large	42 - 45	PTS642510
X Large	46	PTS642511
Ches Cire European sizes		

Easy Slide 'Magnide' On / Off Aid (Closed Toe)



- The easiest donning and doffing aid for closed toe compression stockings and tights
- Reduces the physical strain for • both the caregiver and user

ARION	DVA: AU13	
Size	Shoe Size	Code
Medium / Large	40 - 45	PTS642514

Shoe Size - European sizes

Shoe Size - European sizes

Easy Slide 'Magnide' Applicator (Closed Toe)



- The ultimate donning aid for closed toe compression stockings and tights
- Reduces the physical strain for both the caregiver and user

ARION	DVA: AU13		
Size	Shoe Size	Code	
Small / Medium	40	PTS642512	
Large / X Large	41	PTS642513	
Shoe Size - European sizes			

Medi-Butler Stocking Donner



Shoe Size - European sizes

Medi-Butler Off - Doffing Aid



- This long handled device assists with the removal of compression stockings, socks and normal pantyhose
- It can also double as a shoehorn

	DVA: AU13
Code	Size
PTS642543	Universal

JOBST Stocking Donner



- Application aid for closed-toe garments
- Slide stockings over, then remove from top of stocking
- Suggest to use with rubber gloves

Code	Size	Shoe Size	DVA:
PTS641090	Small / Med	Up to 38	AU13
Note: Shoe size - European sizes			

Sigvaris Doff N' Donner



- Soft device that allows you to quickly guide a compression garment up and down your limb
- The Cone acts as a loading device to make using the Sigvaris Doff N' Donner even easier
- Sold seperately or as a set

D\/A· A1113

Code	Description
PTS642660	Doff N' Donner Device
PTS642661	Doff N' Donner Cone (Black)
PTS642662	Doff N' Donner Device & Cone

Jobst Donning Gloves



- Jobst gloves help create friction, provide grip and help to evenly distribute the fabric on the leg
- This product contains natural rubber Latex

DVA: AU13	
Description	
FBS303490 Jobst Donning Gloves - Small	
Jobst Donning Gloves - Medium	



Stocking Accessories

Jobst Body Adhesive



- Roll-on body adhesive helps Jobst garments stay in place Easy application
- Washes off easily with water without leaving a stain

	DVA: AR26	
Code	Size	
PTS641100	59ml	

Compression Stocking Aid



- Applies support hosiery quickly and evenly
- Grooves catch the compression stocking and prevent it from popping off
- Works like a traditional sock aid - pull the stocking over the shell, insert the foot and pull the cord to bring the hose up onto the leg
- Suitable for Class 1 and 2 stockings

Code	Description	Warranty	DVA:
DLD269110	Compression Stocking Aid - Small	1 Year	AU13
DLD269120	Compression Stocking Aid - Large	1 Year	AU13



- Hypoallergenic skin lubricant makes donning your compression garments an easier process
- Spray bottle for easy application

Code
GMS336050



Compression Stocking Aid



- EUREKA-ON Plastic Fabric Slide is a compression stocking applicator for the leg or arm
- Consists of a triangular shaped durable paper-like slippery material with a central extended pull tab and evelet holes on the sides for threading a thin tubular plastic lacing pin

Code	Description	Warranty	DVA:
DLD269010	Compression Stocking Aid - Eureka	1 Year	AU13



EzyAs Compression Stocking Applicator

- Unique plastic gutter to assist in pulling on compression stockings
- Stocking is fed onto gutter, held in place by side notches, slipped over foot, then pulled into place

Code	Description	Width	Length	Calf/Arm Circ.	Warranty	DVA:
DLD269001	Small - Red	140 mm	140 mm	280 - 350 mm	1 Year	AU13
DLD269002	Med - Yellow	140 mm	140 mm	320 - 400 mm	1 Year	AU13
DLD269003	Large - Blue	140 mm	175 mm	380 - 470 mm	1 Year	AU13
DLD269004	Handle	210 mm	270 mm	-	1 Year	AU13







481

Compression - Armwear



Compression Armsleeve Overview

What is Lymphoedema?

- If the lymphatic system is damaged or blocked, excess fluid, which would normally be drained away by the lymphatics, can build up in the tissues (most co mmonly the arms or legs)
- Lymphoedema in the upper limb is frequently a consequence of surgery or radiotherapy relating to breast cancer, resulting in fluid build-up in the arm, and often the hand

Importance of compression in the treatment of Lymphoedema

Compression is an important adjunct to ongoing professional treatment of the affected limb. Compression plays a major role in the treatment of lymphoedema by encouraging lymph flow.

Your lymphoedema therapist will advise you on the style and compression appropriate for your particular needs.

	INDICATORS					
	20-30 mmHG (class 1)	30-40 mmHg (class 2)				
•	• Treatment of primary and secondary lymphoedema in stages 1 and 2 of the arm with permanent or partial oedema of the hand					
•	Post-operative (eg. mastectomy) and post-traumatic oedema					

Blood poisoning and burn therapy

Styles Available

- The diagrams below depict the most co mmonly used styles of arm sleeves and hand pieces
- Not all of these styles are available in every brand the type of top selected (plain or grip) is a personal preference





Styles Available (contd)



Fitting & Selection Tips

- Never allow the top of the sleeve to roll or double-over. If there seems to be excess material, ease it gently down the • arm until the top is flat against the skin, and ensure no creases or wrinkles remain
- Ensure correct measurements are taken, and that the garment is not too tight or loose
- Advise the therapist of any pain, swelling or colour change of the hand or fingers, or numbness in any part of the limb

	DVA:	DVA Approval Rules	Eligible Prescriber*
	AR22	Limit of 3 every 6 months	S, RN, Physio, OT, LMO
I	OVA Co mments	*For lymphoedema garments, health providers should have speciali order to access, measure, fit and review these garments.	ist post graduate training in oedema/lymphoedema management in



Isotoner Therapeutic Pressure Gloves (pair)

- Washable lycra and spandex blend gives Isotoner gloves a comfortable fit while controlling swelling
- Use this glove anytime day or night and the gentle, even compression will sooth and reduce swollen aching hands
- Seamless support-seams are on the outside and away from pain-sensitive areas like the base of the thumb

		DV,	A: AR22
Size	Palm (cm)	Code (open finger)	Code (full finger)
Small	17 - 19	FBS301159	FBS301142
Medium	19 - 22	FBS301160	FBS301143
Large	22 - 24	FBS301161	FBS301144

Measure (cm) around palm (widest part of hand)

- 2-way stretch provides gentle compression
- Isotoner therapeutic glove compression is considered moderate: 23 - 32mmHg

DVA: AR22



Norco Oedema Glove (each)

			Full F	Open Finger			
		1. Wrist		2. Over Wrist		3. Over Wrist	
Size	Palm (cm)	Left	Right	Left	Right	Left	Right
XS	18	-	-	FBS301151	FBS301152	FBS301162	FBS301163
S	20	FBS301145	FBS301146	FBS301153	FBS301154	FBS301164	FBS301165
Μ	23	FBS301147	FBS301148	FBS301155	FBS301156	FBS301166	FBS301167
L	25	FBS301149	FBS301150	FBS301157	FBS301158	FBS301168	FBS301169

Measure (cm) around palm (widest part of hand).

• Made of nylon/lycra with seams on the outside for smooth oedema control

• Mild compression 15 - 25mmHg

- · Closed finger style available in wrist length or over-the-wrist
- Open finger styles available in over-the-wrist and full arm styles
- Sold individually



For NDIS Item Code



483





Handeze Support Gloves (pair)

- The unique four-way stretch fabric provides three main proven pain relievers: heat, support and massage
- The Handeze glove traps natural body heat, which stimulates blood circulation in wrist, hand and fingers. This reduces swelling and helps maintain flexibility
- Made of modified lycra fabric
- The form-fitting design provides additional support for the muscles, tendons and ligaments of the hand and wrist

	DVA: AR22		
Knuckles (cm)	Colour	Code	
Lin to 7.0	Beige	FBS301170	
Up to 7.0	Slate	FBS301170G0	
6.5 - 8.0	Beige	FBS301170S1	
	Slate	FBS301170G1	
n 7.5 - 9.5	Beige	FBS301170S2	
	Slate	FBS301170G2	
9 5 11 0	Beige	FBS301170S3	
8.5 - 11.0	Slate	FBS301170G3	
	(cm) Up to 7.0 6.5 - 8.0	Knuckles (cm)ColourUp to 7.0Beige6.5 - 8.0Beige7.5 - 9.5Beige8.5 - 11.0Beige	

Measure the width hand across the four knuckles (not finger joints). Do not include the thumb

Jobst - Medical Armwear

Overview

Jobst upper limb garments focus on improving the quality of life for those suffering from lymphatic conditions

To Order

Choose your preferred model, then use the Size Guide to select the appropriate size and code

Jobst Bella Lite

- Soft and air-permeable
- Lightweight knit construction

General Enquiries

• Fits a wide range of hand and arm shapes

- Fits both the left or right arm/hand
- The combination sleeve/glove can be worn separately together

Jobst Bella Lite - Standard Arm Sleeve (20 - 30mmHg)

					DVA: AR22	
	Size Guide			Code		
Size	Wrist (cm)	Mid-Lower Arm (cm)	Mid-Upper Arm (cm)	Regular (length up to 43cm)	Long (length over 43cm)	
Small	14 - 20	18 - 23	22 - 40	FBS303479	FBS303478	
Medium	16 - 22	23 - 29	27 - 45	FBS303481	FBS303480	
Large	18 - 24	29 - 34	33 - 51	FBS303483	FBS303482	

Measurement: Wrist and arm measurements are circumferences. Length is from axilla to wrist

Jobst Bella Lite - Combination Sleeve/Glove (20 - 30mmHg)

					DVA: AR22	
	Size Guide			Code		
Size	Palm (cm)	Mid-Lower Arm (cm)	Mid-Upper Arm (cm)	Regular (length up to 43cm)	Long (length over 43cm)	
Small	15 - 19	18 - 23	22 - 40	FBS303484	FBS303487	
Medium	19 - 22	23 - 29	27 - 45	FBS303485	FBS303488	
Large	22 - 26	29 - 34	33 - 51	FBS303486	FBS303489	

Measurement: Palm and arm measurements are circumferences. Length is from axilla to wrist



Compression - Armwear Pressure Care & Circulation



Jobst Bella Strong

- Ideal for managing moderate upper extremity lymphoedema
- Donning device included with every sleeve for easy donning
- Open knit structure allows for improved air permeability which keeps skin dry
- Two-way stretch makes the arm sleeve comfortable without compromising medical efficacy
- Circular knit construction
- Cotton-like feel; feels soft to the skin
- Dotted, 2" silicone band keeps sleeve firmly in place

Jobst Bella Strong - Standard Arm Sleeve (20 - 30mmHg)

				VA: AR22	
		Size Guide	Co	de	
Size	Wrist (cm)	Mid-Forearm (cm)	Axilla (cm) (Top of Arm)	Regular (length up to 43cm)	Long (length over 43cm)
1	14 - 18	17 - 21	22 - 31	FBS303458	FBS303468
2	14 - 18	17 - 21	31 - 41	FBS303459	FBS303469
3	14 - 18	21 - 25	25 - 34	FBS303460	FBS303470
4	16 - 19	21 - 25	28 - 38	FBS303461	FBS303471
5	17 - 20	21 - 25	34 - 43	FBS303462	FBS303472
6	14 - 18	25 - 29	29 - 38	FBS303463	FBS303473
7	17 - 20	25 - 29	33 - 43	FBS303464	FBS303474
8	18 - 21	25 - 29	38 - 48	FBS303465	FBS303475
9	16 - 19	29 - 33	32 - 42	FBS303466	FBS303476
10	19 - 22	29 - 33	36 - 46	FBS303467	FBS303477



Measurement: Wrist, arm and axilla measurements are circumferences. Length is from axilla to wrist

Jobst Bella Strong - Gauntlet (20 - 30mmHg)

		(DVA: AR22
	Size Guid	le	Code
Size	Wrist (cm)	Palm (cm)	Code
1	14 - 18	15 - 17	PTS642534
2	14 - 18	17 - 19	PTS642535
3	18 - 21	17 - 19	PTS642536
4	14 - 18	19 - 21	PTS642537
5	18 - 21	19 - 21	PTS642538
6	16 - 19	21 - 23	PTS642539
7	19 - 22	21 - 23	PTS642540
8	19 - 22	23 - 26	PTS642541



Measurement: Wrist and palm measurements are circumferences



I





(√)

Jobst Elvarex

- JOBST Elvarex garments feature a special low-stretch, 'double-knit' fabric, which exerts a light massaging effect on the initial lymphatics, and minimises re-oedematisation
- The high cotton content, flat knit construction and short stretch thread means that Elvarex products provide a good anatomical fit with maximum comfort
- A variety of styles available including: Standard Arm Sleeve; Arm Sleeve with Bra Attachment; Arm Sleeve with Shoulder Strap; Hand Piece with and without Fingers
- Available in **beige**

Jobst Elvarex - Standard Arm Sleeve (20 - 30mmHg)



DVA: AR22

JOD21	Jobsi Elvalex - Siandala Anni Sieeve (20 - Summing)						ARZZ
	Size Guide					Code	
Size	Wrist (cm) (Narrowest Part)	Mid-Lower Arm (cm)	Elbow Crease (cm)	Mid-Upper Arm (cm)	Top of Arm (cm) (Below Armpit)	Left	Right
Small	16 - 18	23 - 26	24 - 26	26 - 29	31 - 34	PTS642150	PTS642149
Medium	18 - 20	26 - 28	26 - 29	29 - 33	34 - 38	PTS642152	PTS642151
Large	20 - 23	28 - 30	29 - 32	33 - 37	38 - 44	PTS642154	PTS642153

Measurement: All measurements are circumferences

Jobst Elvarex - Arm Sleeve with Bra Attachment (20 - 30mmHg)

Size				Code			
	Wrist (cm) (Narrowest Part)	Mid-Lower Arm (cm)	Elbow Crease (cm)	Mid-Upper Arm (cm)	Top of Arm (cm) (Below Armpit)	Left	Right
Small	16 - 18	23 - 26	24 - 26	26 - 29	31 - 34	PTS642162	PTS642161
Medium	18 - 20	26 - 28	26 - 29	29 - 33	34 - 38	PTS642164	PTS642163
Large	20 - 23	28 - 30	29 - 32	33 - 37	38 - 44	PTS642166	PTS642165

Measurement: All measurements are circumferences

Jobst Elvarex - Arm Sleeve with Shoulder Strap (20 - 30mmHg) OVA: AR22										
		Code								
Size	Wrist (cm) (Narrowest Part)	Mid-Lower Arm (cm)	Elbow Crease (cm)	Mid-Upper Arm (cm)	Top of Arm (cm) (Below Armpit)	Left	Right			
Small	16 - 18	23 - 26	24 - 26	26 - 29	31 - 34	PTS642156	PTS642155			
Medium	18 - 20	26 - 28	26 - 29	29 - 33	34 - 38	PTS642158	PTS642157			
Large	20 - 23	28 - 30	29 - 32	33 - 37	38 - 44	PTS642160	PTS642159			

Measurement: All measurements are circumferences

Jobst E	Ivarex - Hand	D'	VA: AR22				
		Size (Guide		Code		
Size	Hand (cm) (at knuckles)	Palm (cm) (widest part)	Wrist (cm) (narrowest part)	Thumb (cm) (at base)	Left	Right	
Small	16 - 18	17 - 19	16 - 18	5 - 6	PTS642168	PTS642167	
Medium	17.5 - 19.5	19 - 21	18 - 20	6 - 7	PTS642170	PTS642169	
Large	19.5 - 20.5	21 - 23	20 - 22	7 - 8	PTS642172	PTS642171	

Measurement: All measurements are circumferences

Jobst Elvarex - Handpiece with Fingers (20 - 30mmHg)

• Note: For correct sizing, use the Handpiece Size Guide above and the below Size Guide

001011	0120 00100					
		Size (Code			
Size	Index Finger (cm) (at base)	Middle Finger (cm) (at base)	Ring Finger (cm) (at base)	Little Finger (cm) (at base)	Left	Right
Small	5 - 6	5 - 6	4.5 - 5.5	4 - 5	PTS642174	PTS642173
Medium	6 - 7	6 - 7	5.5 - 6.5	5 - 6	PTS642176	PTS642175
Large	6 - 7	6.5 - 7	6.5 - 7.5	6 - 7	PTS642178	PTS642177

Measurement: All measurements are circumferences

General Enquiries

Mediven

- Mediven Harmony uses a unique fabric that is soft to touch, clinically effective and easy to apply
- It features a new design that incorporates a double ellipse - one in the crease and one on the point of the elbow - for maximum patient comfort
- Seamless, soft fibre offers a comfortable fit
- RAL compression classes 1 and 2

To Order

Choose your preferred model, then use the appropriate Size Guide to select the appropriate size and code

Mediven Harmony - Arm Sleeve (CG)

- CG Arm Sleeves stop at the wrist
- The CG Arm Sleeve is available in 3 styles: Arm Sleeve, Arm Sleeve with Top Band, Arm Sleeve with Shoulder cap
- Colour available: Caramel
- Compression: 18 21 mmHg and 23 32mmHg







		Size Guide										
	Size 1	Size 2	Size 3	Size 4	Size 5	Size 6	Size 7	Size 8				
Top of Arm (cm)	25.5-28	28 - 30.5	30 - 33	32 - 35	33.5 - 36.5	36 - 39.5	38 - 41.5	40.5 - 44.5				
Elbow (cm)	23 - 25	24.5 - 26.5	26 - 28.5	27.5 - 30	29 - 31.5	30.5 - 33.5	32.5 - 35.5	35 - 38.5				
Wrist (cm)	13.5 - 14.5	14.5 - 15.5	16 - 17.5	17.5 - 19	19 - 20.5	20.5 - 22.5	22.5 - 24.5	25 - 27.5				

Measurement: All measurements are circumferences and are in cm

						OVA: AR22	
		Class 1 (18 - 21 mmHg)		Class 2 (23 - 32 mmHg)			
	CG Arm Sleeve	CG Arm Sleeve with Top Band	CG Arm Sleeve with Shoulder cap	CG Arm Sleeve	CG Arm Sleeve with Top Band	CG Arm Sleeve with Shoulder cap	
Size			Co	de			
1	FBS303000	FBS303280	FBS303152	FBS303032	FBS303312	FBS303216	
2	FBS303001	FBS303281	FBS303153	FBS303033	FBS303313	FBS303217	
3	FBS303002	FBS303282	FBS303154	FBS303034	FBS303314	FBS303218	
4	FBS303003	FBS303283	FBS303155	FBS303035	FBS303315	FBS303219	
5	FBS303004	FBS303284	FBS303156	FBS303036	FBS303316	FBS303220	
6	FBS303005	FBS303285	FBS303157	FBS303037	FBS303317	FBS303221	
7	FBS303006	FBS303286	FBS303158	FBS303038	FBS303318	FBS303222	
8	FBS303007	FBS303287	FBS303159	FBS303039	FBS303319	FBS303223	





For NDIS Item Code to page: 500









Mediven Harmony - Arm Sleeve (CG) - Extra Wide

- The Extra Wide CG Arm Sleeve can acco mmodate increased swelling
- CG Arm Sleeves stop at the wrist
- The CG Arm Sleeve is available in 3 styles: Arm Sleeve, Arm Sleeve with Top Band, Arm Sleeve with Shoulder cap
- Colour available: Caramel
- Compression: 18 21mmHg and 23 32mmHg







Ann Sleeve w Shoulder ca

		Size Guide									
	Size 1	Size 2	Size 3	Size 4	Size 5	Size 6	Size 7	Size 8			
Top of Arm (cm)	28.5- 31	30.5 - 33.5	33 - 36	35 - 38.5	37 - 40.5	39.5 - 43	41.5 - 45.5	45 - 49.5			
Elbow (cm)	24 - 26	26.5 - 29	28.5 - 31	30 - 33	31.5 - 34.5	33.5 - 36.5	35.5 - 39	38.5 - 42			
Wrist (cm)	13.5 - 14.5	14.5 - 15.5	16 - 17.5	17.5 - 19	19 - 20.5	20.5 - 22.5	22.5 - 24.5	25 - 27.5			

Measurement: All measurements are circumferences and are in cm

						OVA: AR22	
		Class 1 (18 - 21 mmHg)		Class 2 (23 - 32 mmHg)			
	CG Arm Sleeve (Extra Wide)	CG Arm Sleeve with Top Band (Extra Wide)	CG Arm Sleeve with Shoulder cap (Extra Wide)	CG Arm Sleeve (Extra Wide)	CG Arm Sleeve with Top Band (Extra Wide)	CG Arm Sleeve with Shoulder cap (Extra Wide)	
Size			Co	de			
1	FBS303016	FBS303296	FBS303168	FBS303048	FBS303328	FBS303232	
2	FBS303017	FBS303297	FBS303169	FBS303049	FBS303329	FBS303233	
3	FBS303018	FBS303298	FBS303170	FBS303050	FBS303330	FBS303234	
4	FBS303019	FBS303299	FBS303171	FBS303051	FBS303331	FBS303235	
5	FBS303020	FBS303300	FBS303172	FBS303052	FBS303332	FBS303236	
6	FBS303021	FBS303301	FBS303173	FBS303053	FBS303333	FBS303237	
7	FBS303022	FBS303302	FBS303174	FBS303054	FBS303334	FBS303238	
8	FBS303023	FBS303303	FBS303175	FBS303055	FBS303335	FBS303239	

Mediven Harmony - Combined Arm Sleeve (AG)

- AG Arm Sleeves incorporate a gauntlet
- The AG Arm Sleeve is available in 3 styles: Arm Sleeve, Arm Sleeve with Top Band, Arm Sleeve with Shoulder cap
- Colour available: Caramel
- Compression: 18 21mmHg and 23 32mmHg



		Size Guide									
	Size 1	Size 2	Size 3	Size 4	Size 5	Size 6	Size 7	Size 8			
Top of Arm (cm)	25.5- 28	28 - 30.5	30 - 33	32 - 35	33.5 - 36.5	36 - 39.5	38 - 41.5	40.5 - 44.5			
Elbow (cm)	23 - 25	24.5 - 26.5	26 - 28.5	27.5 - 30	29 - 31.5	30.5 - 33.5	32.5 - 35.5	35 - 38.5			
Wrist (cm)	13.5 - 14.5	14.5 - 15.5	16 - 17.5	17.5 - 19	19 - 20.5	20.5 - 22.5	22.5 - 24.5	25 - 27.5			
Palm (cm)	15- 18	16 - 19	17 - 20	19 - 22	20 - 23	22 - 25	24 - 27	26 - 29			

Measurement: All measurements are circumferences and are in cm



DVA: AR22

Mediven Harmony - Combined Arm Sleeve (AG) (contd)

Class 2 (23 - 32 mmHg) AG Arm Sleeve with Top Band AG Arm Sleeve with Shoulder cap AG Arm Sleeve with Shoulder cap Top Band Code 1 FBS303088 FBS303344 FBS303184 FBS303120 FBS303376 FBS303248 2 FBS303089 FBS303345 FBS303185 FBS303121 FBS303377 FBS303249 3 FBS303090 FBS303346 FBS303186 FBS303122 FBS303378 FBS303250 4 FBS303091 FBS303347 FBS303187 FBS303123 FBS303379 FBS303251 5 FBS303092 FBS303348 FBS303188 FBS303124 FBS303380 FBS303252 6 FBS303093 FBS303349 FBS303189 FBS303125 FBS303381 FBS303253 7 FBS303094 FBS303350 FBS303190 FBS303126 FBS303382 FBS303254 FBS303095 FBS303351 FBS303191 FBS303127 FBS303383 FBS303255 8

Mediven Harmony - Combined Arm Sleeve (AG) - Extra Wide

- The Extra Wide AG Arm Sleeve can acco mmodate increased swelling •
- AG Arm Sleeves incorporate a gauntlet
- The AG Arm Sleeve is available in 3 styles: Arm Sleeve, Arm Sleeve with Top Band, Arm Sleeve with Shoulder cap
- Colour available: Caramel
- Compression: 18 - 21mmHg and 23 - 32mmHg



		Size Guide										
	Size 1	Size 2	Size 3	Size 4	Size 5	Size 6	Size 7	Size 8				
Top of Arm (cm)	28.5- 31	30.5 - 33.5	33 - 36	35 - 38.5	37 - 40.5	39.5 - 43	41.5 - 45.5	45 - 49.5				
Elbow (cm)	24 - 26	26.5 - 29	28.5 - 31	30 - 33	31.5 - 34.5	33.5 - 36.5	35.5 - 39	38.5 - 42				
Wrist (cm)	13.5 - 14.5	14.5 - 15.5	16 - 17.5	17.5 - 19	19 - 20.5	20.5 - 22.5	22.5 - 24.5	25 - 27.5				
Palm (cm)	15- 18	16 - 19	17 - 20	19 - 22	20 - 23	22 - 25	24 - 27	26 - 29				

Measurement: All measurements are circumferences and are in cm

						JVA: ARZZ	
		Class 1 (18 - 21 mmHg)		Class 2 (23 - 32 mmHg)			
	AG Arm Sleeve AG Arm Sleeve with AG (Extra Wide) (Extra Wide)		AG Arm Sleeve with Shoulder cap (Extra Wide)	AG Arm Sleeve (Extra Wide)	AG Arm Sleeve with Top Band (Extra Wide)	AG Arm Sleeve with Shoulder cap (Extra Wide)	
Size			Co	ode			
1	FBS303104	FBS303360	FBS303200	FBS303136	FBS303392	FBS303264	
2	FBS303105	FBS303361	FBS303201	FBS303137	FBS303393	FBS303265	
3	FBS303106	FBS303362	FBS303202	FBS303138	FBS303394	FBS303266	
4	FBS303107	FBS303363	FBS303203	FBS303139	FBS303395	FBS303267	
5	FBS303108	FBS303364	FBS303204	FBS303140	FBS303396	FBS303268	
6	FBS303109	FBS303365	FBS303205	FBS303141	FBS303397	FBS303269	
7	FBS303110	FBS303366	FBS303206	FBS303142	FBS303398	FBS303270	
8	FBS303111	FBS303367	FBS303207	FBS303143	FBS303399	FBS303271	

DVA Contact Details









- The AC1 Hand Part is available with or without compressive fingers •
- Colour available: Caramel

	Size Guide									
	Size 2	Size 3	Size 4	Size 5	Size 6	Size 7				
Wrist	13.5 - 15.5	15 - 18	17 - 20	19 - 21	20 - 23.5	22.5 - 27.5				
Palm	15 - 19	17 - 20	19 - 22	20 - 23	22 - 25	24 - 29				

Measurement: All measurements are circumferences and are in cm

				DVA: AR22	
	Class 1 (18	- 21 mmHg)	Class 2 (23 - 32 mmHg)		
	Hand Part (with fingers)	Hand Part (without fingers)	Hand Part (with fingers)	Hand Part (without fingers)	
Size	Code				
2	FBS303408	FBS303064	FBS303420	FBS303076	
3	FBS303409	FBS303065	FBS303421	FBS303077	
4	FBS303410	FBS303066	FBS303422	FBS303078	
5	FBS303411	FBS303067	FBS303423	FBS303079	
6	FBS303412	FBS303068	FBS303424	FBS303080	
7	FBS303413	FBS303069	FBS303425	FBS303081	

Sigvaris

To Order

Choose your preferred model. Then, order using: Code, Description, Style, Compression (mmHg), Size

Sigvaris - Traditional - Unisex

- The Traditional series is a reliable safe and effective option
- The elasticity and breathability mean that they are very comfortable to wear
- Colour available: Beige

	Size Guide						
Size	Wrist (cm) (Narrowest Part)	Forearm (cm) (Mid-Forearm)	Upper Arm (cm) (Thickest Part)	Palm (cm (Around Knud			
S	13 - 16	20 - 24	24 - 30	16 - 19			
Μ	16 - 19	24 - 28	30 - 36	19 - 22			
L	19 - 22	28 - 32	36 - 42	22 - 26			
Aeasurement: All measurements are circumferences. Palm measurement only required for arm sleeves with mittens							

Compression 23 - 32 mmHg (Class 2) Without Mitten ~ S/M/L PTS642415 Traditional Sleeve PTS642420 Traditional Sleeve With Mitten ~ S/M/L PTS642445 Traditional Mitten S/M/L

1



Pressure Care & Circulation

490









SIGVARIS

PRODUCT AND SERVICES CATALOGUE



PRODUCT AND SERVICES

Sigvaris - Advance - Unisex

- The Advance series has been developed with the help of patients and clinicians to create a garment that meets medical efficacy standards and exceeds expectations
- Constructed with a polyamide microfibre for extreme softness and anti-bacterial properties
- Colour available: Light beige

	Size Guide						
Size	Wrist (cm) (narrowest part)	Mid-Upper Arm (cm) (mid-way between elbow & axilla)	Palm (cm) (around knuckles)				
XS	13 - 14	24 - 32	16 - 20				
S	14 - 16	26 - 34	17 - 21				
Μ	16 - 18	29 - 37	18 - 22				
L	18 - 20	32 - 40	19 - 23				
XL	20 - 23	34 - 42	20 - 24				

Measurement: All measurements are circumferences. Palm measurement only required for arm sleeves with mittens

			Compressio	on (mmHg)		
Code	Description	Style	14 - 18	20 - 25	Size	
PTS642425	Advance Sleeve	Without Mitten	\checkmark	~	XS / S / M / L / XL	
PTS642430	Advance Sleeve	With Mitten	\checkmark	\checkmark	XS / S / M / L / XL	
PTS642435	Advance Sleeve with Grip Top	Without Mitten	\checkmark	\checkmark	XS / S / M / L / XL	
PTS642440	Advance Sleeve with Grip Top	With Mitten	\checkmark	\checkmark	XS/S/M/L/XL	

Solaris

Overview

- The Solaris brand offers a wide range of compression garments to suit users with different needs
- Their upper limb range includes: ReadyWrap, ExoStrong and Caresia

To Order

Choose your preferred model, then use the appropriate Size Guide to select the appropriate size and code

ReadyWrap (Inelastic Compression)

- Created to replace daytime compression bandaging or serve as an alternative/supplement to compression garments
- Inelastic materials produce high working pressure and low resting pressure
- Breathable and wickable materials ensure maximum comfort

Solaris - ReadyWrap - Gauntlet

- Provides coverage from the wrist to the MCP (quarter finger) and overlaps the ReadyWrap arm garment to reduce gapping
- Left / Right specific design contours more naturally to the hand
- Pocket design allows for swell spots or foam pads to provide additional support

					DVA: AR2	2	
		Size Guide			Co	Code	
Size	Length (cm)	Wrist (cm)	Palm (cm)	Colour	Left	Right	
Creatil	up to 10 F	10 5 17	10 00	Beige	FBS303663	FBS303671	
Small up to 1	up to 12.5	13.5 - 17	18 - 20	Black	FBS303664	FBS303672	
Medium	up to 13.5	15 5 20	5.5 - 20 20 - 22	Beige	FBS303661	FBS303669	
Medium	up to 13.5 15.5 - 20	15.5 - 20		Black	FBS303662	FBS303670	
Large	up to 14 5	175 22	22 - 24	Beige	FBS303659	FBS303667	
Large	up to 14.5 17.5 - 23	22 - 24	Black	FBS303660	FBS303668		
XLarge	X Large up to 15.5 19.5 - 26 24 - 26	24 - 26	Beige	FBS303665	FBS303673		
A Large		24 - 20	Black	FBS303666	FBS303674		

Measurement: All measurements are circumferences

- Ch





SOLARIS

DVA: AR22





Solaris - ReadyWrap - Arm

- Double-sewn liner allows for easy, one-handed donning and doffing
- Bilateral design contours to both left and right arms
- Features a padded wrist and elbow to eliminate compression hot spots and provide additional comfort

Size Guide			Short (Arm length up to 43cm)	Average (Arm length up to 46cm)	Tall (Arm length up to 49.5cm)				
Size	Wrist (cm)	Elbow (cm)	Axilla (cm)	Colour		Code			
Care ell	10 5 17	10 00	00 00	Beige	FBS303687	FBS303689	FBS303691		
Small	13.5 - 17 19 - 26	19 - 26	19 - 26 22 - 29	Black	FBS303688	FBS303690	FBS303692		
N/ a all a sea	155 00			04.00	00.00	Beige	FBS303681	FBS303683	FBS303685
Medium	15.5 - 20	24 - 33	24 - 33 26 - 36	Black	FBS303682	FBS303684	FBS303686		
1	17 5 00	00 40	01 40	Beige	FBS303675	FBS303677	FBS303679		
Large	Large 17.5 - 23 29 - 40	29 - 40	31 - 43	Black	FBS303676	FBS303678	FBS303680		
VIerre	17 5 00	20 42		Beige	FBS303693	FBS303695	FBS303697		
X Large 17.5 - 23	30 - 43	36 - 50	Black	FBS303694	FBS303696	FBS303698			

Measurement: All measurements are circumferences

Solaris - ExoStrong (20 - 30 mmHg)

- ExoStrong garments are designed to move with you. With a functional design at the elbow, wrist and web spaces, you can enjoy the freedom to move
- ExoStrong provides **20 30mmHg** of compression
- Available in beige

Solaris - ExoStrong - Arm



Size Guide				Short (Arm length up to 43cm)	Average (Arm length up to 46cm)	Tall (Arm length up to 49cm)	
Size	Wrist (cm)	Elbow (cm)	Upper Arm (cm)	Style		Code	
Carroll	45 47	04 07	07 00	Knit Top	FBS303564	FBS303572	FBS303580
Small	Small 15 - 17 24 - 27	24 - 27	27 - 30	Silicone Top	FBS303565	FBS303573	FBS303581
Madium	47 40	27 - 31	01 05	Knit Top	FBS303566	FBS303574	FBS303582
Medium	17 - 19		31 - 35	Silicone Top	FBS303567	FBS303575	FBS303583
1	00 00	04 05	25 40	Knit Top	FBS303568	FBS303576	FBS303584
Large	Large 20 - 22	31 - 35	35 - 40	Silicone Top	FBS303569	FBS303577	FBS303585
VIerre	22 25			Knit Top	FBS303570	FBS303578	FBS303586
X Large 23 - 25	36 - 41	40 - 45	Silicone Top	FBS303571	FBS303579	FBS303587	

Measurement: All measurements are circumferences

Solaris - ExoStrong - Glove & Gauntlet

	Size Guide		Glove (to DIPs)	Gauntlet (to MCPs)
Size	Palm (cm)	Wrist (cm)	Cc	de
Small	17 - 19	15 - 17	FBS303558S	FBS303559S
Medium	19 - 21	17 - 19	FBS303558M	FBS303559M
Large	22 - 24	20 - 22	FBS303558L	FBS303559L
X Large	25 - 27	23 - 25	FBS303558XL	FBS303559XL



Measurement: All measurements are circumferences

General Enquiries

*DIPs - Half finger

*MCPs - Quarter Finger



Solaris Caresia - Lymphoedema Compression Bandage Liners

- Caresia standard sized bandage liners simplify complex multi-layered bandaging • protocols by replacing and improving upon foam liners, stockinette, and gauze
- Caresia over-wrapping is more comfortable and effective than standard bandaging • during the active phase of therapy
- The unique foam chips used in every Caresia garment help resolve tissue induration, • stimulate interstitial fluid movement, and diffuse irregular wrapping pressures

Glove

(DVA:		
	Size Guio		
Size	Palm (cm)	Wrist (cm)	Code
S	15 - 19	13.5 - 16	FBS303590
Μ	19 - 23	16 - 22	FBS303589
L	23 - 28	22 - 28	FBS303588

Measurement: All measurements are circumferences

Arm - MCPs to Axilla

	DVA: AR22				
	S				
Size	Palm (cm)	Wrist (cm)	Biceps (cm)	Co	de
0126		What (Chi)	Biceps (cm)	Left	Right
S	15 - 19	13.5 - 16	20 - 30	PTS642634	PTS642637
Μ	19 - 23	16 - 22	30 - 40	PTS642633	PTS642636
L	23 - 28	22 - 28	40 - 50	PTS642632	PTS642635

Gauntlet

(DVA:	1	
	Size Guio		
Size	Palm (cm)	Wrist (cm)	Code
S	15 - 19	13.5 - 16	PTS642631
Μ	19 - 23	16 - 22	PTS642630
L	23 - 28	22 - 28	PTS642629

Measurement: All measurements are circumferences

Arm - Wrist to Axilla

(DVA:		
	Size Guio		
Size	Wrist (cm)	Biceps (cm)	Code
S	13.5 - 16	20 - 30	PTS642638
Μ	16 - 22	30 - 40	PTS642640
L	22 - 28	40 - 50	PTS642639

Measurement: All measurements are circumferences *MCPs - Quarter Finger

Measurement: All measurements are circumferences

Solaris Swell Spots - Lymphoedema Pressure Pads

- Designed to be used under compression bandaging, to address areas of • swelling, pad bony prominences, soften scars, or fill crevices
- Tailored (breast, foot, ankle, patella, thigh etc), and multi-function types are • available to be used wherever is appropriate

Solaris Swell Spots	Dimensions	Code	Solaris Swell Spots	Dimensions	Code
Swell Spot - Curved	5 x 12cm	PTS642233	Swell Spot - Thora - Thigh	50 x 22cm	PTS642244
Swell Spot - Dorsal - (S)	9 x 14.5cm	PTS642234S	Swell Spot - Crescent - (S)	54(L) x 6.5(H) x 9(D)cm	PTS642245S
Swell Spot - Dorsal - (L)	9 x 18cm	PTS642234L	Swell Spot - Crescent - (M)	60(L) x 6.5(H) x 11(D)cm	PTS642245M
Swell Spot - Palmar - (S)	8.5 x 7cm	PTS642235S	Swell Spot - Crescent - (L)	71(L) x 6.5(H) x 12(D)cm	PTS642245L
Swell Spot - Palmar - (L)	9.5 x 8cm	PTS642235L	Swell Spot - We'a Toe - (S)	Up to 22cm (MTP Circ)	PTS642246S
Swell Spot - Oval	8 x 19cm	PTS642236	Swell Spot - We'a Toe - (M)	Up to 25cm (MTP Circ)	PTS642246M
Swell Spot - Katina - (S)	10 x 38cm	PTS642237	Swell Spot - We'a Toe - (L)	Up to 32cm (MTP Circ)	PTS642246L
Swell Spot - Katina - (L)	18 x 42cm	PTS642238	Swell Spot - Female Labia	30 x 50cm	PTS642247
Swell Spot - Sobol - (S)	18cm (diam)	PTS642239S	Swell Spot - Male Scrotum	30 X 50cm	PTS642248
Swell Spot - Sobol - (L)	23cm (diam)	PTS642239L	Swell Spot - Lateral Bra	10 x 16cm	PTS642249
Swell Spot - Jo-Ju	5 x 8cm	PTS642240	Swell Spot - LT Knee Shaper	35 x 65cm (Patella Girth)	PTS642250
Swell Spot - Kimbe	25 x 15cm	PTS642241	Swell Spot - Sweet Spot	7 x 14cm	PTS642251
Swell Spot - Breast - (S)	A - B cup	PTS642242S	Swell Spot - Tic Tac - (S)	18 x 18cm	PTS642252S
Swell Spot - Breast - (M)	C - D cup	PTS642242M	Swell Spot - Tic Tac - (M)	28 x 28cm	PTS642252M
Swell Spot - Breast - (L)	D - DD cup	PTS642242L	Swell Spot - Tic Tac - (L)	38 x 38cm	PTS642252L
Swell Spot - Versi - Face	26 x 17cm	PTS642243			

DVA: AR22

DVA Contact Details





Solaris Tribute - Night

- Tribute custom-designed garments are essential for a successful lymphoedema home management program
- Worn nightly, Tribute garments eliminate the need for complex bandaging, and offer an effective solution to help regulate fluid fluctuations, decrease tissue turgor, assist a weakened lymphatic system, and maintain gains made in therapy
- Foam-filled channels in the garments create localised tissue stretch and pressure differential, opening the initial lymphatic gaps to encourage re-absorption of interstitial fluids and particles
- **Note:** Custom garments require DVA prior approval on a case-by-case basis

Flowtron 'Hydroven 3' Compression Pump

Application

- The 'Hydroven 3' intermittent pneumatic compression system treats a wide variety of vascular and lymphatic conditions, including venous ulceration, dependent and traumatic oedema and chronic venous insufficiency
- Client should be assessed by a qualified medical professional before using this product.

Features

- Clinically effective, non-invasive dynamic compression
- Efficient use in both the clinic or home environment
- Three chamber leg or arm compression
- Adjustable pressure control provides tailored therapy to suit individual patients needs
- Pressure range of 30-110 mmHg
- Inflation/deflation cycle time of 90 seconds each
- Secure snap-lock connection of all garments to the pump promotes easy set up and prevents accidental disconnection
- Single tube attachment to each garment ensures less clutter

Garments

- Hydroven's triple-chamber arm and leg garments provide gradient sequential compression, with 10% reduction in pressure in each chamber distally to proximally
- All garments can be used on either left or right limb, or two units can be used concurrently on both limbs

VenaPro Post-Op DVT Leg Pump

Application

- A completely portable device which helps prevent postoperative DVT for moderate to high risk patients
- Client should be assessed by a qualified medical professional before using this product.

Features

- Product comprises of two wrap-around leg units, one for each leg
- Each leg unit has an attached re-chargeable battery (5-6 hours)
- No tubes or cords simple one button operation weighs less than 0.5kg
- Graduated and asymmetric compression applies 50mmHg once per minute





Solaris Tribute - Night	Code
Fingertips to Axilla (20 mmHg)	PTS642223
Toes to Knee (20 mmHg)	PTS642224

hydroven

	DVA: A	R23		
Product	Code			
'Hydroven 3' Pump	MDC402100	MDC402100		
	DVA: A	R26		
Garments	Limb Garment	Inserts		
Half leg 50cm	MDC402118	MDC402100		
Full leg 66cm	MDC402119	MDC402110		
Full leg 71cm	MDC402120	MDC402111		
Full leg 76cm	MDC402121	MDC402112		
Full leg 84cm	MDC402122	MDC402113		
Full leg 92cm	MDC402123	MDC402114		
Full arm 68cm	MDC402124	MDC402115		
Full arm 78cm	MDC402125	MDC402116		







DVA: AR	23
Product	Code
VenaPro Post-Op DVT Leg Pump	PTS642655

Pressure Care & Circula