Distributed by:



SOLUTIONS FOR HUNTINGTON'S DISEASE

www.aidacare.com.au 1300 133 120



HUNTINGTON'S BED - ELECTRIC

BFB046010

Made in Germany, the Careline Huntingtons bed is of the highest quality. The beds are sturdy in construction and built to resist constant movements by the user. The mattress is also waterproof and designed to reduce the risk of shearing.

The bed offers high low functions, backrest and leg support in addition to being fully padded for the protection of the user. This padding extends to the head and footboards and also the rails are padded.

The padded rails are electrically operated to minimize risk to users and also carers. The castors are lockable by a foot control.

Suitable For:

- Huntington's Disease
- **Cerebral Palsy**
- **Multiple Sclerosis**



- **Restless Legs** People needing padded
- protection from thrashing







430 - 610 mm

900 - 1100mm



0 - 40°

2000 x 900 mm



1300 133 120



HUNTINGTON'S BED - MANUAL

BEB046000H

The Manual Huntington's Bed is a fully adjustable hospital bed with manually controlled padded protective side rails.

Suitable For:

- Huntington's Disease
- Cerebral Palsy
- Multiple Sclerosis
- Restless Legs
- People needing padded protection from thrashing





FEATURES:



www.aidacare.com.au





ATLANTA CHAIR

CHP197820











A robust and durable chair which provides a safe and comfortable seat for those with vigorous involuntary movements while reducing the risk of falling from the chair.

The Atlanta[™] is necessary to comfortably accommodate individuals who may be challenging to seat. It is a unique in that it provides security without compromising function, pressure or postural management.

The integrated tilt and recline mechanism is fundamental in positioning and will reduce the risk of sliding, pressure damage and hazardous postural tendencies like posterior pelvic tilt. The resultant deep seat and high arms is beneficial for those who may be at risk of falling from the chair.

The Atlanta[™] is successfully used in the management of lower limb oedema. The variable angle leg rest used in combination with the integrated tilt and recline helps to elevate the feet higher than the heart which is the optimum position for fluid drainage. Individuals with tight hamstrings and contracted legs can also be accommodated due to the negative angle clearance underneath the seat which allows the feet to retract behind the knee without compromising the persons position in the chair.

The Atlanta[™] is effective for adults with the following disabilities and conditions:

- Huntington's disease •
- Alzheimer's disease
- Lower limb oedema
- Cellulitis
- Advanced dementia
- Parkinson's disease
- **Brain Injury**

The Atlanta[™] has been proven effective fin seating clients with the following postural and medical conditions:

- Lower limb oedema and swelling
- **Contracted legs** •
- Risk of pressure uclers •
- Vigorous involuntary movement
- Agitation
- Risk of falling or sliding from chairs •

SEAT DEPTH 480 mm



OVERALL HEIGHT



750 mn

SEAT TO FLOOR HEIGHT



SAFE WORKING LOAD

ABR45.5 210518



380 - 500 mm

1135 mn

540 mm

1300 133 120

www.aidacare.com.au